

DRAFT Guidance Notes on the Adaptation of the In Vitro micronucleus assay (OECD TG 487) for Testing of Manufactured Nanomaterials

Study Report reflecting State of the Art Knowledge for the in vitro micronucleus assay

Enter any logistical information related to the meeting e.g. meeting date, time and location.

DRAFT circulated to the WNT for inputs by **31st December 2021**

Mar Gonzalez: mar.gonzalez@oecd.org

3 Table of contents

4	Guidance Notes on the Adaptation of the In Vitro micronucleus assay (OECD TG 487) for	
5	Testing of Manufactured Nanomaterials	4
6	Aim and scope of the project	4
7	Phase 2: Methods	7
8	Preparation of ENMs.	7
9	Phase 2: Results	12
10	Phase 2: SU TK6 Results	13
11	SU HepG2 Results	15
12	BASF Buffy coat Results	17
13	BASF whole blood Results	18
14	BASF V79 Results	19
15	SU V79 Results following BASF SOP	20
16	Exposure assessment	21
17	Recommendations	23
18	Limitations	25
19	Additional Supporting Interlaboratory Trial	27
20	Annex 1 – Harmonised SOP generated by BASF & SU	28
21	Experiments for characterization of nanoparticles, recommended by both BASF & Swansea	
22	University	28
23	<i>In vitro</i> micronucleus test on primary human lymphocytes; Cytochalasin B method with	
24	nanoparticles	BASF
25	Standard Operating Procedure 2020	30
26	Preparation of the cells	31
27	Stimulation of the cells with phytohemagglutinin	31
28	Test substance preparation and preparation of the particulate positive control	31
29	The test substance is tested with the following concentrations:	32
30	Preparation and treatment of test cultures	32
31	Checking or determining further parameters	33
32	Test substance removal	34
33	Cytochalasin B treatment	34
34	Cell harvest and preparation of slides	34

35	Micronucleus analysis	36
36	The following criteria must be evaluated for the assessment as a micronucleus	36
37	Materials and chemicals	37
38	Recipe	38

39	<i>In vitro</i> cytokinesis-blocked micronucleus (CBMN) test protocol for use with TK6 cells and engineered nanomaterials (ENMs)	Swansea	
40		University	
41		Standard Operating Procedure 2020	39
42	Biological Setting & Test System		40
43	Chemicals & Reagents		40
44	Cell culture preparation		40
45	Procedure		40
46	MMC CAS 50-07-7 (Fisher Scientific; BP2531-2), dose for 1 – 1.5 cell cycles at 0.01µg/ml		41
47	Preparation of ENMs		41
48	Cytotoxicity Calculations		45
49	References		47

50

51 FIGURES

52	Figure 1. TK6 cytotoxicity & genotoxicity data following 1 cell-cycle exposure to 5nm AuNPs (A), 30nm AuNPs (B) and 22nm SiO ₂ NPs (C). MMC (0.01µg/ml) was used as the positive chemical control, WC-Co was used as the particulate control at 20 and 100µg/ml. Results were considered statistically significant (*) when p<0.05.	13
53		
54	Figure 2. HepG2 cytotoxicity & genotoxicity data following 1 cell-cycle exposure to 5nm AuNPs (A), 30nm AuNPs (B) and 22nm SiO ₂ NPs (C). MMC (0.01µg/ml) was used as the positive chemical control, WC-Co was used as the particulate control at 20 and 100µg/ml. Results were considered statistically significant (*) when p<0.05	15
55		
56	Figure 3. Buffy coat CBPI and frequency of micronuclei determined by the <i>in vitro</i> CBMN assay following 1 cell-cycle exposure to 5nm AuNPs (A), 30nm AuNPs (B), CeO ₂ NPs (C), 22nm SiO ₂ NPs (D). Positive controls used were MMC at 0.04µg/ml, Colcemide at 0.05µg/ml and WC/Co at 10, 30, 60 & 100µg/ml. Results were considered statistically significant (*) when p≤0.05 or (**) when p≤0.01.	17
57		
58	Figure 4. Whole blood CBPI and frequency of micronuclei determined by the <i>in vitro</i> CBMN assay following 1 cell-cycle exposure to 5nm AuNPs (A), 30nm AuNPs (B), BaSO ₄ NPs (C), 22nm SiO ₂ NPs (D). Positive controls used were MMC at 0.04 µg/ml, Colcemide at 0.05 µg/mL and WC/Co at 10, 30, 60 & 100µg/ml. Results were considered statistically significant (*) when p≤0.05 or (**) when p≤0.01.	18
59		
60	Figure 5. V79 CBPI and frequency of micronuclei determined by the <i>in vitro</i> CBMN assay following 1 cell-cycle exposure to 5nm AuNPs (A), 30nm AuNPs (B), CeO ₂ NPs (C), 22nm SiO ₂ NPs (D), BaSO ₄ NPs (E). Positive controls used were EMS 600 µg/ml and WC/Co at 60 & 100µg/ml. Results were considered statistically significant (*) when p≤0.05 or (**) when p≤0.01.	19
61		
62	Figure 6. V79 cytostasis & genotoxicity data following 1 cell-cycle exposure to 5nm AuNPs (A), 30nm AuNPs (B) and 22nm SiO ₂ NPs (C). EMS (600µg/ml) was used as the positive chemical control, WC-Co was used as the particulate control at 30 and 100µg/ml. Results were considered statistically significant (*) when p<0.05.	20
63		
64		
65		
66		
67		
68		
69		
70		
71		
72		

73

74 No table of figures entries found.

75

Guidance Notes on the Adaptation of the In Vitro micronucleus assay (OECD TG 487) for Testing of Manufactured Nanomaterials

76 Aim and scope of the project

77 1. The OECD Working Party on Manufactured Nanomaterials addresses all the components needed for
78 thorough hazard assessment of nanomaterials (NM) for human health and the environment, including making
79 available reliable test methods for the purpose of NMs safety assessment. At an OECD expert meeting held in
80 Paris in October 2014 it was concluded that the scientific information available was not sufficient to fully support
81 development of a harmonised version of the *in vitro* micronucleus test protocol because the existing or published
82 data were not generated in a harmonised way and did not provide answers to some key questions (e.g. the
83 uptake on NMs, dosing regimens, necessary length of exposure to NMs or scheme of treatment to cytochalasin
84 B). The OECD WPMN approved a JRC led project with the aim to identify needed adaptations of the *in vitro*
85 micronucleus test (Test Guideline 487) for the assessment of manufactured NMs and eventually provide
86 appropriate guidance. The same project was subsequently included by the OECD-WNT in the work plan of the
87 Test Guidelines Programme (Project 4.95). Therefore, work aiming to develop a Guidance Document (GD)
88 addressing the necessary adaptations of current genotoxicity TGs for nanomaterials testing was initiated in
89 January 2019 under the coordination of JRC. In order to share the lessons learnt from work carried out so far
90 under this project and guide the users of the Test Guideline 487 when testing nanomaterials, it was agreed to
91 publish this document, mindful that it will be refined in future with the testing of additional nanomaterials.

92 2. The whole study was divided into 3 phases:

93 **Phase 1: Clarification of main technical issues**

- 94 • Physical-chemical characterization of gold, silver and silica Nanoparticles.
- 95 • Cytotoxicity assessments of nanoparticles in A549, Caco-2, V79, TK6 cells
- 96 • Uptake of Nanoparticles into the various cell systems
- 97 • Define the doubling times

98 3. JRC coordinated this phase of work, which was completed and culminated in the following reports:

99 4. The Report on Physico-Chemical Characterisation of the nanomaterials chosen for the study was
100 published in 2018 by JRC.

- 101 ○ Drewes, C., Ojea Jimenez, I., Mèhn, D., Colpo, P., Gioria, S., Bogni, A., Ponti, J., Kinsner-
 102 Ovaskainen, A., Gilliland, D. and Riego Sintes, J., Physicochemical characterisation of
 103 gold, silica and silver nanoparticles in water and in serum-containing cell culture media,
 104 EUR 29054 EN, Publications Office of the European Union, Luxembourg, 2018, ISBN 978-
 105 92-79-77705-9 (online),978-92-79-77704-2 (print), doi:10.2760/818663
 106 (online),10.2760/58721 (print), JRC110379.

107 5. The Report on In vitro toxicity and uptake with the selected cell lines and selected nanomaterials was
 108 published on 2020 by JRC

- 109 ○ Bogni A., Ponti J., Drewes C., Kinsner-Ovaskainen A., Bremer-Hoffmann S., In vitro
 110 cytotoxicity and cellular uptake evaluation of gold, silica and silver nanoparticles in five
 111 different cell lines: Caco-2, A549, CHO, V79 and TK6, European Commission, Ispra, 2020,
 112 JRC120791

113 **Phase 2: Interlaboratory comparison**

- 114 • Selection, purchase and characterization of material
- 115 • Protocol development by lead laboratory
- 116 • Protocol transfer to other laboratories
- 117 • Interlaboratory study
- 118 • Study evaluation

119 6. The JRC coordinated this phase of work, while the experimental work was performed by BASF
 120 (Germany) and Swansea University (SU; UK). The outcomes of the experimental work conducted under this
 121 phase of the work programme are detailed within this document. Additionally, a publication detailing this phase
 122 of the study is in preparation.

123 7. The selection of cell lines to evaluate during this phase of study were agreed upon during an expert
 124 group meeting and included TK6 cells, V79 and human lymphocytes. These cells were selected as they are
 125 cited as appropriate for use with the *in vitro* micronucleus assay in OECD TG487, while cells such as A549 and
 126 Caco-2 (which were considered in Phase 1) have not been evaluated for use with the assay.

127 **Phase 3: Guidance development**

- 128 • Set up of expert drafting group
- 129 • Guidance development
- 130 • Presentation to OECD

131 8. The JRC was not able to coordinate this phase of work and so the leadership was transferred to the UK
 132 and Germany.

133 9. The OECD guideline 487 (adopted 29. July, 2016) is the current guideline under the Mutual Acceptance
 134 of Data (MAD) agreement that instructs how an *in vitro* micronucleus assay should be performed. The
 135 procedures described in the guideline are, however, directed at testing chemicals and not nanomaterials. In the
 136 guideline this deficiency is noted, yet a recommendation of a suitable adaptation of the protocol for the
 137 nanomaterial testing is not given. Thus, the current regulatory situation for the mutagenicity testing of
 138 nanomaterials is officially to follow the current procedure of the guideline, which is not necessarily the optimal
 139 procedure from a scientific point of view. Parameters such as the treatment interval or selection of the test
 140 concentration are a few examples of aspects, which need specific consideration when testing nanomaterials; a
 141 pulse treatment of cultures with nanomaterials do not provide sufficient time for the internalization process of a

142 particular material and nanomaterials can per se only be tested at precipitating levels. Since there is no official
143 document on how to adapt the OECD 487 protocol for nanomaterial testing, there is a certain level of urgency
144 in providing a guidance on the necessary adaptations, which would allow a more relevant assessment of the
145 mutagenic potential of nanomaterials and avoid repetition of inadequately performed studies. This preliminary
146 guidance document is therefore aimed at addressing this urgent deficiency in the guidance available.

Phase 2: Methods

147 Preparation of ENMs.

148 10. 5nm AuNPs (#OECD-RLS-03) & 30nm AuNPs (#OECD-RLS-10)), 22nm SiO₂ NPs (#OECD-RLS-0102)
149 were supplied by JRC. JRC stock ENMs were sonicated in a 90W Ultrasonic Bath (Fisher Scientific #FB15046)
150 for 20 minutes at 37°C to encourage destabilisation of agglomerate material. Following sonication, working
151 concentrations of particles were prepared in cell culture media using a 1:1 dilution for the highest concentration
152 and using serial dilutions to prepare the rest of the dose range (Supplementary Information Table 1). Given
153 each JRC ENM stock were at different concentrations, this process was ENM-specific. CeO₂ (Frauenhofer,
154 Lot#NM-212 Reference Nanomaterial), BaSO₄ (Solvay, Lot#V106) and Tungsten-Carbide/Cobalt, (8wt%
155 WC/Co <200nm, 99.5% LOT#5561-072018, Nanostructured & Amorphous Materials Inc., USA), was used as
156 a trial positive particle control and was weighed, suspended and sonicated according to the NANoREG protocol
157 (Alstrup Jensen et al., 2014).

158 **Cell culture (University of Swansea protocol).**

160 11. The human lymphoblastoid (TK6) cells were purchased from Public Health England (PHE). The cells
161 were cultured in RPMI 1640 supplemented with 10% Horse serum and 1% L-glutamine. TK6 cells were routinely
162 sub-cultured for 2 weeks prior to ENM testing, cells were regularly checked for potential changes to morphology
163 and density by light microscopy. The human liver cancer (HepG2) cell line was purchased from American Type
164 Culture Collection (ATCC). The cells were cultured in DMEM containing 10% FBS, 1% Penicillin/Streptomycin
165 (P/S) and 1%L-glutamine. HepG2 cells were grown to 80% confluency before being routinely sub-cultured. The
166 V79 cells (AbbVie GmbH) were cultured in MEM Eagle media (Pan-Biotech, UK) containing 1% L-glutamine,
167 1% Amphotericin, 10%FBS and 1% P/S. V79 cells could grow to 80% confluency before being routinely sub-
168 cultured.

169 **Semi-automated in vitro cytokinesis-blocked micronucleus assay.**

171 12. Two slightly different approaches were taken with the CBMN assay depending on whether suspension
172 or adherent cell lines were used. TK6 (suspension cells), HepG2 and V79 (adherent) cells were seeded at
173 1.0×10^5 cells/ml in T25 flasks along with satellite flasks per concentration to be counted for cytotoxicity and were
174 then treated with ENMs for 1 cell-cycle (12 hours). Mitomycin-C (MMC) at 0.01µg/ml was used as the chemical
175 positive control, and WC/Co at 20µg/ml & 100µg/ml was tested as a potential positive particle control. The *in*
176 *vitro* cytokinesis-blocked micronucleus assay was performed as described previously by Evans *et al* & Burgum
177 and colleagues (Evans et al., 2019, Burgum et al., 2020). On the day of exposures, cells are counted at least 2
178 hours before exposures (with suspension cells this can be done immediately prior to exposures). When using
179 the HepG2 cells time is required for them to re-adhere to the flasks before exposures. Cells were dosed (both
180 the satellite flasks and the CBMN-flasks) with ENMs prepared in cell-culture media and the negative control
181 being media only. Exposures were performed as close to sonication times as possible to avoid sedimentation.
182 Following the exposure period, the satellite flask cells were counted for the calculation of relative population
183 doubling (RPD), a measure of cytotoxicity (Equation 1). Cell counts were performed using a Beckman coulter

184 counter by adding 100µl of cells to 10ml of diluent. Following ENM exposure for 1 – 1.5 cell cycle, TK6 cells
185 were centrifuged and washed with PBS in triplicate (removal of as much ENM as possible), satellite flasks were
186 counted to facilitate calculation of cytotoxicity data. V79 and HepG2 cells had exposure media aspirated and
187 were then washed twice with PBS before being counted in the same manner. TK6 cells were re-seeded into
188 clean T25 flasks in cell-culture media containing 3µg/ml of cytochalasin-B (cyto-B), the cells were returned to
189 the incubator for a further 1.5 cell cycles. V79 and HepG2 cells remained in the same flasks and were
190 supplemented with fresh media containing the same concentration of cyto-B.

191 **Cell culture (BASF SE Protocol):**

192 13. **V79 cell cultures:** V79 (adherent) cells were seeded at 5.0×10^5 cells/ml in T25 flasks (attachment
193 period 24-28 hours) and were treated with ENMs along with satellite flasks per concentration to be counted for
194 cytotoxicity and LA ICP MS assessment for 1 cell-cycle (14 hours). Ethyl methanesulfonate (EMS) at 500µg/ml
195 and 600µg/mL was used as the chemical positive control, and WC/Co at 30 µg/mL & 100µg/ml was tested as a
196 potential positive particle control. After 1 cell-cycle test substance treatment the cultures were rinsed twice with
197 HBSS (Hanks Balanced Salt Solution). The cultures intended for mutagenicity assessment were incubated in
198 MEM (incl. 10% [v/v] FCS) supplemented with CytB (final concentration: 3 µg/mL; stock: 0.6 mg/mL in DMSO;
199 AppliChem, Cat.No. A7657) for 24 hours. Cultures used for LA ICP MS assessments were trypsinized, fixed
200 twice in methanol : acetic acid (19:1; -20 °C) and spread on slides.

201 14. **Whole blood cultures:** After 48 hours the activated cell cultures were pooled and centrifugated in 10
202 mL aliquotes at 900 g for 5 minutes. After centrifugation the supernatant (culture medium) was removed and
203 the cells suspended in ENM dilutions in culture medium along with satellite tubes for each concentration for LA
204 ICP MS assessment. All tubes were transferred into cell culture flasks and incubated for 20 hours. Mitomycin C
205 (MMC; Roche Diagnostics) at 0.04 µg/mL and Colchicine (Col; Roche Diagnostics) at 0.05 µg/mL were used as
206 the chemical positive controls, and WC/Co at 10, 30, 60 and 100µg/mL was tested as potential positive particle
207 control. At the end of the exposure period, the cells were transferred in tubes, centrifuged for 5 minutes at 900
208 g and resuspended in HBSS (Hanks Balanced Salt Solution). Washing of the cells was repeated at least once.
209 Then the cells were centrifuged at (900 g, 5 min) and resuspended in DMEM/F12 medium with 10% [v/v] FCS
210 and transferred into 25 cm² cell culture flasks. Cyt B (see 3.2.2) was added to the cultures intended for
211 mutagenicity assessment and incubated at 37°C, 5% (v/v) CO₂ and ≥ 90% relative humidity for 20 hours. To
212 prepare the cells for the LA ICP MS assessments, they were separated from the nanoparticles using density
213 centrifugation (Ficoll paque) and washed once at 900 g for 5 minutes. The obtained cells were fixed twice with
214 methanol:acetic acid (19:1; -20°C) and spread on slides.

215 15. **Buffy coat cell cultures:** After 48 hours the activated cell cultures were pooled and centrifugated in 10
216 mL aliquotes at 900 g for 5 minutes. After centrifugation the supernatant (culture medium) was removed and
217 the cells suspended in ENMs dilution in culture medium along with satellite tube for each concentration for LA
218 ICP MS assessment. All tubes were transferred into cell culture flasks and incubated under agitation for 20
219 hours (corresponding to 1 cell cycle) . Mitomycin C (MMC; Roche Diagnostics) at 0.04 µg/mL and Colchicine
220 (Col; Roche Diagnostics) at 0.05 µg/mL were used as the chemical positive controls, and WC/Co at 10, 25 30,
221 60 and 100µg/mL was tested as potential positive particle control. At the end of the exposure period, the cells
222 were transferred in tubes, centrifuged for 5 minutes at 900 g and resuspended in HBSS (Hanks Balanced Salt
223 Solution). Washing of the cells was repeated at least once. Then the cells were centrifuged at (900 g, 5 min)
224 and resuspended in RPMI medium with 20% [v/v] FCS and transferred into 25 cm² cell culture flasks. Cyt B was
225 added to the cultures intended for mutagenicity assessment and incubated at 37°C, 5% (v/v) CO₂ and ≥ 90%
226 relative humidity for 20 hours. In order to prepare the cells for the LA ICP MS assessments, they were separated
227 from the nanoparticles using density centrifugation (Ficoll paque) and washed once at 900 g for 5 minutes. The
228 obtained cells were fixed twice with methanol:acetic acid (19:1; -20 °C) and spread on slides.

229

230 **Cell harvesting for semi-automated Metafer scoring (University of Swansea Protocol).**

231 16. The TK6, HepG2 and V79 cells were harvested by centrifugation (230g for 5 min), resuspended in 5ml
 232 of pre-warmed PBS and centrifuged at 230g for 10 minutes. The supernatant was discarded, this was repeated
 233 a second time. The cells were then resuspended in hypotonic solution (KCl 0.56%), before being centrifuged
 234 immediately at 230g for 10 minutes. The cells were resuspended in Fixative 1 (methanol: acetic acid: NaCl
 235 (0.09%) (5:1:6 parts)) and incubated at 4°C for 10 minutes before centrifugation at 230g for 10 minutes. Cells
 236 were resuspended in Fixative 2 (methanol: acetic acid (5:1 parts)) and incubated at 4°C for 10 minutes before
 237 centrifugation at 230g for 10 minutes, this was repeated 3 times. Cells can be maintained overnight in Fixative
 238 2 at 4°C, tubes covered by foil. The day before making slides, freshly opened microscope slides were placed in
 239 a glass tank of Fixative 2 at 4°C, for 2 hours before slide preparation (ideally this is done overnight). On the day
 240 of preparing slides, the fixative was replaced with ddH₂O. On the day of slide preparation, the fixed cell
 241 suspensions were centrifuged at 230g for 10 minutes and thoroughly re-suspended in ~1ml of Fixative 2. Slides
 242 were removed from the ddH₂O and wiped dry with slide tissue. A total of 100µl of the cell suspension was evenly
 243 pipetted onto the slide held at an angle. The slides were then stood near-vertical on tissue paper to dry. The
 244 cell density was checked to ensure cells were evenly distributed, without clumping. Once dried, the slides were
 245 stained with 30µl of Vectashield mounting medium with DAPI, coverslip applied and incubated in the dark for
 246 15 minutes. Slides were scored using the Zeiss AxioCam HRc (Carl Zeiss Microscopy and Imaging, UK) semi-
 247 automated Metafer system. The details for the classifier used to support the analysis can be found in Table 1
 248 of the Supplementary Information. All experiments were performed in triplicate ($n=3$) and 2000 binucleated (BN)
 249 cells per replicate were scored per concentration (6000 BN cells in total).

250

251 **Cell harvesting for manual scoring approach (University of Swansea Protocol).**

252 17. Cells were harvested by centrifugation (230g for 5 min), resuspended in 5 ml of pre-warmed PBS and
 253 centrifuged at 230g for 5 minutes. The supernatant was discarded, and this step was repeated once more before
 254 the cells were resuspended in 5ml of PBS. The cell density was checked at this stage using a practice slide.
 255 Labelled slides were 90% methanol-cleaned, dried and placed onto cyto-clips with 5mm cyto-dot filter cards
 256 and cyto-funnel attached. A total 100 µl of cell suspension was pipetted into each funnel and centrifuged at
 257 200g for 10 minutes in a Shandon Cytospin. Slides were then examined for density. Cells were then fixed by
 258 immersing slides in 90% methanol at -20°C for 10 minutes. Once dried the slides were stained in 20% Giemsa
 259 solution Gurr buffer (pH 6.8), filtered for 10 minutes. The slides were then rinsed in pH 6.8 Gurr buffer and
 260 soaked in pH 6.8 phosphatase buffer for 1-2 minutes. Slides were dried by standing them vertically. The
 261 following stages were performed in a fume hood. Slides were dipped in xylene for 10 seconds and dried. One
 262 drop of DPX was placed on the cytodot and a 22mm x 22mm coverslip applied, the slides were then left to dry
 263 for 24 hours. Slides were viewed under a light microscope, 1000BN cells per replicate were scored per
 264 concentration (3000 BN total). Cells were scored to calculate Cytokinesis Black Proliferation Index (CBPI) as a
 265 measure of cytostasis (Equation 2).

266

Equation 1)

267

No. of population doublings in treated cultures

268

RPD= -----x 100

269

No. of population doublings in control cultures

270

Where population doubling = [log (Post-treatment cell number/ Initial cell number)]/ log 2

271

272

273

274

Equation 2)

275

% Cytotoxicity = $100 - 100 \left(\frac{\text{CBPI}_T - 1}{\text{CBPI}_C - 1} \right)$; where T is treatment, C is control

276

CBPI = No mononucleated cells + 2x BN cells + 3 x MNCs / N; N is the total number of cells scored

277

278 Cell harvesting for manual scoring approach (BASF SE Protocol).

279

18. After the CytB treatment the adherent cells were rinsed with pre-warmed HBSS and trypsinized. The obtained cell suspensions from all different cultures were treated in the same manner, namely the cells were centrifuged and the cell pellet treated with a hypotonic solution (10 minutes with 0.4% KCl for V79 and buffy coat cells; 20 minutes with 0.28% KCl for whole blood cultures). After the hypotonic treatment, the cells were fixed twice by adding fixative (19 parts methanol and 1-part acetic acid; -20°C). Slides were prepared by immersing in deionized water followed by pipetting the fixed cells on the slide. The cells were stained with May-Grünwald (3min) and 10%[v/v] Giemsa (in Titrisol, pH7.2, 20min), mounted and scored as described above.

286

287 Kinetochore staining of micronuclei.

288

19. To determine if the DNA damage was a consequence of an aneugenic or clastogenic response, kinetochore staining of micronuclei in BN cells was performed. Following the exposure period, the cells were washed three times with PBS before being cytocentrifuged (500g, 5 minutes) onto microscope slides and fixed in ice-cold 90% methanol at -20°C. Immunofluorescent staining of kinetochore proteins was performed as previously described (Singh et al., 2012, Burgum et al., 2020). Kinetochore scoring was performed on a Zeiss AxioCam HRc (Carl Zeiss Microscopy and Imaging, UK). Only the lowest and highest significant NP concentrations were used for each particle type, with WC/Co included at 100µg/ml, of which 50 micronuclei were scored for the presence or absence of FITC fluorescence in the micronucleus indicating the presence of a whole chromosome (K+) or chromosome fragment (K-) respectively (n=1 biological replicate).

297

298 Exposure assessment.

299

20. For LA-ICP-MS analysis, a 193 nm ArF Excimer Laser (NWR193 Excimer Laser Ablation System, Elemental Scientific Lasers, Bozeman, MT USA) equipped with a two-volume cell (TwoVol² Ablation Cell, Elemental Scientific Lasers) was coupled to an ICP-MS Triple Quadrupole (8900 ICP-MS Triple Quad, Agilent Technologies, Santa Clara, CA, USA). The samples were transported in a carrier gas flow (He, 800 mL/min) and introduced *via* a Dual Concentric Injector (DCI, Elemental Scientific Lasers) to the ICP-MS. An additional gas flow (Ar, 1 L/min) was added, and the sample transferred *via* a quartz injector pipe (inner diameter: 3.5 mm) into the plasma. The ICP-MS was equipped with platinum sampler and skimmer. To resolve the issue of polyatomic interferences especially for low masses (e.g. ³¹P), it was operated in TQ modus with O₂ as reaction gas. The optimum RF was set to 1300 V and the set-up was tuned daily for maximum signal intensity and an oxide ratio (*m/z* 232/248) below 1.5 % with a NIST Glass standard (NIST SRM 612, National Institute of Standards and Technology, Gaithersburg, MD, USA). Each cell was ablated separately with 50 bursts and a spot size of 25 µm. A laser pulse frequency of 100 Hz and a laser energy of 0.5 J/cm² ensured a fully ablation avoiding ablation of the glass slide. The isotope ³¹P¹⁶O⁺ was monitored with a dwell time of 50 ms regardless of the type of experiment and chosen as internal standard. For Au NPs the isotope ¹⁹⁷Au, for CeO₂ NPs the isotope ¹⁴⁰Ce¹⁶O⁺ and for WC-Co NPs the isotopes ⁵⁹Co⁺ and ¹⁸⁴W¹⁶O⁺ were detected each with a dwell time of 50 ms. To distinguish a cell event from the continuous background, a threshold three times higher than the mean signal of the background was applied. For an identified cell, an association rate was calculated by dividing the summed signal intensities of the NPs by the summed phosphorus signal intensities. This results in an association rate to what extent the NPs are allocated to the cells.

317

318

319

Statistics.

320

321

322

323

324

21. All data is presented as the mean +/- the standard deviation (SD). For data generated at the Swansea University lab: statistical analysis was performed in SPSS statistics software (v.20 IBM, UK) where all data sets were firstly analysed for normality (Shapiro-Wilk test, $p \leq 0.05$) and for equal variance $p \leq 0.05$). A one-way analysis of variance (ANOVA) was performed with post hoc Dunnett's multiple comparisons applied to evaluate pairwise statistical significance between control and concentrations; the alpha value was set to 0.05.

325

326

327

328

329

330

331

332

333

334

22. For data generated at the BASF lab: An appropriate statistical analysis was performed. The proportion of cells containing micronuclei was calculated for each test group. A comparison of the micronucleus rates of each test group with the concurrent vehicle control group was carried out for the hypothesis of equal proportions (i.e. one-sided Fisher's exact test, BASF SE). If the results of this test were statistically significant compared with the respective vehicle control ($p \leq 0.05$), labels (s) were printed in the tables. In addition, a statistical trend test (SAS procedure REG) was performed to assess a possible dose-related increase of micronucleated cells. The used model is one of the proposed models of the International Workshop on Genotoxicity Test procedures Workgroup Report. The dependent variable was the number of micronucleated cells and the independent variable was the concentration. The trend was judged as statistically significant whenever the one-sided p value (probability value) was below 0.05.

335

Phase 2: Results

336

337 23. The data generated by SU with TK6 cells revealed no cytotoxic effects following 1 cell-cycle exposure
338 to either of the AuNPs. Some significant (*) cytotoxicity was observed following the SiO₂ exposure but only at
339 the highest test concentration of 100µg/ml. Chromosomal breakage, measured by the *in vitro* CBMN assay
340 revealed some significant (*) genotoxicity following exposure to each JRC particle type, the significant
341 responses were most apparent at the highest test concentration of each material where 2-fold increases over
342 background levels (of %Mn/BN) were observed. The WC-Co particle control was observed to be positive
343 following 1 cell-cycle exposures to TK6 cells, with 2-fold increases in micronuclei induced at the two test
344 concentrations of 20 and 100µg/ml (**Figure 1**).

345 24. HepG2 cell line data generated by SU showed a more muted response to the JRC materials as
346 compared to the TK6 suspension cells. Firstly, no significant cytotoxicity was observed across any test material.
347 Secondly, the HepG2 cells showed no variation in %Mn/BN following 1 cell-cycle exposure to 30nm AuNPs.
348 Whilst the 5nm AuNPs and the SiO₂ NPs did induce significant (*) %Mn/BN, this had limited biological
349 significance as there was no increase in micronuclei that was greater than 1-fold elevation. Unlike the TK6
350 response to WC-Co particles, the HepG2 cell line only showed significant (*) %Mn/BN at the highest test
351 concentration of the material, 100µg/ml (**Figure 2**).

352 25. The BASF-generated data on buffy coat cells showed no significant cytotoxic response to the JRC
353 materials. This muted response to the materials was mirrored in the genotoxicity data whereby the only
354 significant (*) data was reported by the positive chemical controls and the WC-Co at 10, 60, and 100µg/ml
355 (**Figure 3**). Whole blood exposures to JRC nanomaterials showed largely corresponding data sets with no
356 significant (*) cytotoxic or genotoxic damage with exception of the positive chemical controls and the WC-Co
357 particle control, however only at the highest test concentration of 100µg/ml (**Figure 4**). BASF-generated data
358 on V79 cells using the JRC reference materials showed no significant (*) cytotoxicity or genotoxicity following 1
359 cell cycle exposures. The exception being the positive chemical control, MMC (**Figure 5**).

360 26. While in the above data sets, SU and BASF evaluated the response of differing cell lines to the test
361 nanoparticles, one data set was generated using the same cell line, evaluated in both labs. This final data set
362 involved exposure of the test nanoparticles to the V79 cells by SU, utilising the BASF protocol to determine the
363 transferability of the protocol and ability to generate harmonised data sets. Following the 1 cell-cycle exposures
364 to the JRC reference materials the V79 cells showed a muted response with no significant (*) cytotoxicity or
365 genotoxicity with the only exception being the positive chemical control, EMS (**Figure 6**). Thus, there was no
366 difference in the data generated by both the SU and BASF labs using this cell line; demonstrating the
367 transferability and reproducibility of this SOP.



368


369 Phase 2: SU TK6 Results

370 Figure 1. TK6 cytotoxicity & genotoxicity data following 1 cell-cycle exposure to 5nm AuNPs (A), 30nm AuNPs (B) and 22nm SiO₂NPs (C). MMC
 371 (0.01µg/ml) was used as the positive chemical control, WC-Co was used as the particulate control at 20 and 100µg/ml. Results were considered
 372 statistically significant (*) when p<0.05.

373

374





375


C


376

377 SU HepG2 Results

378 Figure 2. HepG2 cytotoxicity & genotoxicity data following 1 cell-cycle exposure to 5nm AuNPs (A), 30nm AuNPs (B) and 22nm SiO₂NPs (C). MMC
 379 (0.01µg/ml) was used as the positive chemical control, WC-Co was used as the particulate control at 20 and 100µg/ml. Results were considered
 380 statistically significant (*) when p<0.05

381






383

384

C


385 **BASF Buffy coat Results**

386 **Figure 3.** Buffy coat CBPI and frequency of micronuclei determined by the *in vitro* CBMN assay following 1 cell-cycle exposure to 5nm AuNPs (A), 30nm AuNPs (B), CeO₂ NPs (C), 22nm SiO₂ NPs (D). Positive controls used were MMC at 0.04µg/ml, Colcemide at 0.05µg/ml and WC/Co at 10, 30, 60 & 100µg/ml. Results were considered statistically significant (*) when $p \leq 0.05$ or (**) when $p \leq 0.01$.




391 **BASF whole blood Results**

392 **Figure 4. Whole blood CBPI and frequency of micronuclei determined by the in vitro CBMN assay following 1 cell-cycle exposure to 5nm AuNPs (A),**
 393 **30nm AuNPs (B), BaSO₄ NPs (C), 22nm SiO₂ NPs (D). Positive controls used were MMC at 0.04 µg/ml, Colcemide at 0.05 µg/mL and WC/Co at 10, 30,**
 394 **60 & 100µg/ml. Results were considered statistically significant (*) when $p \leq 0.05$ or (**) when $p \leq 0.01$.**



396 **BASF V79 Results**

397 **Figure 5. V79 CBPI and frequency of micronuclei determined by the in vitro CBMN assay following 1 cell-cycle exposure to 5nm AuNPs (A), 30nm**
 398 **AuNPs (B), CeO₂ NPs (C), 22nm SiO₂ NPs (D), BaSO₄ NPs (E). Positive controls used were EMS 600 µg/ml and WC/Co at 60 & 100µg/ml. Results**
 399 **were considered statistically significant (*) when p≤0.05 or (**) when p≤0.01.**






401 **SU V79 Results following BASF SOP**

402 **Figure 6. V79 cytostasis & genotoxicity data following 1 cell-cycle exposure to 5nm AuNPs (A), 30nm AuNPs (B) and 22nm SiO₂NPs (C). EMS**
 403 **(600µg/ml) was used as the positive chemical control, WC-Co was used as the particulate control at 30 and 100µg/ml. Results were considered**
 404 **statistically significant (*) when p<0.05.**

405

406





408 Exposure assessment

409 27. In addition to the cytotoxicity and genotoxicity analysis, BASF evaluated the association of the test nanoparticles with the buffy coat cells, whole blood
 410 cultures and V79 cells by LA-ICP-MS analysis. The ratio of the respective elements Au, Ce, Ba and Co to Phosphorus is given as a semi-quantitative measure
 411 for the association of the tested nanoparticles with the cells. The data given in **Table 1** show that ratio of the individual element to the intracellular phosphorus
 412 increased concentration relatedly in each individual experiment. A comparison of the ratios between the elements is not possible. However, the ratios between
 413 the cell culture type used can be compared. Thus, it can be observed that V79 cells in all cases have a higher ratio of the element to their intracellular phosphorus
 414 levels as compared to buffy coat or whole blood cells. Consequently, despite the lack of genotoxicity in V79 cells, this data demonstrates that cell association
 415 did occur following exposure to the test nanoparticles.

416

417 **Table 1.** The ratio of the indicated elements to Phosphorus in buffy coat, whole blood and V79 cells after treatment with the indicated
 418 concentrations of the given nanoparticles.

$\mu\text{g/mL}$	Buffy coat cells				Whole blood cultures				V79				
	Au _{5nm}	Au _{30nm}	CeO ₂	WC/Co	Au _{5nm}	Au _{30nm}	BaSO ₄	WC/Co	Au _{5nm}	Au _{30nm}	CeO ₂	BaSO ₄	WC/Co
	Au/P	Au/P	Ce/P	Co/P	Au/P	Au/P	Ba/P	Co/P	Au/P	Au/P	Ce/P	Ba/P	Co/P
0	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	0.02	<LOD	<LOD	<LOD	<LOD	0.03	<LOD
10	0.01	0.003	0.04	--	0.008	0.01	--	--	0.23	0.46	57.53	--	--
25	--	--	--	--	--	--	--	--	--	--	--	--	--
30	--	--	--	0.08	--	--	1.24	0.019	--	--	--	27.32	0.45
50	0.02	0.010	0.36	--	0.035	0.03	--	--	0.48	2.01	219.86	--	--
60	--	--	--	0.1	--	--	1.61	0.045	--	--	--	32.53	0.95
100	0.04	0.058	3.81	0.16	0.089	0.16	1.97	0.072	0.76	4.56	334.19	35.71	1.26

419

Recommendations

420 28. From the Phase 2 experimental work presented here, the following recommendations can be made
421 to answer the original questions posed during the expert meeting held in Paris in October 2014:

422 1. *Appropriate cell lines to use in the in vitro micronucleus assay when testing*
423 *nanomaterials:*

424 29. Suspension cells are recommended (e.g. TK6 cells, buffy coat cells). Within the Phase 2 work, it
425 was evident that these cells were more sensitive as compared to adherent V79 and HepG2 cells.

426 30. It is important to note that the most sensitive cell system, preferably of human origin (regardless
427 the expected route of exposure to the nanomaterial in the final consumer product) and p53 competent
428 should be used. Intracellular uptake of the NMs should be demonstrated, and only cells that can efficiently
429 internalise NMs should be used.

430 31. As this Phase 2 study only included a small number of nanomaterials, it is recommended that
431 expanding the sensitivity comparison of TK6 versus V79 cells is continued through further inter-laboratory
432 trials using a wider range of nanomaterials.

433 2. *Positive controls:*

434 32. Currently, there are no suitable particle positive controls that have been identified for the *in vitro*
435 micronucleus assay. Thus, in the absence of particle positive controls, it is recommended that chemical
436 positive controls continue to be used (e.g. MMC, EMS).

437 33. Previously, there has been a suggestion in the scientific literature that tungsten carbide
438 nanoparticles could be considered as a particle positive control for the *in vitro* micronucleus assay with
439 human lymphocytes¹. Thus, within the Phase 2 study, the nanoparticle tungsten carbide was included for
440 further evaluation. This study demonstrated that tungsten carbide nanoparticles have the potential to be a
441 suitable particle positive control, but only when used with TK6 or buffy coat cells (however, the *in vivo*
442 relevance and the biological significance, i.e. germ cell mutagenesis or somatic cell carcinogenicity have
443 not yet been demonstrated).

444 3. *Top dose:*

445 34. It is recommended that the top dose should be restricted to 100 µg/mL or 100 µg/cm². Doses
446 higher than this are not physiologically relevant, but more importantly, can result in interference with scoring
447 due to high deposition on cells. Lower concentrations can be considered if they are justified by *in vivo*
448 organ burdens. Furthermore, agglomeration dynamics need to be carefully considered. Typically, at high
449 doses, agglomeration is more predominant with lower dispersion stability than at lower doses when a
450 stable dispersion of the particles can be more readily achieved.

451

452

¹ Moche H, Chevalier D, Barois N, Lorge E, Claude N, Nesslany F. Tungsten carbide-cobalt as a nanoparticulate reference positive control in in vitro genotoxicity assays. *Toxicol Sci.* 2014 Jan;137(1):125-34.

453 4. *Treatment time:*

454 35. One cell cycle treatment time, followed by 1.5 cell cycle exposure to cytochalasin B is
455 recommended for the cytokinesis blocked micronucleus assay.

456 36. It should be noted that this experimental time can be challenging to conduct within normal working
457 hours for the TK6 cells, which have a cell cycle time of 12-14h. Swansea University have therefore
458 repeated the Phase 2 *in vitro* micronucleus assay experiments using 1.5 cell cycle exposures. The data
459 demonstrated that there was no significant differences in the data generated when using either 1 or 1.5
460 cell cycle exposure with the TK6 cells. Thus, to enable the assay to be conducted during standard working
461 hours, for TK6 cells, a 1.5 cell cycle exposure is recommended and is detailed within the SOP in Annex 1.

462 5. *Dispersion:*

463 37. Several methods can be utilised to disperse nanomaterials, however within the Phase 2 trial
464 conducted, the NANoREG protocol (Alstrup Jensen et al., 2014) was implemented as this is a well-
465 established method.

466 6. *BASF and Swansea University have generated a harmonised SOP for the in vitro*
467 *micronucleus assay for testing of manufactured nanomaterials (Annex 1).*

468

Limitations

469

470 38. The described project has addressed several key issues regarding the testing of nanomaterials in
471 the *in vitro* micronucleus assay, nevertheless there are some issues which have either not been addressed
472 in this study or not resolved.

473 39. The recommendations in this guidance document are derived from data obtained from a small
474 number of compounds performed in two laboratories using different cellular test systems. A cross validation
475 of the protocols was only partly performed by Swansea University using V79 cells. A larger scale inter-
476 laboratory validation comprising of several independent laboratories would definitely add value to the
477 obtained results. It should be noted that if future ring-trials are to be conducted, it is critical to ensure that
478 all laboratories use an identical approach to dispersing their materials; if harmonisation in dispersion
479 approach across the participating laboratories cannot be achieved then the data may not align. For
480 example, if sonication is being applied to disperse nanomaterials, then delivery of specific energy levels is
481 required.

482 40. The number and type of nanomaterials used for this project were more or less predefined by the
483 previous work carryout out by JRC². A follow up trial should focus on a broader range of nanomaterials
484 representing a larger portion of the existing nanomaterials. In this context, the assessment of the *in vivo*
485 relevance of the obtained results with the tested nanomaterials is very limited. In this study CeO₂ was the
486 only compound with relevant *in vivo* follow up data (Ma-Hock et al³). Thus, a selection of further
487 nanomaterials with solid *in vivo* data would provide a better opportunity for interpretation of *in vitro* data.

488 41. During the meeting where the structure of the protocol was decided upon by the expert team the
489 subject of the highest concentration to be tested was thoroughly discussed. Thus, 100 µg/mL was selected
490 as the optimal highest concentration to be used in the experiments. Despite the unanimous decision by
491 the expert team, it is recommended to discuss this topic in a broader community to establish a range finding
492 workflow for the selection of the top concentration to be used for mutagenicity assessment.

493 42. Inclusion of S9 mix in the study protocol was not considered as relevant, since all selected
494 nanomaterials were inorganic and not likely subject to metabolization for induction of mutagenicity.
495 However, ECHA recommends the use of S9 mix for organic nanomaterials. The issues not addressed in
496 this document relate to firstly identify the plausibility of use of S9 for organic nanomaterials, the protocol
497 on how to incorporate an S9 Treatment within a test substance treatment period corresponding to one cell
498 cycle as well as the identification of a corresponding positive nanomaterial requiring metabolization. Future

² Drewes, C., Ojea Jimenez, I., Mèhn, D., Colpo, P., Gioria, S., Bogni, A., Ponti, J., Kinsner-Ovaskainen, A., Gilliland, D. and Riego Sintes, J., Physicochemical characterisation of gold, silica and silver nanoparticles in water and in serum-containing cell culture media, EUR 29054 EN, Publications Office of the European Union, Luxembourg, 2018, ISBN 978-92-79-77705-9 (online),978-92-79-77704-2 (print), doi:10.2760/818663 (online),10.2760/58721 (print), JRC110379.

³ Ma-Hock, L., Gröters, S., Strauss, V., Keller, J., Wiench, K., van Ravenzwaay B. and Landsiedel, R., Long-term inhalation study with CeO₂ and BaSO₄ nanomaterials. Manuscript in preparation.

499 metabolism considerations with respect to the replacement of animal derived S9 mix could include in vitro
500 metabolism assays such as TGP project 4.150.

501

Additional Supporting Interlaboratory Trial

502

503 43. Since completion of the Phase 2 experimental work, the Horizon 2020 project RiskGONE has
504 agreed to undertake an interlaboratory trial using the harmonised SOP generated through this Phase 2
505 project. This interlaboratory trial will include three laboratories: Swansea University (UK), Norwegian
506 Institute for Air Research (NILU; Norway) and the French Agency for Food, Environmental and
507 Occupational Health & Safety (ANSES; France). The interlaboratory trial will involve each partner exposing
508 the TK6 cell line to the engineered nanomaterials:

- 509 • Zinc Oxide (ZnO) – Sigma #MKCJ4155
- 510 • Titanium dioxide (TiO₂) – JRC #JRCNM01005a990582
- 511 • Tungsten carbide-cobalt (WC/Co) – NanoAmor #5561HW
- 512 • Chemical positive control: Mitomycin C (MMC) – Sigma, #M4287-2MG. Suspend in double
513 distilled water at 1mg/ml, keep at 4°C.

514 44. The ring trial will consist of each partner testing the TK6 cell line with the selected engineered
515 nanomaterials (ENMs), following the SOP in **Annex 1**. Experiments are to be performed in biological
516 duplicate for an $n=2$ (Two independent experiments). The ENM exposure period will be for 1.5 cell cycles
517 of TK6 cells and a cyto B time of 1.5 cell cycles.

518 45. This trial started in September 2021 and the data collation is planned for completion by January
519 2022. *[note from the Secretariat: analysis of the data and refinement to the considerations provided in this
520 document will be considered in a future version of this document].*

521

522

523

524

525

526

527

Annex 1 – Harmonised SOP generated by BASF & SU

528

529 46. This annex includes the finalised SOPs for testing engineered nanomaterials with the *in vitro*
530 micronucleus assay using TK6 and buffy coat cells, as these were the most sensitive cell lines identified.

531 Experiments for characterization of nanoparticles, recommended by both BASF 532 & Swansea University

533 47. The interpretation of *in vitro* genotoxicity requires a minimum physical-chemical characterisation
534 of the test item, including both intrinsic properties (size, shape, specific surface area) and extrinsic
535 properties (agglomeration and solubility in the genotoxicity test medium). The following (non-GLP)
536 technologies can be applied to generate this information:

537 1. Transmission Electron Microscopy (TEM) to determine the particle size distribution and shape
538 distribution. These are ECHA requirements to identify the Nanoform (or non-nano-form) of the
539 test item. (European_Chemicals_Agency_(ECHA), 2019) The method is selected accordingly
540 to the NanoDefine Method Manual, published by JRC in 2020 (Mech et al., 2020). Specifically,
541 a tip of a spatula of material is dispersed in 5ml of ethanol. The dispersion is treated for 5min
542 in an ultrasonic bath. 5 drops of the dispersion are applied between two glass slides to create
543 a thin liquid film. A carbon coated TEM grid is dipped on the film to transfer the particles to the
544 TEM substrate. The TEM samples are analysed using a Tecnai Osiris machine (Thermo-
545 Fisher) operated at 200 keV under bright-field (BF) conditions (Spotsize; 100µm condenser
546 aperture; 20µm objective aperture). About 10 Images per sample are acquired using a Gatan
547 1000XP 2K CCD camera with an acquisition time of 2s at magnifications of 13500x and
548 26500x. The particle size distribution and aspect ratio distribution are evaluated using the free
549 ParticleSizer plugin for ImageJ developed and validated within the EC FP7 project
550 NanoDefine. Automated evaluation is faster, reproducible and limits operator bias to a
551 minimum compared to manual image analysis. The report includes distributions, median D50
552 and statistics on particle size and aspect ratio.

553 2. Brunauer-Emmett-Teller (BET) method to determine the specific surface area. This is an
554 ECHA requirement to identify the Nanoform (or non-nano-form) of the test item
555 (European_Chemicals_Agency_(ECHA), 2019). The method is selected accordingly to the
556 NanoDefine Method Manual, published by JRC in 2020 (Mech et al., 2020). Specifically, the
557 method conforms with ISO_9277:2014 “Determination of the specific surface area of solids by
558 gas adsorption – BET method. If the test item is organic, then the degassing temperature is
559 limited to 80°C, or at least 50°C below melt temperature.

560 3. Analytical Ultracentrifugation (AUC) method to determine the agglomeration in genotoxicity
561 test medium. Characterise two concentrations (0.03 mg/mL and 0.1 mg/mL) and repeat one

562 concentration at 20h (conforming with the duration of the genotoxicity testing) to assess the
563 stability throughout the incubation time. AUC implements an optical system that is
564 synchronized with the centrifugal frequency to monitor the radial concentration profile during
565 centrifugal separation. The samples are prepared by the identical protocol as for the
566 genotoxicity testing in the identical medium. Suspensions are then directly loaded in AUC cells
567 without further preparation. All measurements are done on a Beckman Coulter XLI Proteome
568 Lab (Beckman Coulter Inc., Brea, USA). The cells are mounted on a 4-sample or 8-sample
569 titanium rotor, and data is acquired during a gravitational sweep from 73g to 116000g over 2h.
570 The samples are measured using either interference or absorbance optics to detect
571 sedimentation as described by the NanoDefine project (Mehn et al., 2018). Both types of
572 sedimentation data are evaluated with sedfit v15.0 software and are represented in particle
573 size distributions in mass metrics and in number metrics.

574 4. Incubation and filtration with analysis by Inductively-Coupled-Plasma Mass Spectrometry
575 (ICPMS) or by UV Vis-Spectrometry, to determine the static solubility in genotoxicity test
576 medium. Characterise two concentrations (0.03 mg/mL and 0.1 mg/mL). The samples are
577 prepared by the identical protocol as for the genotoxicity testing in the identical medium.
578 Suspensions are then incubated at 37°C for 20h, conforming with the conditions of the
579 genotoxicity testing. After 20 hours, the suspensions are filtered through a 1 µm glass filter
580 directly followed by a 0.02 µm (= 20 nm) alumina membrane. Thereafter, the filtrates are
581 analysed by UV-Vis (for organic pigments) or ICP-MS (for metal-based particles). If the test
582 item consists of Si or Al, polymeric filters are used instead.

583 a) Calibration of UV-Vis at peak pigment wavelength for quantification of the dissolved
584 material is performed by dissolving the pigment in concentrated sulfuric acid to get the
585 mass attenuation coefficient. The UV-Vis instrument Agilent Cary is used with 50mm
586 quartz cells to optimize detection down to 0.001 mg/L. Blank controls (no test item in the
587 dissolution setup) and medium controls (pure medium) are conducted. The general
588 procedure is following "UV-VIS Absorption Spectra (spectrophotometric method)", OECD
589 guideline for Testing of Chemicals, guideline 101, adopted 12th May 1981. For pigments
590 containing metals, the dissolved ions are quantified by ICP-MS (Perkin Elmer Nexion
591 2000b). The samples are diluted by a factor of 10 to 100 in 1% HNO₃, then measured with
592 50ms integration time on each mass. External calibration uses concentrations of
593 0/0.01/0.1/1/10 µg/L with matrix-matched standards. The ultra-high purity quartz
594 nebulizer has an argon flow of 0.94 mL/min and a sample flow of 0.29 mL/min at a pump
595 rate of 35rpm.

596 *In vitro* micronucleus test on
597 primary human lymphocytes;
598 Cytochalasin B method with
599 nanoparticles
600 BASF
601 Standard Operating Procedure
602 2020

603 Preparation of the cells

604 48. The blood used for testing comes from volunteer test persons, the administration of personal data
605 as well as the sampling procedure are described in the SOP "Organization of blood sampling for
606 genotoxicological tests". The selection of possible test persons is based on the criteria specified in OECD
607 guideline 487.

608 • Included are; Men and women aged 18 to 35, non-smokers, do not depend on regular
609 medication, and do not use drugs.

610 • Excluded are; pregnant women, people with known diseases or exposure to genotoxic
611 substances (chemicals / ionizing radiation).

612 • The cells will be incubated at 37°C, 5% (v/v) CO₂ and ≥ 90% relative humidity in the
613 incubator. The culture medium RPMI contains stable glutamine and is supplemented with
614 20% [v/v] FCS and 1% (v/v) Penicillin / Streptomycin. The culture medium is preheated to
615 37°C.

616 Stimulation of the cells with phytohemagglutinin

617 49. Erythrocytes are first isolated from the whole blood. For this purpose, the blood is first diluted 1:1
618 with culture medium. The erythrocytes are then isolated from the remaining nucleated cells (buffy coat) by
619 means of a Ficoll separation.

620 50. For the Ficoll separation, 25 mL Ficoll-Paque™ (GE Healthcare Bio-Sciences AB) are placed in
621 50 mL centrifuge tubes and carefully covered with 20 mL of the diluted blood. The centrifuge tubes are
622 centrifuged at 1055g for 20 minutes. The centrifuge's braking function should be switched off. The "buffy
623 coat" layer (the layer between the medium and Ficoll) is carefully collected and centrifuged again at 491g
624 (15 minutes with braking function). The pellet is then washed twice with medium and centrifuged (5 minutes
625 at 900g). The pellet is taken up in a volume of cell culture medium (RPMI; 20% FCS, 0.6 mg / mL PHA)
626 that corresponds to 5 times the blood volume (e.g. for 20 mL whole blood, the cells are taken up in 100
627 mL).

628 51. 50 mL of each batch are placed in a 175 cm² cell culture flask and cultivated in the incubator for
629 48 hours.

630 Test substance preparation and preparation of the particulate positive control

631 52. Tungsten carbide-Cobalt will be used as particulate positive control. The nanoparticle test
632 substance and particulate positive control will be prepared as described by the NANoREG protocol (Alstrup
633 Jensen et al., 2014).

634 53. In summary, a certain amount of the respective nanoparticle is weighed and mixed with ethanol
635 (pre-wetting) so that the final ethanol concentration is 0.5 vol%. Then the dispersion medium (0.05% w / v
636 BSA water) is added so that a particle concentration of 2.56 mg / mL is reached. The dispersion is treated
637 with a probe ultrasonicator (Branson Sonifier 550 W) for 16 minutes at 10% amplitude (approx. 7 W). This
638 stock dispersion is used to create the different testing concentrations.

639 **The test substance is tested with the following concentrations:**

- 640 • 1, 3, 10, 30, 60, 100 µg/mL

641 **Tungsten-Carbide/Cobalt (8wt% WC/Co <200nm, 99.5% LOT#5561-072018, Nanostructured &**
 642 **Amorphous Materials Inc., USA) particulate positive control is tested with the following**
 643 **concentrations:**

- 644 • 30, 60, 100 µg/mL

645 **Preparation of the chemical positive control**

646 54. The positive controls are prepared fresh on the day of the test shortly before the substance is
 647 administered or prepared in larger quantities at an earlier point in time, stored in ready-to-use portions at
 648 80°C, and thawed and used shortly before use.

649 55. Mitomycin C (0.04µg/ml) and colchicine (0.05µg/ml) are used as chemical positive controls.

650 56. Several concentrations of a positive control substance can be used, it being enough if only one
 651 concentration used is evaluated if it shows a sufficiently strong genotoxic effect in the experiment.

652 **Preparation and treatment of test cultures**

653 57. The substance preparations and the preparations of the particulate positive control are diluted to
 654 the final concentration in culture with culture medium, the total volume (for two cultures (duplicate
 655 determination) is 30mL. The pipetting scheme for dilution is shown in Table 1.

656 **Table 1.** Pipetting scheme of test substance preparations

Dose [µg/mL]	dilution in BSA			dilution in culture medium RPMI incl. 20% FCS	
	[mL]	from Dosis	+ vehicle [mL]	[mL] Ansatz	aus + vehicle [mL]
100*	5	2560	7.8	3	27
60*	6	100	4	3	27
30*	5	60	5	3	27
10	0.5	100	4.5	3	27
3	0.5	30	4.5	3	27
1	0.5	10	4.5	3	27

657 *: applies to both test substance and particulate positive control

658 58. The blood cultures are removed from the incubator, suspended and pooled. For each culture (two
 659 cultures per test group = double determination), 10 mL are pipetted into a prepared 15 mL centrifuge tube.
 660 The cells are centrifuged at 900g for 5 min and the supernatant is removed.

661 59. For each culture (double determination) 10 mL of the prepared substance dispersion are added to
662 the cells. The cells are suspended with a 10 mL serological pipette and transferred into prepared 25 cm²
663 cell culture flasks. The cells are incubated on a shaker (approx. 150 rpm) in the incubator for 20 hours (1-
664 cell cycle).

665 **Checking or determining further parameters**

666 60. Changes in the pH value can be recorded by the colour change of the indicator dye of the culture
667 medium (phenol red: no colour change from pH 6.7 - 8.3).

668 61. The pH value is measured at the start of treatment at least for the highest dose and for the
669 untreated control or the vehicle control (pH 6.7-8.3).

670

671 **Test substance removal**

- 672 • The exposure phase is ended by washing the cells several times with HBSS
- 673 • Transfer the cultures into prepared 15 mL centrifuge tubes
- 674 • Centrifugation (5 min at 900g)
- 675 • Remove the supernatant
- 676 • Suspend the cells in 5 mL HBSS
- 677 • Centrifugation (5 min at 900g)
- 678 • Remove the supernatant
- 679 • Suspend the cells in 5 mL HBSS
- 680 • Centrifugation (5 min at 900g)
- 681 • Remove the supernatant
- 682 • Suspend the cells in 10 mL fresh medium
- 683 • Transfer of the cultures into prepared 25 cm² cell culture flasks

684 **Cytochalasin B treatment**

- 685 • Following a 1-cell cycle exposure to the test nanoparticles, each culture is treated with 30 µL initial
- 686 solution (final concentration 6 µg/mL culture medium). The cultures are incubated for an additional
- 687 20h (1-cell cycle) in the incubator.

688 **Cell harvest and preparation of slides**

- 689 • The preparation takes place under a fume hood under non-sterile conditions.
- 690 • The fixative (19 parts methanol: 1 part acetic acid) is freshly prepared every working day, during
- 691 the incubation steps the fixative is stored at -20°C.
- 692 • The hypotonic solution (0.65% KCl) is pre-warmed to 37°C in a water bath, it can be prepared the
- 693 day before and stored in the refrigerator until use.

694 **Hypotonic treatment and fixation of the cells:**

- 695 62. The cells will be transferred into 15mL tubes, centrifugated at 900g for 5 min (room temperature)
- 696 and washed with HBSS. Then the cells will be treated with a hypotonic solution and fixative:
- 697 • Centrifugation (5 min, 900g, at room temperature)
 - 698 • Suspend in 5 mL KCl, cover with Parafilm
 - 699 • Incubate for 10 min at 37°C
 - 700 • Fix the cells by adding 1 mL of cold fixative and suspending them in the solution
 - 701 • Centrifugation (5 min at 900g and 4°C)
 - 702 • Remove the supernatant

- 703 • Suspend in 5 mL cold fixative
 704 • Incubate for 20 min at 4°C
 705 • Repeat the fixation step with 5 mL fixative (incubation for 20 min at 4°C)
 706 • Centrifugation (5 min at 900g and 4°C)
 707 • Suspend in 3 mL fixative
 708 • After the last fixation step, the cells can be stored for up to six hours at 4°C or applied directly to
 709 the slide

710 Preparation of the cells on microscope slides

- 711 • Three preparations are made per culture as follows:
 712 • Centrifuge the cells (5 min at 900g and 4°C) and suspend in 1-2 mL cold fixative. Depending on
 713 the pellet size / cloudiness of the cell suspension, additional fixatives can be added.
 714 • Moisten the slide by immersing it in deionized water, allow it to drain briefly
 715 • Apply 250 µL of the cell suspension evenly to the microscope slides
 716 • Pull the slide through the flame of the Bunsen burner
 717 • Extinguish the flame after approx. 1 second
 718 • The cell density can already be assessed microscopically on the native preparation and, if
 719 necessary, adjusted (by varying the amount of fixative). The cells should be isolated and not
 720 overlapping on the slide.
 721 • Let the preparations dry at least overnight.

722 Staining of the slides

- 723 63. The cells were stained with May-Grünwald (3 min) and 10% [v/v] Giemsa (in Titrisol, pH 7.2, 10
 724 min) and mounted.

725 Cytotoxicity – proliferation index (CBPI)

- 726 64. Micronuclei can only arise in proliferating cells. The cytokinesis block proliferation index (CBPI) is
 727 a direct measure for determining the cell division activity. The number of mononuclear, binuclear and
 728 multinucleated cells is determined in 500 cells per culture (culture = slide), in two cultures per test group =
 729 1000 cells. The CBPI can be evaluated without coding.

- 730 65. If only isolated cells are found, no analysis is carried out. This is noted on the evaluation form (n.s.
 731 = not scorable).

- 732 66. The CBPI calculated from this distribution indicates the average number of cell cycles that a culture
 733 has undergone in the presence of the actin polymerization inhibitor Cyto B. The CBPI is given as an
 734 absolute value without a unit.

$$735 \text{ CBPI} = \frac{((\text{No. mononucleate cells}) + (2 \times \text{No. binucleate cells}) + (3 \times \text{No. multinucleate cells}))}{(\text{Total number of cells})}$$

- 736 67. The CBPI was used to calculate the % cytostasis (relative inhibition of cell growth compared to the
 737 respective vehicle control group) - a CBPI of 1 (all cells are mononucleate) is equivalent to 100% cytostasis.

$$738 \% \text{ Cytostasis} = 100 - 100 \{(\text{CBPIT} - 1) / (\text{CBPIC} - 1)\}$$

739 T = test substance treated culture C = vehicle control culture

740 **Micronucleus analysis**

741 68. The selection of the dose groups that are used to record the micronucleus frequency depends on
742 the quality of the preparations or cells, the cytotoxicity and possible other physiological parameters (e.g.
743 pH value). This dose selection is made considering the applicable international test regulations. The
744 evaluation of the cells for micronucleus frequency is carried out manually with a transmitted light
745 microscope (e.g. Axio Scope.A1, Zeiss) on coded preparations. A magnification of 400 times is used. The
746 findings are recorded manually on an evaluation form.

747 69. At least 1000 binuclear cells are analysed per culture (slide) and the number of micronucleated
748 cells is recorded. A total of 2 slides per test group, i.e. at least 2000 binuclear cells evaluated. The criteria
749 for evaluating a cell are:

- 750 • The plasma border must be complete (otherwise main cores or MN may be missing)
- 751 • The plasma must not be too dark in colour
- 752 • The main nucleus should not differ too much in size
- 753 • The cell must be clearly distinguishable from the neighbouring cell
- 754 • Apoptotic cells with a fragmented nucleus are not assessed

755 **The following criteria must be evaluated for the assessment as a micronucleus**

- 756 • The micronucleus must be round (no angular / crystalline structures)
- 757 • The diameter of the micronucleus was less than 1/3 of the main nucleus
- 758 • The micronucleus was not linked to the main nucleus and was located within the cytoplasm of the
759 cell
- 760 • Only binucleated cells were scored
- 761 • The micronucleus must lie completely in the plasma

762 **Appendix 1 for BASF SOP**763 **Materials and chemicals**

Materials	Manufacture	Order-no.
Cell culture flasks, 25 cm ²	TPP	P 90025
Cell culture flasks, 75 cm ²	Greiner	661160
Centrifuge tubes 15 mL	TPP	TPPA 91015
Centrifuge tubes 50 mL	TPP	TPPA 91050
microscope slide, free of grease	Superfrost	631-1320

764

Reagents und medium	Hersteller	Best.-Nr.
RPMI 1640, w: Stab. Glutamine, w: NaHCO ₃ , ohne Phenolrot	PAN biotech	P04-16520
Hanks' Balanced Salt Solution <u>ohne</u> Ca ²⁺ , Mg ²⁺ (HBSS)	PAN biotech	P04-34500
Fetal Calf Serum (FCS)	Biowest	-
Penicillin/Streptomycin (10 000 E/10 000 µg/mL)	PAN biotech	P06-07100
Mitomycin C (MMC)	Roche	10107409001
Colchicin	Roche	10295892001
KCl	Bernd Kraft	4016056
methanol	Bernd Kraft	34940
glacial acetic acid	Bernd Kraft	4146307
Mounting Pertex®	Medite	MEDT 41-4011-00

765 70. A product or supplier of comparable quality can be used for all articles at any time. Reagents not
766 listed here are described in more detail in the test plan if required. Other chemicals used are customary in
767 the laboratory and are obtained from common manufacturers or are available in the laboratory.

768

769

770 **Appendix 2 for BASF SOP**771 **Recipe**772 **Culture medium**

773 71. RPMI 1640, w: Stab. Glutamine, w: NaHCO₃, ohne Phenolrot unter Zusatz von:

- 774 • fetal calf serum (FCS), ca. 20% (v/v)
- 775 • Penicillin/Streptomycin (10 000 E/10 000 µg/mL), ca. 1% (v/v)

776 72. All components are mixed under sterile conditions in the medium bottle and marked with the
777 ingredients, date and stored for a maximum of 6 weeks at approx. +4°C.

778 **Phytohämagglutinin**

779 73. The PHA (1.2 mg) is dissolved in 2 mL of ultrapure water in a container and can be stored in the
780 refrigerator for up to 6 weeks. The documentation is carried out in the raw data across studies.

781 **Cytochalasin B – solution**

782 74. The cytochalasin B (10 mg) is dissolved in the container in 5 mL DMSO (final concentration 2
783 mg/mL) and can be stored in the refrigerator for up to 6 weeks.

784 **Hypotonic treatment**

785 75. The hypotonic solution (0.65% KCl) can be prepared on the day before and stored in the
786 refrigerator.

787 76. Approach: Weigh of 0.65g KCl/ 100 mL with ultrapure water

788 **Fixative**

789 77. The fixative is freshly applied on a working day and stored at -20°C during the fixation steps.

790 78. Approach: 19 parts methanol (cooled overnight at -20°C) + 1-part glacial acetic acid

791 *In vitro* cytokinesis-blocked
792 micronucleus (CBMN) test
793 protocol for use with TK6 cells
794 and engineered nanomaterials
795 (ENMs)
796 Swansea University
797 Standard Operating Procedure
798 2020

799

Abbreviations

- Cyto B – Cytochalasin B
- EDTA – Ethylenediaminetetraacetic Acid
- ENM – Engineered Nanomaterial
- FBS – Foetal Bovine Serum
- HS – Horse Serum
- KCl – Potassium Chloride
- MMC – Mitomycin C
- MMS – Methyl Methanesulphate
- MNU - N-nitro-N-methylurea

- NaCl – Sodium Chloride
- P/S – Penicillin/Streptomycin
- PBS – Phosphate Buffered Saline

Biological Setting & Test System

79. This SOP should be carried out under strict laboratory conditions, with all work performed under sterile conditions and in a Class II Laminar Tissue Culture Hood.

Chemicals & Reagents

Table 1. List of chemicals and reagents required for *in vitro* CBMN assay.

Reagent	Supplier
PBS	ThermoFisher; 10010023
Cell culture medium	ThermoFisher; 21870076
Cytochalasin B	Sigma; C6762
Methanol	Sigma; 34860
Beckman Coulter Diluent	Beckman; 628017
Giemsa Stain	VWR; 350864
Gurr pH 6.8 Buffer tablets	VWR; 331542
Glutamine	ThermoFisher; 25030081
Horse Serum	ThermoFisher; 16050122
Trypsin-EDTA (0.25%)	ThermoFisher; 25200056
Potassium Chloride (KCl)	Sigma; P9333
Sodium Chloride (NaCl)	Sigma; S7653
Acetic acid	Sigma; 1005706
VECTASHEILD Mounting Medium with DAPI	VECTOR Laboratories; H-1200
DPX Mountant	VWR; 13510
Xylene	VWR; 214736

Cell culture preparation

80. In preparation for this SOP, it is advised that cell culture medium is prepared and pre-warmed at 37°C for 30 minutes prior to use. TK6 complete cell culture medium is prepared by adding 50ml of horse serum to 500ml of media (RPMI 1640), 5ml of L-glutamine is also added to complete the media. The full cell culture medium should then be mixed prior to use by inverting the bottle.

Procedure

81. To thaw cells from Liquid Nitrogen:

- Place a bottle of complete medium in the water bath at 37°C (25-30 minutes before use). Remove your vial of cells from the Liquid N₂ and place in a foam float in the water bath. Do not allow the

vial to be immersed. Cells will thaw in ~1 minute, transfer them to a flask of the warm medium slowly and carefully using a Pasteur pipette. (The amount of medium depends on cell type and on rate of cell division; 40-50ml for suspension cell lines and 10-15ml for adherent).

- Label the flask with the name of the cells, the passage number, your name, the date of resuscitation and the date that the cells were previously frozen. (This is helpful because if there is a problem with them, the other vials from the batch can be identified.)
- Place the flask in the CO₂ incubator.
- Check them after 24 hours for growth and contamination. They can be counted to determine when they will need splitting. Cells should not be allowed to become too confluent (split once 70-80% confluent). If required change media every 2 days.
- Monitor passage number in your experiments; do not let this rise beyond a reasonable range. TK6 cells should not be grown continuously for more than 4 weeks.
- At day 4 the SOP can be performed 2 different ways, such that scoring can be conducted using a manual approach, or alternatively, a semi-automated scoring approach using the Metafer microscope.

DAY 1 – seed cells

- Count cells using a Beckmann Coulter Counter and seed at 1×10^5 cells/ml in 10ml of media in 25cm³ flasks per treatment. Two sets of flasks are required per treatment dose: one for genotoxicity (micronucleus frequency analysis) and a second additional satellite flask per treatment for cytotoxicity calculations.
- Incubate overnight at 37°C/ 5% CO₂

DAY 2 – cytotoxicity count 1 & dosing of cells with test agent

Concurrent positive and solvent/vehicle controls should be included in each replicate. For TK6 cells, the recommended chemical positive control is:

MMC CAS 50-07-7 (Fisher Scientific; BP2531-2), dose for 1 – 1.5 cell cycles at 0.01µg/ml

This procedure must be conducted in a Class II Laminar Tissue Culture Toxic Hood, with the user wearing double gloves to ensure safety when dosing with chemicals and / or ENMs.

Preparation of ENMs

ENMs will require preparation, suspension and sonication prior to use in toxicology

testing, allow time to prepare ENMs fresh on the day of dosing. User's should refer to a recognised dispersion procedure which generates a stable ENM suspension and where the final quantity of ENM suspension added to the cells does not exceed 1:100 as to avoid disturbing cell culture conditions (OECD TG487). For this protocol we use the NANoREG guidance document for ENM suspension and sonication (Alstrup Jensen et al., 2014).

- Count satellite flasks for cytotoxicity (initial cell number) 1h before dosing. 100µl of cell suspension is added to a cuvette containing 10ml of diluent for cell concentration determination using a Beckmann Coulter Counter.
- Dose cells in both sets of flasks (for genotoxicity and cytotoxicity evaluation) for 1 – 1.5 cell cycles (for TK6 cells, 1 cell cycle time is typically 12-15 hours) with the ENM at the desired final concentration; this is to be performed in triplicate per concentration, including triplicate negative control (media only, or vehicle used to suspend ENM) and positive controls.

DAY 3 – cytotoxicity count 2 & addition of cytochalasin B

- Count satellite flasks (post-treatment cell number) to calculate cytotoxicity. Add 100µl of cell suspension directly from satellite flasks to a cuvette containing 10ml of diluent for cell concentration determination using a Beckmann Coulter Counter.
- Transfer the content of the genotoxicity flasks to a 15 ml centrifuge tube, centrifuge at 230g for 5 min and discard the supernatant. Re-suspend the pellet in pre-warmed PBS and re-centrifuge at 230g for 5 min, repeat this wash step a second time and discard the supernatant.
- Resuspend cell pellets in fresh media containing 3-6 µg/ml cytochalasin B and place into new T25 flasks.
- Incubate for 1.5 cell cycles (18 hours)

DAY 4 – harvesting cells for automated scoring protocol using the Metafer system

- Harvest cell pellets by centrifugation (230g for 5 min), resuspend in 5 ml pre-warmed PBS and centrifuge cells at 230g for 10 min. Discard the supernatant and repeat this wash step a second time.
- Re-suspend cell pellets with hypotonic solution (KCl 0.56%), then centrifuge immediately at 230g for 10min
- Re-suspend pellets in Fixative 1 (methanol: acetic acid: 0.9% NaCl (5:1:6 parts)) and incubate at 4°C for 10min before centrifugation at 230g for 10min.
- Re-suspend pellets in Fixative 2 (methanol: acetic acid (5:1 parts)) and incubate at 4°C for 10 min before centrifugation at 230g for 10 min. Repeat this wash step a further 3 times. Maintain cells in the last fix wash overnight at 4°C.
- Place freshly opened microscope slides in a glass tank of Fixative 2 at 4°C, at least a 2 hours before slide preparation (ideally overnight). On the day of preparing slides, replace the fix with ddH₂O.
- Centrifuge (230g for 10 minutes) the fixed cell suspensions and thoroughly re-suspend in ~1ml Fixative 2
- Take a slide out of the ddH₂O and wipe the water off the upper side with slide tissue, ideally with one movement (the surface should be dry or with only a faint film of water remaining). Pipette 100µl of the cell suspension evenly onto the slide
- Wait a few seconds until the suspension is evenly spread over the slide, and then put it in a vertical position for drying
- Check cell density of binuclear cells and if required adjust the final re-suspension volume by either lowering or increasing the volume of Fixative 2 added. Cells should not be overlapping, densely packed or too sparse.
- Stain slides with 30µl (3 dots of 10µl) of Vectashield mounting medium with DAPI, apply coverslip and incubate in the dark for 15 min
- Score slides on a TK6 cell classifier using the automated Metafer microscope (Axio-imager Z2 fluorescent microscope, Carl Zeiss UK) Metafer 4 software version 3.5 (MetaSystems, Germany). Not all cell lines will be suited to this classifier – it is specific to TK6 cells. Score 2000 BN cells per ENM

concentration/per replicate, (6000 BN cells in total per ENM concentration). Classifier information can be found in the work by Seager and colleagues (Seager et al., 2014). We have also provided nuclei and micronucleus classifier settings in the Appendix.

DAY 4 – harvesting cells for manual scoring

- Harvest cell pellets by centrifugation (230g for 5 min), resuspend in 5 ml pre-warmed PBS and centrifuge cells at 230g for 5 min. Discard the supernatant and repeat this wash step a second time.
- Resuspend cells in 10ml of PBS (if the cell pellet looks small adjust the volume as necessary)
- Place labelled slides (90% methanol-cleaned) into cyto-clips, place a filter card (with 5 mm cyto-dot hole) on top of the slide and clip in place a cyto-funnel. Load 100µl of cell suspension into each funnel and centrifuge at 200g for 10 min in a Shandon Cytospin.
- Examine slides for correct cell density and adjust the volume of cell suspension as required. If the cells are too sparse, centrifuge the cell suspension at 230g for 5 min and resuspend in a smaller volume of PBS.
- Fix slides for 10min in ice cold 90% methanol and leave to air dry (at this point slides can be stored at -20°C)
- Stain slides in 20% Giemsa solution (5ml of Giemsa + 20ml pH 6.8 Gurr buffer, filtered) for 10 minutes
- Rinse slides in pH 6.8 Gurr buffer, then soak in pH 6.8 phosphatase buffer for 1-2min. Note: two slide tanks can be set up simultaneously to improve efficiency
- Leave slides to air dry standing vertically
- Dip slides in xylene for 10 seconds using tweezers and drain off the excess
- Drop DPX over the area of cells
- Place 22x22 mm coverslip over the DPX, ensuring there are no air bubbles by pressing lightly on the coverslip
- Leave slides to dry for 24h in the fume hood

- View slides under a light microscope, evaluate 1000 BN cells per replicate per concentration of ENM for the presence of micronuclei (3000 BN cells in total).

Cytotoxicity Calculations

As described in OECD TG487, CBPI, RPD or RICC are the most suitable for quantifying cytostasis and cytotoxicity (paused growth and cell death respectively).

Relative Population Doubling (RPD):

$$\text{RPD} = \frac{\text{No. of population doublings in treated cultures}}{\text{No. of population doublings in control cultures}} \times 100$$

Where population doubling = $[\log (\text{Post-treatment cell number} / \text{Initial cell number})] / \log 2$

Relative Increase in Cell Counts (RICC):

$$\text{RICC} = \frac{\text{Increase in number of cells in treated cultures (final – starting)}}{\text{Increase in number of cells in control cultures (final – starting)}} \times 100$$

Cytokinesis Block Proliferation Index (CBPI):

% Cytotoxicity = $100 - 100 (\text{CBPI}_T - 1 / \text{CBPI}_C - 1)$; where T is treatment, C is control

CBPI = $(\text{No mononucleated cells} + 2 \times \text{BN cells} + 3 \times \text{MNCs}) / N$; N is the total number of cells scored

Appendices for Swansea University SOP

TK6 classifier settings on the Metafer microscope (Axio-imager Z2 fluorescent microscope, Carl Zeiss UK) Metafer 4 software version 3.5 (MetaSystems, Germany)

Nuclei

- Image Processing Operations: Sharpen (3,2) Median V (3) Median H (3)
- Object Threshold: 30%
- Minimum Area: $20.0\mu\text{m}^2$
- Maximum Area: $400.0\mu\text{m}^2$

- Maximum Relative Concavity Depth: 0.900
- Maximum Aspect Ratio: 1.500
- Maximum Distance: 30.0 μm^2
- Maximum Area Asymmetry: 70%
- Region of Interest Radius: 40.0 μm^2
- Maximum Object Area in ROI: 50.0 μm^2

Micronuclei

- Image Processing Operations: Median V (3) Median H (3) Sharpen (5,3) SB Histomax
- Object Threshold: 8%
- Minimum Area: 1.50 μm^2
- Maximum Area: 55.0 μm^2
- Maximum Relative Concavity Depth: 0.900
- Maximum Aspect Ratio: 4.000
- Maximum Distance: 25.0 μm^2

References

ALSTRUP JENSEN, K., KEMBOUCHE, Y., LOESCHUNER, K. & CORREIA, M. 2014. SOP for probe-sonicator calibration of delivered acoustic power and de-agglomeration efficiency for in vitro and in vivo toxicological testing, version 1.0. NANoREG.

EUROPEAN_CHEMICALS_AGENCY_(ECHA) 2019. Appendix for nanoforms applicable to the Guidance on Registration and Substance Identification. In: ECHA (ed.) ECHA-19-H-14-EN.

MECH, A., RAUSCHER, H., RASMUSSEN, K., BABICK, F., HODOROABA, V.-D., GHANEM, A., WOHLLEBEN, W., MARVIN, H., BRÜNGEL, R. & FRIEDRICH, C. M. 2020. The NanoDefine Methods Manual - Part 2: Evaluation of methods. The NanoDefine Methods Manual. Luxembourg: Publications Office of the European Union.

MEHN, D., RIO-ECHEVARRIA, I. M., GILLILAND, D., KAISER, M., VILSMEIER, K., SCHUCK, P. & WOHLLEBEN, W. 2018. Identification of nanomaterials: A validation report of two laboratories using analytical ultracentrifugation with fixed and ramped speed options. *NanoImpact*, 10, 87-96.

SEAGER, A. L., SHAH, U. K., BRÜSEHAFFER, K., WILLS, J., MANSHIAN, B., CHAPMAN, K. E., THOMAS, A. D., SCOTT, A. D., DOHERTY, A. T., DOAK, S. H., JOHNSON, G. E. & JENKINS, G. J. 2014. Recommendations, evaluation and validation of a semi-automated, fluorescent-based scoring protocol for micronucleus testing in human cells. *Mutagenesis*, 29, 155-64.