

1
2 **Draft GUIDANCE DOCUMENT ON GOOD *IN VITRO* METHOD PRACTICES (GIVIMP)**
3 **FOR THE DEVELOPMENT AND IMPLEMENTATION OF *IN VITRO* METHODS FOR**
4 **REGULATORY USE IN HUMAN SAFETY ASSESSMENT**

5
6 **FOREWORD**
7

8 A guidance document on Good *In Vitro* Method Practices (GIVIMP) for the development and
9 implementation of *in vitro* methods for regulatory use in human safety assessment was
10 identified as a high priority requirement. The aim is to reduce the uncertainties in cell and
11 tissue-based *in vitro* method derived predictions by applying all necessary good scientific,
12 technical and quality practices from *in vitro* method development to *in vitro* method
13 implementation for regulatory use.

14 The draft guidance is coordinated by the European validation body EURL ECVAM and has
15 been accepted on the work plan of the OECD test guideline programme since April 2015 as a
16 joint activity between the Working Group on Good Laboratory Practice (GLP) and the
17 Working Group of the National Coordinators of the Test Guidelines Programme (WNT).

18 The draft document prepared by the principal co-authors has been sent in September 2016 to
19 all 37 members of the European Union Network of Laboratories for the Validation of
20 Alternative Methods (EU-NETVAL¹) and has been subsequently discussed at the EU-
21 NETVAL meeting on the 10th of October 2016.

22 By November/December 2016 the comments of the OECD Working Group on GLP and
23 nominated experts of the OECD WNT will be forwarded to EURL ECVAM who will
24 incorporate these and prepare an updated version. A second round of commenting shall be
25 concluded in the beginning of 2017. EURL ECVAM shall then prepare the final GIVIMP
26 version, which will be submitted to OECD for proposed adoption at the OECD Joint Meeting
27 in April 2017.

28

¹ <https://eurl-ecvam.jrc.ec.europa.eu/eu-netval>

Field Code Changed

68	4.6	Dedicated media to particular cell lines	48	Field Code Changed
69	5	Test Systems	49	Field Code Changed
70	5.1	GCCP.....	49	Field Code Changed
71	5.2	Cell and tissue sourcing	49	Field Code Changed
72	5.3	Handling and maintenance of the test system.....	50	Field Code Changed
73	5.4	Cryopreservation	52	Field Code Changed
74	5.5	Cell line identity and genetic aberrations	52	Field Code Changed
75	5.6	Contaminants screening: sterility, mycoplasma, virus.....	53	Field Code Changed
76	5.7	Quality Control	54	Field Code Changed
77	5.8	Biomarkers and functional tests to confirm the required cell function state	55	Field Code Changed
78	5.9	Special issues for microbial strains	55	Field Code Changed
79	5.10	Qualification of reference strains.....	55	Field Code Changed
80	6	Test and reference items	57	Field Code Changed
81	6.1	Test item	57	Field Code Changed
82	6.1.1	Considerations during the development of the method	57	Field Code Changed
83	6.1.2	Considerations for the final user of the validated method	61	Field Code Changed
84	6.2	Interaction between test item and test system	61	Field Code Changed
85	6.2.1	Interference with the test system	62	Field Code Changed
86	6.2.2	Interference with <i>in vitro</i> method endpoint.....	67	Field Code Changed
87	6.2.3	Interference with the analytical endpoint	67	Field Code Changed
88	6.2.4	Consideration of interferences not coming from the active ingredient	68	Field Code Changed
89	6.3	Biokinetics / dose extrapolation	68	Field Code Changed
90	6.3.1	Kinetic processes	69	Field Code Changed
91	6.3.2	Measurement of free concentration /passive dosing	71	Field Code Changed
92	6.4	Reference and control items	73	Field Code Changed
93	6.5	Use of proficiency chemicals.....	74	Field Code Changed
94	7	Standard Operating Procedures (SOPs)	76	Field Code Changed
95	7.1	<i>In vitro</i> method SOP(s) development.....	76	Field Code Changed
96	7.2	Preparing Standard Operating Procedures	77	Field Code Changed
97	8	Performance of the method.....	80	Field Code Changed
98	8.1	Plate layout	80	Field Code Changed
99	8.2	Data analysis.....	82	Field Code Changed
100	8.3	Data-intensive <i>in vitro</i> methods	82	Field Code Changed
101	8.4	Acceptance criteria.....	83	Field Code Changed
102	8.5	Dynamic range/range of application	86	Field Code Changed
103	8.6	Signal intensity	87	Field Code Changed
104	8.7	Signal variability and plate uniformity assessment.....	88	Field Code Changed
105	8.8	Reliability of endpoint calculations.....	88	Field Code Changed

106	8.8.1	Confidence intervals for EC ₅₀ s.....	88	Field Code Changed
107	8.8.2	Minimum significant ratio.....	89	Field Code Changed
108	8.8.3	Variability and outliers.....	89	Field Code Changed
109	8.8.4	Transferability: in between lab transfer.....	89	Field Code Changed
110	8.9	Accuracy, reliability and uncertainty.....	89	Field Code Changed
111	9	Reporting of Results.....	92	Field Code Changed
112	9.1	Publishing.....	93	Field Code Changed
113	9.2	Reporting of method validation.....	93	Field Code Changed
114	9.2.1	Reporting Templates.....	94	Field Code Changed
115	9.3	Data reporting for regulatory purposes.....	95	Field Code Changed
116	10	Storage and retention of records and materials.....	97	Field Code Changed
117	10.1	Archiving of materials.....	97	Field Code Changed
118	10.2	Document and record management.....	98	Field Code Changed
119	10.2.1	Documents and records to be retained.....	99	Field Code Changed
120	10.3	Archiving and retention.....	99	Field Code Changed
121	10.3.1	Retrieval.....	100	Field Code Changed
122	10.3.2	Backup and restore.....	100	Field Code Changed
123	11	References.....	101	Field Code Changed
124		List of Figures.....	108	Field Code Changed
125		List of Tables.....	109	Field Code Changed
126		Annex 1 Good Cell Culture Practice (GCCP).....	110	Field Code Changed
127		Annex 2 Good Cell Culture Practice for stem cells and stem-cell-derived models.....	111	Field Code Changed
128		Annex 3 Experts participating at the GIVIMP meeting 24-25 February 2015.....	112	Field Code Changed
129				

130

Glossary of important terms used in the Guidance Document²

131	Acceptance criteria	Criteria for when study results can be accepted, i.e. a set of well-defined parameters describing aspects of the <i>in vitro</i> method such as control and reference item output, acceptable range for positive and negative controls, etc. These should primarily be established based on information from existing data on the finalised <i>in vitro</i> method or described in relevant bibliographic data. However, relevant development data obtained on a version of the <i>in vitro</i> method equal to the one used for generating the study results, as well as historical data, may need to be taken into account where available.
132		
133		
134		
135		
136		
137		
138		
139		
140	Adverse outcome pathway (AOP)	An AOP represents the existing knowledge concerning the causal linkages between the Molecular Initiating Events (MIEs) and the cascade of intermediate or key events at subcellular, cellular, tissue, organ, whole animal, and population levels that lead to a specific adverse outcome.
141		
142		
143		
144		
145	Apoptosis	Process of programmed cell death generally characterised by distinct morphological characteristics and energy-dependent biochemical mechanisms. Apoptosis is considered a vital component of various processes including normal cell turnover, proper development and functioning of the immune system, hormone-dependent atrophy, embryonic development and chemical-induced cell death.
146		
147		
148		
149		
150		
151		
152	Archive	A designated area or facility (e.g. cabinet, room, building or computerised system) for the secure storage and long term, permanent retention of raw data, completed data and relevant metadata in its final form and records for the purposes of reconstruction of the activity or study. It also enables audits.
153		
154		
155		
156		
157	Assay	A defined laboratory procedure that produces results. Also known as testing method. An assay can be considered as a technical operation that consists of determination of one or more characteristics of a given product, process or service according to a specified procedure. Often an assay is part of an experiment. The testing result can be qualitative (yes/no), categorical, or quantitative (a measured value). It can be a personal observation or the output of a precision measuring instrument. Usually the testing result is the dependent variable, the measured response based on the particular conditions of the testing method or the level of the independent variable. Some testing methods, however, involve changing the independent variable to determine the level at which a certain response occurs: in this case, the testing result is the independent variable.
158		
159		
160		
161		
162		
163		
164		
165		
166		
167		
168		
169	Batch	A specific quantity or lot of a test item or reference item, test system, assay reagent or other consumable, produced during a defined cycle of manufacture in such a way that it could be expected to be of a uniform character and should be designated as such.
170		
171		
172		
173		
174		
175	Benchmark dose (BMD) or concentration (BMC)	A point of reference by which something can be measured. The benchmark dose/concentration approach was developed as an alternative to the use of No Observed Adverse Effect Level (NOAEL) and Lowest Observed Adverse Effect Level (LOAEL).
176		
177		
178		
179		
180		
181	Best practice	A method or technique that has consistently shown results superior to those achieved with other means, and that is used as a benchmark. The term is also used to describe the process of developing and following a standard way of doing things that multiple organisations can use.
182		
183		
184		
185		
186		
187	Between-laboratory assessment	Phase in which different operators from different laboratories perform (or run) the <i>in vitro</i> method independently to establish whether or not an <i>in vitro</i> method can be successfully established in different laboratories.
188		
189		
190		

² All terms and their descriptions should be considered as working definitions for the purpose of this Guidance Document only.

191	Biokinetics	Time-course of a chemical (substance and mixture) and its metabolites in a living organism, i.e., increase or decrease of substance concentration at the site of measurement due to transport or due to formation or breakdown. The term "toxicokinetics" is also often used synonymously.
192		
193		
194		
195		
196	Biological pathway	A number of biochemical steps, linked together in a certain order, with a start and an end. Some typical types of biological pathways are metabolic pathways and signalling pathways.
197		
198		
199		
200	Blank control/untreated control	Separate untreated part of a test system that is kept under the original culture conditions; the untreated control provides baseline data of the test system under the conditions of the <i>in vitro</i> method. It provides the background response from the test system, obtained by treatment with only the buffer or media used for the administration of test, control and reference items. The primary purpose of a blank or untreated control is to trace sources of artificially introduced variation on the <i>in vitro</i> method results.
201		
202		
203		
204		
205		
206		
207		
208	Bovine spongiform encephalopathy (BSE)	Bovine spongiform encephalopathy, commonly known as mad cow disease, is a fatal neurodegenerative disease (encephalopathy) in cattle that causes a spongy degeneration of the brain and spinal cord.
209		
210		
211		
212		
213	Carcinogenicity	The property of any agent (chemical, physical or biological agent) directly involved in causing cancer (carcinogen). Carcinogenicity results in an increased incidence of tumours, increased proportion of malignant tumours or a reduction in the time to appearance of tumours, compared with concurrent control groups. The process of carcinogenesis involves the transition of normal cells into cancer cells via a sequence of stages that entail both genetic alterations (i.e. mutations) and non-genetic events.
214		
215		
216		
217		
218		
219		
220		
221	Coefficient of variation (CV)	A measure of spread that describes the amount of variability relative to the mean. Because the coefficient of variation is per definition unrelated to the magnitude of the mean and also unitless, it can be used instead of the standard deviation to compare the spread of data sets that have different units or different means.
222		
223		
224		
225		
226	Comparative genomic hybridisation analysis (aCGH)	A molecular cytogenetic method for analysing copy number variations relative to ploidy level in the DNA of a test sample compared to a reference sample, without the need for culturing cells. The aim of this technique is to quickly and efficiently compare two genomic DNA samples arising from two sources, which are most often closely related, because it is suspected that they contain differences in terms of either gains or losses of either whole chromosomes or subchromosomal regions (a portion of a whole chromosome).
227		
228		
229		
230		
231		
232		
233		
234		
235	Computerised systems	A group of hardware components and associated software designed and assembled to perform a specific function or group of functions.
236		
237		
238	Cytotoxicity	General cytotoxicity (or basal cytotoxicity) is the result of toxic effects on structures and functions common to all cells of the body, such as DNA, chromosomes, mitochondria, the cytoskeleton and various membranes. A large number of general (basal) cytotoxicity tests have been developed for a wide range of purposes. Tissue-specific cytotoxicity involves adverse effects on particular types of differentiated cells, in terms of either their tissue-type specific properties (e.g. hormone production, conductivity, contractility) or their modulation of general cytotoxic effects (e.g. via the metabolic activation or detoxification of xenobiotics). For this purpose target organ(system)toxicity tests have been developed.
239		
240		
241		
242		
243		
244		
245		
246		
247		
248		
249		
250	Data	Information derived or obtained from raw data (e.g. a reported analytical result)
251		
252	Data Governance	The total sum of arrangements to ensure that data, irrespective of the format in which it is generated, is recorded, processed, retained and used to ensure a complete, consistent and accurate record throughout the data lifecycle
253		
254		
255		
256	Data Integrity	The extent to which all data are complete, consistent and accurate throughout the data lifecycle.
257		

258		
259	Data Lifecycle	All phases in the life of the data (including raw data) from initial generation and recording through processing (including transformation or migration), use, data retention, archive / retrieval and destruction.
260		
261		
262		
263	Design qualification (DQ), installation qualification (IQ), operational qualification (OQ), performance qualification (PQ)	Verification of <i>in vitro</i> method equipment usually consists of design qualification installation qualification, operational qualification and performance qualification. When acquiring a new piece of equipment in an <i>in vitro</i> environment, design specifications are needed for the intended use of the equipment. Installation qualification is the first step in the process to ensure that the equipment will live up to its expectations. Operational qualification verifies that the equipment is achieving its operational requirements. For the performance qualification, the last phase, the equipment will run several times under normal operating conditions and its functions will be challenged.
264		
265		
266		
267		
268		
269		
270		
271		
272		
273		
274	ECVAM DataBase service on ALternative Methods (DB-ALM)	A database aimed at implementing the communication and dissemination strategy of the Joint Research Centre on animal alternatives, as requested by the European Commission and the European Parliament.
275		
276		
277		
278		
279	Effective concentration 50 (EC₅₀) and Inhibition concentration 50 (IC₅₀)	For <i>in vitro</i> cell and tissue culture work the terms effective concentration 50 (EC ₅₀) and inhibition concentration 50 (IC ₅₀) are used, in analogy to median effective dose (ED ₅₀) and median lethal dose (LD ₅₀) used in animal experiments. IC ₅₀ is used in case of an <i>in vitro</i> method where there is a decline in read-out. IC ₅₀ is therefore the test item concentration causing 50% inhibition of the desired activity. EC ₅₀ is used for read-outs that increase with concentration. EC ₅₀ is therefore the concentration causing 50% of maximum effect for any measured biological effect of interest.
280		
281		
282		
283		
284		
285		
286		
287		
288		
289	Emulsion	A stable dispersion of liquid droplets in another liquid, where the two are immiscible.
290		
291		
292	Engelbreth-Holm-Swarm (EHS)	A mouse sarcoma which is a rich source of both individual basement membrane components and Matrigel often used in cell and tissue culture work.
293		
294		
295	European Chemicals Agency (ECHA)	Agency of the European Union (EU) that manages technical, scientific and administrative aspects of EU chemicals legislation, notably the regulation on the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH).
296		
297		
298		
299	European Directorate for the Quality of Medicines & HealthCare (EDQM)	Organisation that is responsible for the European Pharmacopoeia and the European biological standardisation programme.
300		
301		
302		
303		
304	European Food Safety Authority (EFSA)	Agency of the European Union that provides independent scientific advice in the fields of food and feed safety, animal health and welfare, plant protection and plant health and communicates on existing and emerging risks associated with the food chain.
305		
306		
307		
308		
309		
310	European Medicines Agency (EMA)	Agency of the European Union that is responsible for the protection of public and animal health through the scientific evaluation and supervision of medicines.
311		
312		
313	European Union Network of Laboratories for the Validation of Alternative Methods (EU-NETVAL)	A network of highly qualified laboratories to (1) respond to some of the provisions of Directive 2010/63/EU, (2) generate <i>in vitro</i> method information that is reliable, relevant and based on current best quality and scientific practices, (3) increase the European Commission's validation capacity of <i>in vitro</i> methods and (4) provide a laboratory network knowledgeable on the routine implementation of good <i>in vitro</i> method practices for regulatory use in human safety assessment.
314		
315		
316		
317		
318		
319		
320		
321		
322	Foetal bovine serum (FBS)	Fetal bovine serum derived from clotted blood is the most widely used growth supplement for cell and tissue culture media because of its high content of embryonic growth promoting factors. When used at appropriate concentrations it
323		
324		

325		may supply many defined and undefined components that have been shown to
326		satisfy specific metabolic requirements for the culture of cells and tissues.
327		
328	Genetically modified micro-organisms	
329	(GMMs),	A micro-organism in which the genetic material has been altered in a way that
330		does not occur naturally by mating and/or natural recombination.
331		
332		
333	Good cell culture practice (GCCP) and	
334	GCCP2.0	Guidelines developed in 2005 to define minimum standards in cell and tissue
335		culture work. This GCCP guidance lists a set of six principles intended to support
336		best practice in all aspects of the use of cells and tissues <i>in vitro</i> , and to
337		complement, but not to replace, any existing guidance, guidelines or regulations.
338		GCCP2.0 is the updated version which is currently being drafted.
339		
340	Good Laboratory Practice (GLP)	A quality system applied to the conduct of non-clinical health and environmental
341		safety testing that is intended for submission to regulatory authorities in support
342		of the registration, licensing or regulation of chemical and related products. It
343		concerns the organisational process and the conditions under which these studies
344		are planned, performed, monitored, recorded, archived and reported. It ensures
345		uniformity, consistency, reproducibility, quality and integrity of chemical non-
346		clinical safety tests.
347		
348	Hazard	An intrinsic feature of a stressor (e.g. chemical or physical in nature) to cause
349		harm or adverse effects to human health and to the environment. It is a qualitative
350		(for example in the case of classifications) or quantitative expression of the
351		adverse effects elicited by a test item under defined conditions of exposure.
352		
353		
354	High performance liquid chromatography	
355	(HPLC)	High performance liquid chromatography (or high-pressure liquid
356		chromatography) is a chromatographic technique that can separate a mixture of
357		compounds when in solution and is used in biochemistry and analytical chemistry
358		to identify, quantify and purify the individual components of the mixture.
359		
360	High-efficiency particulate arrestance	
361	(HEPA)	High-efficiency particulate arrestance, also sometimes called high-efficiency
362		particulate arresting or high-efficiency particulate air, is a type of air filter used to
363		create an aseptic environment, through retention of a certain number of particles,
364		depending of the category of the filter. Aseptic conditions are required for <i>in vitro</i>
365		cell and tissue culture work to ensure reliability and reproducibility of results
366		obtained.
367		
368	High-throughput screening (HTS)	A high throughput method for scientific experimentation, relevant to the fields of
369		biology and chemistry, and especially used in drug discovery. HTS involves an
370		automated-operation platform, data processing and control software. HTS allows
371		a researcher to quickly conduct many biochemical, genetic or pharmacological
372		tests and obtain a large number of information from a single experiment.
373		
374	<i>In silico</i>	The term <i>in silico</i> refers to the technique of performing experiments via computer
375		simulations. Examples are structure-activity relationships (SAR) and quantitative
376		structure-activity relationships (QSAR)
377		
378	<i>In vitro</i>	The term <i>in vitro</i> (Latin for "in the glass") refers to the technique of performing a
379		given experiment in a test tube, or, more generally, in a controlled environment
380		outside of a living organism.
381		
382	<i>In vitro</i> method endpoint	Defined as quantitative measurable characteristics that serve as indicators of a
383		pathologic process or related biochemical or molecular events, e.g. measured
384		absorbance in a cytotoxicity assay or a skin irritation <i>in vitro</i> method.
385		
386	<i>In vitro</i> to <i>in vivo</i> extrapolation (IVIVE)	<i>In vitro</i> to <i>in vivo</i> extrapolation refers to the qualitative or quantitative
387		transposition of experimental results or observations made <i>in vitro</i> to predict
388		phenomena <i>in vivo</i> , i.e. in whole organisms.
389		
390	<i>In vivo</i>	
391		

392		The term <i>in vivo</i> ("within the living") refers to experimentation using a whole, living organism as opposed to a partial or dead organism, or an <i>in vitro</i> controlled environment. Animal testing and clinical trials are two forms of <i>in vivo</i> research.
393		
394		
395		
396	Inhibitor or spiked up control	Mix of test item and positive control to assess any effect of inhibition of the test item on the test system endpoint measurements.
397		
398		
399	Integrated testing strategies (ITS)	Integrated testing strategies enable to significantly increase the use of non-animal testing information for regulatory decision making, and thus to minimise the need for animal testing. To this end, operational procedures are developed, tested and disseminated that guide a transparent and scientifically sound evaluation of test items in a risk-driven, context-specific and substance-tailored manner. The envisaged decision theory framework includes alternative methods such as chemical and biological read-across, <i>in vitro</i> results, <i>in vivo</i> information on analogues, qualitative and quantitative structure-activity relationships, thresholds of toxicological concern and exposure-based waiving.
400		
401		
402		
403		
404		
405		
406		
407		
408		
409	Intellectual property rights (IPR)	Intellectual property refers to creations of the mind: inventions; literary and artistic works; and symbols, names and images used in commerce. Industrial Property includes patents for inventions, trademarks, industrial designs and geographical indications.
410		
411		
412		
413		
414	International Uniform Chemical Information Database (IUCLID)	A software application designed to capture, store, maintain and exchange data on intrinsic and hazard properties of chemicals (substances and mixtures). It is essential for chemical industry to comply with the new legislation which entered into force on 1 June 2007. The freely downloadable tool will assist chemical companies globally in fulfilling their obligation to submit data to the Agency under the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH) legislation from 1 June 2008.
415		
416		
417		
418		
419		
420		
421		
422		
423		
424	Limit of detection (LOD), Lower limit of quantification (LLOQ) and Upper limit of quantification (ULOQ)	The Limit of detection and the upper and lower limits of quantification are important parameters that need to be determined during <i>in vitro</i> method development. The LOD is the lowest amount of test item in a sample that can be detected with (stated) probability, although not quantified as an exact value. The LLOQ and ULOQ are the lowest and highest amounts of test item in a sample that can be quantitatively determined with a stated acceptable precision and accuracy, under stated experimental conditions. The usefulness and optimal throughput of an <i>in vitro</i> method may depend on the appropriate determination of the LOD and the ULOQ and LLOQ.
425		
426		
427		
428		
429		
430		
431		
432		
433		
434		
435	Lipophilicity	The ability of a chemical (substance and mixture) to dissolve in non-polar environments such as oils, lipid membranes, and non-polar solvents such as hexane or toluene.
436		
437		
438		
439	Mass spectrometry (MS)	Analytical technique that measures the mass-to-charge ratio of charged particles. It is used for determining masses of particles, for determining the elemental composition of a sample or molecule, and for elucidating the chemical structures of molecules such as peptides and other chemical compounds.
440		
441		
442		
443	Maximum average score (MAS)	The maximum average Draize score was utilised as the primary quantitative measurement of eye irritation potential in rabbits.
444		
445		
446	Metadata	Metadata is data that describe the attributes of other data, and provide context and meaning. Typically, these are data that describe the structure, data elements, inter-relationships and other characteristics of data. It also permits data to be attributable to an individual.
447		
448		
449		
450		
451	Micro-organism	Any microbiological entity, cellular or non-cellular, capable of replication or of transferring genetic material, including viruses, viroids, animal and plant cells in culture.
452		
453		
454		
455	Minimal essential medium (MEM)	Developed by Harry Eagle, is one of the most widely used of all synthetic cell culture media for <i>in vitro</i> cell and tissue culture work.
456		
457		

458	Minimum significant ratio (MSR)	Parameter that can be used to quantify assay reproducibility and resolution (the smallest ratio between compound potency which can be detected in the <i>in vitro</i> method).
459		
460		
461	Mixture²	A combination of two or more chemicals (liquid or solid) that do not react with each other
462		
463		
464		
465	Multi-component test chemicals³	Mixtures comprising a complex mix of individual test chemicals with different solubility and physical-chemical properties. In most cases, they can be characterised as a homologous series of test chemicals with a certain range of carbon chain length/number or degree of substitution. These materials are frequently referred to as “complex mixtures”. However, in this guidance document, these are referred to as “multi-component test chemicals”.
466		
467		
468		
469		
470		
471		
472	Mutual Acceptance of Data (MAD)	The OECD Mutual Acceptance of Data is a multilateral agreement which states that test data generated in any member country in accordance with OECD Test Guidelines and GLP shall be accepted in other member countries for assessment purposes and other uses relating to the protection of human health and the environment. The application of MAD avoids unnecessary and costly duplication of testing as well as non-tariff barriers to trade. In addition, it saves laboratory animals used for <i>in vivo</i> testing.
473		
474		
475		
476		
477		
478		
479		
480	Nanomaterials	A natural, incidental or manufactured material containing particles, in an unbound state or as an aggregate or as an agglomerate and where, for 50 % or more of the particles in the number size distribution, one or more external dimensions is in the size range 1 nm - 100 nm.
481		
482		
483		
484		
485	Negative control	Separate part of a test system treated with an item for which it is known that the test system should not respond; the negative control provides evidence that the test system is not responsive under the actual conditions of the assay.
486		
487		
488		
489	Omics	Omics is a general term for a broad discipline of science and engineering for analysing the interactions of biological information objects in various omes (these include genome, transcriptome, proteome, metabolome, expressome, and interactome). Some examples of ‘Omics’ technologies: - genomics - proteomics - metabolomics - transcriptomics
490		
491		
492		
493		
494		
495		
496		
497		
498		
499	Particulates	Also known as particulate matter (PM), fine particles and soot, are tiny subdivisions of solid matter suspended in a gas or liquid. In contrast, aerosol refers to particles and/or liquid droplets and the gas together. Sources of particulate matter can be man-made or natural. Air pollution and water pollution can take the form of solid particulate matter or be dissolved. Salt is an example of a dissolved contaminant in water, while sand is generally a solid particulate.
500		
501		
502		
503		
504		
505		
506	Physiologically based pharmacokinetic, physiologically based toxicokinetic, physiologically based biokinetic (PBPK/PBTK/PBBK)	Physiologically based toxicokinetic, or alternatively referred to as physiologically based pharmacokinetic or biokinetic models, are quantitative descriptions of absorption, distribution, metabolism, and excretion (ADME) of synthetic or natural chemical substances in humans and other animal species. PBTK models are increasingly being used as an effective tool for designing toxicology experiments and for conducting extrapolations essential for risk assessments (e.g. in pharmaceutical research and drug development, and in health risk assessment for cosmetics or general chemicals).
507		
508		
509		
510		
511		
512		
513		
514		
515		
516	Polymerase chain reaction (PCR)	Polymerase chain reaction is a molecular biology <i>in vitro</i> technique. Using the natural ability of DNA polymerase to synthesise a new strand of DNA
517		

³ Consideration is given to the consistency with the definition of “multi-component test chemicals” (or “complex test chemicals”) in Chapter 1 of a draft “Guidance Document on the Use of the Harmonized System for the Classification of Chemicals which are Hazardous for the Aquatic Environment” (ENV/JM/HCL(2000)11).

518		complementary to the offered template strand, millions of copies of a specific DNA sequence are made.
519		
520		
521	Positive control	Separate part of the test system treated with an item the response to which is known for the test system; the positive control provides evidence that the test system is responsive under the actual conditions of the assay.
522		
523		
524		
525	Prediction Model	The method by which the <i>in vitro</i> endpoint value(s) is used to predict the <i>in vivo</i> equivalent activity (i.e., degree of toxicity).
526		
527		
528	Proficiency chemicals	A panel of chemicals with known and statistically well-defined responses in a particular <i>in vitro</i> method. These are items used e.g. to verify that a laboratory carries out the <i>in vitro</i> method correctly, or to validate alternative newly developed test systems intended for use with the particular <i>in vitro</i> method.
529		
530		
531		
532		
533	Quality assurance	A system of procedures, checks, audits, and corrective actions to ensure that all research, testing, monitoring, sampling, analysis, and other technical and reporting activities are of the highest achievable quality.
534		
535		
536		
537		
538	Quality assurance programme	A defined system, including personnel, which is independent of study conduct and is designed to assure test facility management of compliance with GLP.
539		
540		
541	Quality control	Documented activity which seeks to confirm that starting materials for assays (including cell lines) and key stages of individual assays and the final results, meet prescribed specifications. It should be clear that this is not the same as quality assurance. ISO 9000 defines quality control as "part of quality management focused on providing confidence that quality requirements will be fulfilled".
542		
543		
544		
545		
546		
547		
548	Quality management system (QMS)	Can be expressed as the organisational structure, procedures, processes and resources needed to implement quality management. GLP specifically refers to a quality system of management controls for test facilities and organisations to try to ensure the uniformity, consistency, reliability, reproducibility, quality, and integrity of test item non-clinical safety tests. Of all QMS regimes, the ISO 9000 family of standards is probably the most widely implemented worldwide.
549		
550		
551		
552		
553		
554		
555	Raw data	Original records and documentation, retained in the format in which they were originally generated (i.e. paper or electronic), or as a 'true copy'. Raw data must be contemporaneously and accurately recorded by permanent means. In the case of basic electronic equipment which does not store electronic data, or provides only a printed data output (e.g. balance or pH meter), the printout constitutes the raw data.
556		
557		
558		
559		
560		
561		
562	Reagent	Term used for media additives, compounds added to a system to induce a chemical reaction, anything added to get the <i>in vitro</i> method or related assays to work etc.
563		
564		
565		
566	Reference item	Any chemical (substance and mixture) or product used to provide a basis for comparison with the test item reference items are used e.g. during method development to verify that the method classifies test items correctly, or in every experimental run of a finalised method to verify system performance according to the acceptance criteria.
567		
568		
569		
570		
571		
572	Relevance	The term "Relevance" describes whether a procedure is meaningful and useful for a particular purpose.
573		
574		
575	Reliability	The term "Reliability" describes whether a procedure can be performed reproducibly within and between laboratories and over time.
576		
577		
578	Replace, Reduce, Refine (3Rs)	3Rs is the short for "Replace, Reduce, Refine". A term describing current internationally accepted strategies for minimising suffering of laboratory animals used in experimental research. The optimal solution is to replace the test method requiring animal experiments with one or several <i>in vitro</i> methods; if this is not possible at least it might be possible to modify the methods so that it is possible to Reduce the number of animals being used in each study without compromising data quality; if this is also not possible it might at least be possible to Refine the
579		
580		
581		
582		
583		
584		

585		test method so that experiments are conducted in a way minimising stress and
586		other impact on the animals.
587		
588	Robustness	The insensitivity of test results to departures from the specified test conditions
589		when conducted in different laboratories or over a range of conditions under
590		which the test method might normally be used. If a test is not robust, it will be
591		difficult to use in a reproducible manner within and between laboratories.
592		
593	Safe Harbour	The international Safe Harbour Privacy Principles or Safe Harbour Privacy
594		Principles were principles which were overturned on October 24, 2015 by the
595		European Court of Justice, which enabled some US companies to comply with
596		privacy laws protecting European Union and Swiss citizens. US companies
597		storing customer data would self-certify that they adhere to 7 principles, to
598		comply with the EU Data Protection Directive and with Swiss requirements. The
599		US Department of Commerce developed privacy frameworks in conjunction with
600		both the European Union and the Federal Data Protection and Information
601		Commissioner of Switzerland.
602		Within the context of a series of decisions on the adequacy of the protection of
603		personal data transferred to other countries, the European Commission made a
604		decision in 2000 that the United States' principles complied with the EU
605		Directive - the so-called "Safe Harbour decision". However, after a customer
606		complained that his Facebook data were insufficiently protected, the European
607		Court of Justice declared in October 2015 that the Safe Harbour Decision was
608		invalid, leading to further talks being held by the Commission with the US
609		authorities towards "a renewed and sound framework for transatlantic data
610		flows".
611		The European Commission and the United States agreed to establish a new
612		framework for transatlantic data flows on 2nd February 2016, known as the "EU-
613		US Privacy Shield".
614		
615	Saturation concentration	The maximum dissolved concentration of a test chemical that can be achieved
616		under the test conditions.
617		
618	Sensitivity	A measure of <i>in vitro</i> method performance that describes the proportion of all
619		evaluated test items that are classified as positive for a particular toxicological
620		endpoint, which are predicted as positive by the actual <i>in vitro</i> method.
621		
622	Service level agreement (SLA)	A contract between a service provider (either internal or external) and the end
623		user that defines the level of service expected from the service provider.
624		
625	Short tandem repeat (STR)	Short Tandem Repeat (STR) DNA profiling, is used to uniquely identify human
626		cell lines derived from the tissue of a single individual allowing researchers to
627		ascertain if their cultures were misidentified or cross-contaminated.
628		
629	Signal windows (SW)	A measure of the separation between the sample (positive control) and the blank
630		(negative control/solvent) including the variability of both measurements.
631		
632	Single nucleotide polymorphism analysis	
633	(aSNP)	Single nucleotide polymorphism or SNP (pronounced snip) analysis is a
634		technique to detect a DNA sequence variation occurring when a single nucleotide
635		- A, T, C, or G - in the genome (or other shared sequence) differs between
636		members of a species (or between paired chromosomes in an individual). For
637		example, two sequenced DNA fragments from different individuals, AAGCCTA
638		to AAGCTTA, contain a difference in a single nucleotide.
639		
640	Solid phase microextraction (SPME)	Is a solid phase extraction sampling technique that involves the use of a fiber
641		coated with an extracting phase, that can be a liquid (polymer) or a solid
642		(sorbent), which extracts different kinds of analytes (including both volatile and
643		non-volatile) from different kinds of media.
644		
645	Solubility limit in water	The maximum attainable concentration or concentration at thermodynamic
646		equilibrium between aqueous pure phase and solid (or liquid or gaseous) pure
647		phase.
648		
649	Specificity	A measure of <i>in vitro</i> method performance that describes the proportion of all
650		evaluated test items that are classified as negative for a particular toxicological
651		endpoint, which are predicted as negative by the actual <i>in vitro</i> method.

652		
653	Standard deviation (SD)	A measure that is used to quantify the amount of variation or dispersion of a set of data values.
654		
655		
656	Standard operating procedure (SOP)	A documented procedure which describes how to perform tests or activities normally not specified in detail in study plans or test guidelines.
657		
658		
659	Structure-activity relationships and quantitative structure-activity relationships (SAR/QSAR)	Structure-activity relationships and quantitative structure-activity relationships, collectively referred to as (Q)SARs, are simplified mathematical representations of complex chemical-biological interactions that can be used to predict the physicochemical and biological properties of molecules.
660		
661		
662		
663		
664		
665	Study plan	A document which defines the objectives and experimental design for the conduct of the study, and includes amendments (i.e. an intended change to the study plan after the study initiation date).
666		
667		
668		
669	Suspension	A stable dispersion of solid particles in a liquid.
670		
671	Test item	A chemical (substance and mixture) or product that is the subject of a study.
672		
673	Test pre-submission form (TPF) and Test submission template (TST)	For the evaluation of the readiness of an <i>in vitro</i> method to enter the EURL ECVAM validation process, the method needs to be officially submitted to EURL ECVAM by compiling as a first step the electronic version of the test pre-submission form to allow a preliminary assessment of the status of development, optimisation and/or validation of an <i>in vitro</i> method and its potential relevance with regard to the 3Rs (replacement, reduction, refinement of animal testing). If this step is satisfactory a complete submission is requested which requires the compilation of a detailed test submission template.
674		
675		
676		
677		
678		
679		
680		
681		
682		
683	Test system	A test system means any biological, chemical or physical system or a combination thereof used in a study. <i>In vitro</i> test systems are mainly biological systems (e.g. cells or tissues), although some of the more recent developments in alternatives to conventional <i>in vivo</i> testing (e.g., gene arrays for toxicogenomics) may also exhibit some attributes of physical-chemical test systems, and still others, e.g., toxicometabonomics, may mainly rely on analytical methodology. Test kits, including proprietary test kits, should also be considered as test systems.
684		
685		
686		
687		
688		
689		
690		
691		
692	Testing method	The term testing method is used to describe both a published, well-evaluated test guideline method (e.g. OECD) and a not-fully-developed method soon ready for submission to validation bodies. A testing method is a process or procedure used to obtain information on the characteristic of a substance or agent. Toxicological testing methods generate information regarding the ability of a substance or agent to produce a specific biological effect under specified conditions.
693		
694		
695		
696		
697		
698		
699	Toxicological endpoint	A direct marker of progression to an adverse outcome - e.g. morphological or physiological changes, functional impairments, disease symptoms or death - used to describe an adverse health effect (or a probability of that adverse effect) resulting from exposure to a test item. The test system response to an exposure of a test item may be measured by a series of endpoints. The most sensitive endpoint (critical endpoint) is the one that occurs at the lowest exposure level and associated with an adverse response (committed step).
700		
701		
702		
703		
704		
705		
706		
707	Training Set	The set of test items used to develop the prediction model for an assay. The training set items should have strong reference data (i.e., values from a recognised regulatory assay) against which the <i>in vitro</i> assay endpoint values can be compared.
708		
709		
710		
711		
712	Untreated control	Separate untreated part of a test system that is kept under the original culture conditions; the untreated control provides baseline data of the test system under the conditions of the assay.
713		
714		
715		
716	Validation	Validation is the process by which the reliability and relevance of a procedure are established for a specific purpose.
717		
718		

719	Validation set (test items)	The set of test items used to assess the predictive capacity of an <i>in vitro</i> method based on the performance of the endpoint values by the reference test results.
720		
721		Testing of the validation test items set is a principal part of <i>in vitro</i> method
722		validation.
723		
724	Vehicle or solvent control	The separate part of a test system to which the vehicle (i.e. solvent) for the test
725		item is added without the test item; the vehicle control provides evidence for a
726		lack of influence of the chosen vehicle on the test system under the actual
727		conditions of the <i>in vitro</i> method.
728		
729	Within-laboratory assessment	Phase in which different operators from the same laboratory perform (or run) the
730		<i>in vitro</i> method independently and at different times to establish whether or not
731		an <i>in vitro</i> method can be successfully established in one laboratory.
732		
733	Xenobiotic	A chemical foreign to the biological system, structurally distinct from
734		endogenous compounds present within the biological system.
735		
736	Z-factor	A measure of the separation between control and sample signal which takes into
737		account the dynamic range of the <i>in vitro</i> method and the data variation
738		associated with the signal and control measurements. It is suitable for <i>in vitro</i>
739		method quality assessment.

DRAFT

740 **Introduction**

741 The availability of a guidance document on Good *In vitro* Method Practices (GIVIMP) for
742 "The development and implementation of *in vitro* methods for regulatory use in human safety
743 assessment" was identified by the scientific and regulatory community (OECD) as a high
744 priority requirement to reduce the uncertainties in cell and tissue-based *in vitro* method
745 derived predictions.

746 An Expert Group was therefore established to develop such a guidance document. The first
747 draft guidance document was prepared following a GIVIMP meeting on the 24th and 25th of
748 February 2015 in Ispra, Italy (see Annex 3) with additional input from experts who could not
749 be present at the meeting. For this first draft version expert input was received from EURL
750 ECVAM, European receiving authorities (European Food Safety Authority EFSA, European
751 Medicine Agency EMA, the European Chemicals Agency ECHA), from the European Union
752 Network of Laboratories for the Validation of Alternative Methods (EU-NETVAL, e.g. from
753 the Belgian, Dutch, Italian, Spanish and Swedish laboratories), from ECVAM's Stakeholder
754 Forum (ESTAF, e.g. the European Society of *In vitro* Toxicology), from the EU and OECD
755 Working Group on GLP (e.g. delegates from Belgium, The Netherlands, The United
756 Kingdom, Poland, Italy, France, Singapore), from 3Rs Centres (Centre for Alternatives to
757 Animal Testing, CAAT), from regulatory agencies (e.g. RIVM), from scientists from large
758 industries and SMEs and from international scientists with expertise in stem cells, cell
759 biology, Good Laboratory Practice (GLP) and *in vitro* methods.

760
761 The main authors drafting this first version are:
762

763 Maria Antonietta Antonelli¹, Maria Astridou², Bohumila Bichlmaier Suchanová³, Camilla
764 Bernasconi⁴, Ann-Charlotte Bostroem⁴, Gerard Bowe⁴, Sandra Coecke⁴, Thomas Cole⁴, Hui
765 Da Chong⁵, Rita Cortvrindt⁶, Laureano Cuevas⁷, Agnieszka Dudra⁸, Eugene Elmore⁹,
766 Chantra Eskes¹⁰, Kristina Fant¹¹, Andrew Gray¹², John Harbell¹³, Minne Heringa¹⁴, Delyan
767 Ivanov¹⁵, Guido Jacobs¹⁶, Rob Jaspers¹⁷, Elisabeth Joossens⁴, Ingrid Langezaal⁴, Marcel
768 Leist¹⁸, Thomas Lucotte¹⁹, Richard McFarland²⁰, Marisa Meloni²¹, Anne Milcamps⁴, Emma
769 Pedersen¹¹, Laura Pioppo¹, Maurits-Jan Prinz²², Annamaria Rossi², Laura Rossi³, Glyn N
770 Stacey²³, Emanuela Testai²⁴, Frederique van Acker²⁵, Sander van der Linden⁴, Maurice
771 Whelan⁴, Iwona Wilk-Zasadna²⁶

772
773 ¹European Medicines Agency (EMA), London, United Kingdom.

774 ²European Food Safety Authority (EFSA), Parma, Italy.

775 ³European Chemicals Agency (ECHA), Helsinki, Finland.

776 ⁴European Commission, Directorate General Joint Research Centre, Directorate F-Health,
777 Consumers and Reference Materials, F.3 Chemicals Safety and alternative Methods, The
778 European Union Reference Laboratory for Alternatives to Animal Testing (EURL ECVAM),
779 Ispra, Italy.

780 ⁵Singapore GLP Monitoring Authority, SPRING Singapore — Singapore Accreditation
781 Council.

782 ⁶Vito, Mol, Belgium.

783 ⁷Centro Nacional de Microbiologia, Majadahonda, Spain.

784 ⁸Bureau for Chemical Substances, Lodz, Poland.

785 ⁹University of California Irvine, Trabuco Canyon, Unites States.

786 ¹⁰SeCAM, Magliaso, Switzerland.

787 ¹¹SP Sveriges Tekniska Forskningsinstitut, Borås, Sweden
788 ¹²United Kingdom GLP Monitoring Authority, The Medicines and Healthcare products
789 Regulatory Agency (MHRA), London, United Kingdom
790 ¹³JHarbell Consulting LCC, Dallas, United States.
791 ¹⁴RIVM, The National Institute for Public Health and the Environment, Centre for Safety of
792 Substances and Products, Bilthoven, The Netherlands.
793 ¹⁵University of Nottingham, Cancer Biology Queens Medical Centre, Nottingham, United
794 Kingdom.
795 ¹⁶WIV-ISP, Brussels, Belgium.
796 ¹⁷Health Care Inspectorate, Ministry of Health, Welfare and Sport , Utrecht, The Netherlands.
797 ¹⁸University of Konstanz, Konstanz, Germany.
798 ¹⁹ANSM, Saint-Denis, France.
799 ²⁰Policy Office of Cellular, Tissue, and Gene Therapies Center for Biologics Evaluation and
800 Research United States Food and Drug Administration, USA
801 ²¹Vitroscreen SRL, Milano, Italy.
802 ²²European Commission, Directorate General Internal Market, Industry, Entrepreneurship
803 and SMEs, Brussels, Belgium.
804 ²³United Kingdom Stem Cell Bank, Advanced Therapies Division, NIBSC, Medicines and
805 Healthcare products Agency, United Kingdom.
806 ²⁴Istituto Superiore di Sanita, Rome, Italy.
807 ²⁵Triskelion, Zeist, The Netherlands.
808 ²⁶Novartis Vaccines, Basel, Switzerland.
809

810 Acknowledgements are also made to the EURL ECVAM colleagues that reviewed this first
811 version and gave useful additional input⁴.

812 This first draft, currently for revision by the OECD Working Group on GLP and the
813 nominated experts from the WNT, has been be circulated in September for review by all 37
814 members of the European Union Network of Laboratories for the Validation of Alternative
815 Methods (EU-NETVAL). An updated version revised draft will be submitted at the beginning
816 of January for a second OECD commenting round. EURL ECVAM will then prepare the
817 final version for proposed adoption at the OECD meeting in April 2017.

⁴ Susanne Belz, Mounir Bouhifd, Laura Gribaldo, Tomislav Horvat, Tracey Holley, Annett Janusch Roi, Roman Liska, Alfonso Lostia, Agnieszka Lidia Swiatek, Francesca Pistollata, Anna Price, Clemens Wittwehr, Andrew Worth.

818

Scope

819 There is a community desire for non-animal methods, but regulators demand validated and
820 internationally accepted *in vitro* methods (i.e. OECD test guidelines or ISO standards). To
821 accommodate the desires of regulatory authorities, a number of *in vitro* methods, often based
822 on the use of human cells and tissues, were submitted to international validation bodies
823 during the last two decennia. However, the experience gained during these validations
824 revealed that many *in vitro* methods need serious improvements in design, robustness and
825 reliability before they can be successfully implemented in a routine laboratory environment
826 and generate data sets which can be used to support regulatory decisions. Therefore, OECD
827 approached EURL ECVAM to coordinate the issuing of a guidance on Good *In Vitro* Method
828 Practices (GIVIMP) for the development and implementation of *in vitro* methods for
829 regulatory use in human safety assessment. The major goal of GIVIMP consists of improving
830 the reliability and robustness of *in vitro* methods, reducing the uncertainties of *in vitro* based
831 predictions and therefore increasing the acceptance of the *in vitro* estimated safety measures
832 by regulatory agencies. The scope of the GIVIMP guidance is taking into account good
833 scientific, technical and quality practices, to ensure that the overall process, starting from *in*
834 *vitro* method development up to the final *in vitro* method implementation for regulatory use
835 becomes more efficient and effective.

836 This guidance document targets all players involved in the process, e.g. *in vitro* method
837 developers, *in vitro* test system producers, validation bodies, producers of equipment,
838 materials and reagents, *in vitro* method end-users such as EU-NETVAL test facilities, testing
839 laboratories, large industries and small to medium enterprises as well as receiving authorities,
840 monitoring authorities, accreditation bodies and OECD. The guidance aims to further
841 facilitate the application of the OECD Mutual Acceptance of Data agreement for data
842 generated by *in vitro* methods and as such contribute to avoidance of unnecessary duplicate
843 testing. This guidance describes the areas related to *in vitro* method development,
844 standardisation, harmonisation, and international acceptance that would benefit from more
845 detailed scientific, technical and quality guidance.

846 The GIVIMP document has been written with different end users in mind, including GLP
847 routine test facilities but also research laboratories developing new *in vitro* methods. In the
848 latter case it is obvious that the laboratory will not be able to fully comply with this document
849 on all points. However, it is necessary to comply with a set of good practices within the *in*
850 *vitro* method life cycle so as not to jeopardise the acceptance and routine use of the *in vitro*
851 method in a regulatory environment.

852 This guidance is not intended to duplicate or replace existing OECD guidance or advisory
853 documents but rather it is complementary, addresses specific gaps and aims to collect
854 available references and information on best scientific, technical and quality practices in one
855 document. GIVIMP takes into account the requirements of the existing OECD guidelines and
856 advisory documents to ensure that the guidance is complementary and fully in line with these
857 documents.

858

859 This document is divided into 10 sections covering:

- 860 1 Responsibilities
- 861 2 Quality considerations
- 862 3 Facilities
- 863 4 Apparatus, material and reagents
- 864 5 Test systems

865	6	Test and reference items
866	7	Standard operating procedures (SOPs)
867	8	Performance of the method
868	9	Reporting of results
869	10	Storage and retention of records and materials
870		

DRAFT

871 1 Responsibilities

872 1.1 *In vitro* method developers

873 *In vitro* methods are often developed without the primary aim of being used for regulatory
874 purposes, but are rather focused on the discovery of disease pathways or investigation of
875 mechanisms of action induced by external factors causing cell disturbance. However, these *in*
876 *vitro* methods in development can form the basis for *in vitro* methods for specific toxicity
877 endpoints, during drug and/or other chemicals safety assessment or for toxicity screening
878 during product quality control processes.

879 Researchers aiming to develop *in vitro* methods suitable for regulatory testing purposes must
880 be aware that beyond the ‘short-term’ repeatability, that is consonant with the good scientific
881 work required in discovery, the quality principles for test acceptance by receiving authorities
882 have additional requirements (OECD, 2005a). To harmonise and speed up the validation
883 process and accelerate the acceptance of new *in vitro* methods by receiving authorities, the *in*
884 *vitro* method developer should keep in mind that the quality of historical data and
885 documentation regarding the *in vitro* method submitted will have a significant impact on the
886 validation process.

887 Briefly, the *in vitro* method developer is responsible for providing a clearly written and well
888 documented *in vitro* method description, and related standard operating procedure(s)
889 (SOP(s)), considering all aspects described in the present guidance document.

890 The developer's knowledge and understanding of the *in vitro* method is the basis for
891 establishing an approach to control the *in vitro* method and to set for instance adequate
892 acceptance criteria for the results obtained when running an *in vitro* method.

893 *In vitro* method developers should:

- 894 • Understand the sources of variation of the *in vitro* method
- 895 • Detect the presence and degree of variation in the results
- 896 • Understand the impact of variation of the *in vitro* results on the related predictions
- 897 • Control the variation in a manner to make a sound, relevant and reliable *in vitro*
898 method

899 Each developer should judge whether he or she has gained sufficient understanding of the *in*
900 *vitro* method to provide a high degree of assurance to successfully propose the *in vitro*
901 method for regulatory applications.

902 *In vitro* method developers should also take into account the Intellectual Property (IP)
903 guidelines regarding test systems as set out on the OECD website (<http://www.oecd.org/>).
904 Proposals for projects aiming at the development of new test guidelines should provide
905 information on Intellectual Property Rights (IPR) aspects, as transparently as possible. In
906 particular, the following information is expected to be provided: "Describe if the *in vitro*
907 method includes components, equipment or other scientific procedures that are covered (or
908 pending) by IPR (e.g., patents, patent applications, industrial designs and trademarks) and/or
909 intended to remain confidential. Information should be provided on the overall availability of
910 the IPR-protected components including whether they are commercially available or require a
911 Material Transfer Agreement (MTA) or other licensing agreements. In addition, the
912 possibility of providing a generic description of the IPR-covered component/test system as
913 well as any other element intended to remain confidential should be disclosed and whether
914 Performance Standards have been developed for the *in vitro* method."

915 *In vitro* methods proposed for regulatory use should not contain elements that are confidential
916 to the extent that this impedes adequate scientific validation of the mechanistic relevance of
917 the method.

918 When *in vitro* method developers conclude that their *in vitro* method is sufficiently
919 developed, they can then proceed to an in-house validation. When such internal assessment is
920 successful, they can submit the *in vitro* method to a validation body for the formal validation
921 of the method, or, can organise the validation by themselves. In order to have the *in vitro*
922 method considered for regulatory acceptance, *in vitro* method developers will have to contact
923 their national coordinator to develop a project proposal for a new Test Guideline. Project
924 proposals for new Test Guidelines need the active support of regulatory authorities in at least
925 one member country, and have to meet a regulatory need in member countries.

926 In conclusion, at the end of the test development process, the *in vitro* method developers
927 should be able to prove that the *in vitro* method they offer to the validation body is robust,
928 reliable, relevant, and supported by high quality data as described in the present guidance.

929 **1.2 Test system providers**

930 *In vitro* test systems are mainly biological systems, quite often consisting of tissues or cell
931 lines. It is difficult, if not impossible, to identify cell lines from different origins based solely
932 on morphology and or culture characteristics. Infection or contamination of a cell line with an
933 adventitious virus or mycoplasma may significantly change the characteristics of the cells but
934 again such contamination may not be evident. Cell lines will also change with time in culture,
935 and to add to all these natural hazards it is all too easy to incorrectly label or cross-
936 contaminate different cell lines in a busy cell culture laboratory (ECACC Handbook –
937 Fundamental Techniques for ECACC Cell Lines⁵).

938 The OECD GLP recommends that test system providers should adhere to a formal quality
939 system, such as International Standard ISO 9001, and particularly Part 1 - Specification for
940 Design/Development, Production, Installation and Servicing.

941 The test system provider must also provide all relevant safety information, in compliance
942 with national and international regulations, for the transport, use and disposal, including
943 containment in the case of an accident.

944 The Guidance on Good Cell Culture Practice: A Report of the Second ECVAM Task Force
945 on Good Cell Culture Practice (GCCP) (Coecke *et al.*, 2005) provides a minimal set of
946 information that is essential when working with cells or tissues of animal or human origin
947 (Table 1), while the OECD GLP No 14 (The Application of the Principles of GLP to *in vitro*
948 Studies (OECD, 2004a)) states that the characterisation of test systems can be directly
949 fulfilled by information from the supplier.

950

⁵ <http://www.sigmaaldrich.com/life-science/cell-culture/learning-center/ecacc-handbook.html>

Field Code Changed

951 **Table 1: Examples of requirements for documentation concerning the origins of cells and**
 952 **tissues**

	Isolated organs and tissues of animal origin (e.g. rat brain tissue)	Primary cultures of animal origin (e.g. rat hepatocytes)	All materials of human origin (e.g. cord blood)	Cell lines (e.g. Balb/c, 3T3)
Ethical and safety issues	+	+	+	Applicable, if human or involving recombinant DNA or pathogens
Species/strain	+	+	+	+
Source	+	+	+	+
Sex	+	+	+	+
Age	+	+	+	+
Number of donors	+	+	If applicable	na
Health status	+	+	+	+
Any special pre-treatment	+	+	+	+
Organ/tissue of origin	+	+	+	+
Cell type(s) isolated	+	+	+	+
Isolation technique	+	+	+	+
Date of isolation	+	+	+	+
Operator	+	+	+	+
Supplier	+	+	+	+
Informed consent	na	na	+	If human, may be applicable
Material transfer agreement	na	na	+	+
Medical history of donor	na	na	+(if available)	If human, may be applicable (if available)
Pathogen testing	If applicable ^a	If applicable ^a	+ ^a	+ ^a
Shipping conditions	+	+	+	+
State of material on arrival	+	+	+	+
Cell line identification and authentication	na	na	na	+
Mycoplasma testing	na	na ^b	na ^b	+

^aScreening tests for animal colonies or donors of cells and tissue may be appropriate.

^bMay be important if material is preserved for longer term use (e.g. as feeder layers for other cultures).

na = not applicable.

953

954 1.3 Validation bodies

955 In general, it is the mandate of national and international governmental validation agencies
 956 (ECVAM, ICCVAM, JACVAM, etc.) to promote and facilitate *in vitro* method validation for
 957 regulatory acceptance to replace or reduce *in vivo* animal testing.

958 The validation body's responsibility is to contribute to both an effective validation process
 959 and to *in vitro* method quality. The basic principle of validation is that an *in vitro* method
 960 should be produced that is fit for its intended use. The validation process consists of
 961 collection and evaluation of data, from the *in vitro* method design stage to the availability of a
 962 routine method, which establishes scientific evidence that the *in vitro* method is capable of
 963 consistently delivering quality and scientifically relevant data for the specific purpose it is
 964 designed for. Ideally the following conditions exist:

- 965
- Quality, safety, and efficacy are designed or built into the *in vitro* method.
- 966
- Quality and a sound scientific basis of the *in vitro* method are assured during the
- 967
- entire *in vitro* method life cycle towards regulatory acceptance.
- 968
- Each step towards a routine running of the *in vitro* method is controlled, and
- 969
- documented to assure that the *in vitro* method meets all scientific and quality
- 970
- attributes.

971 The validation process involves a series of activities taking place retrospectively and/or
972 prospectively over the lifecycle of the *in vitro* method once submitted to a validation body.

973 Taking into account the above conditions, test developers can submit their *in vitro* methods to
974 a validation body and ask for support in the validation process.

975 For example, the entire EURL ECVAM test submission process⁶ follows 2 mandatory steps:

976 *Step 1: Pre-submission - is based on the electronic submission of the EURL ECVAM*
977 *Test Pre-submission Form (TPF). The TPF will allow EURL ECVAM to perform a*
978 *preliminary assessment of the status of development of the in vitro method and the*
979 *potential relevance to the 3Rs (replacement, reduction, refinement of animal testing).*

980 *Step 2: Complete Submission - requires the compilation of a detailed Test Submission*
981 *Template (TST). The TST will be made available by EURL ECVAM after successful*
982 *conclusion of step 1 and the test submitter will be formally invited to complete it and*
983 *providing also a SOP in the ECVAM DataBase service on ALternative Methods (DB-*
984 *ALM) format.*

985 On the basis of a test submission, a validation body can make a final decision as to whether
986 the submitted *in vitro* method qualifies for entering the validation process.

987 In "Practical aspects of designing and conducting validation studies involving multi-study
988 trials" (Coecke *et al.* 2016), details are given that focus on practical aspects of conducting
989 prospective *in vitro* validation studies by laboratories that are EU-NETVAL laboratories.
990 Prospective validation studies within EU-NETVAL comprise multi-study trials involving
991 several laboratories or "test facilities" and typically consist of two main steps:

992 (1) The design of the validation study by EURL ECVAM and

993 (2) The execution of the multi-study trial by a number of qualified laboratories within
994 EU-NETVAL coordinated and supported by EURL ECVAM.

995 The approach adopted in the conduct of these validation studies adheres to the principles
996 described in the OECD guidance document on the Validation and International Acceptance
997 of new or updated *in vitro* methods for Hazard Assessment No. 34 (OECD, 2005a). The
998 (Coecke *et al.*, 2016) paper, mainly focuses on the processes followed to carry out a
999 prospective validation of *in vitro* methods involving different laboratories with the ultimate
1000 aim of generating a dataset that can support the development of an international test guideline
1001 (e.g. by the OECD) or the establishment of performance standards of *in vitro* methods.

1002 Upon successful validation by validation bodies or other entities, an *in vitro* method can be
1003 presented to the OECD for regulatory acceptance. Once *in vitro* methods are consolidated
1004 within an OECD test guideline (TG), data produced by using those methods are mutually
1005 accepted by all receiving authorities.

1006

⁶ <https://eurl-ecvam.jrc.ec.europa.eu/test-submission>

1007 **1.4 Producers of equipment, materials and reagents**

1008 When performing established testing methods (validated or not), the test results can only be
1009 accepted if the equipment, materials and reagents (test and reference items, media additives,
1010 compounds added to a system to induce a chemical reaction, etc.) used, are of proven quality
1011 as established by formal testing or evaluation procedures. To be able to prove this, test end
1012 users will need to work with preferred suppliers who are selected on predefined criteria (e.g.
1013 controlled transport, technical support, assured delivery, batch selection allowed, ISO
1014 certification present, etc.). The characteristics of the supplied materials and reagents should
1015 be appropriately documented in adequate quality documents such as a certificate of analysis,
1016 batch release certificate or similar. Suppliers who cannot fulfil these criteria have to be
1017 excluded from delivering products used to run *in vitro* methods, as the quality of test data
1018 cannot be guaranteed nor be acceptable for receiving bodies. It is the responsibility of the test
1019 end-user to verify and assure the quality of the products and materials (see chapter 2.4).

1020 **1.5 *In vitro* method end-users**

1021 As the aim of this document is to provide guidance for *in vitro* methods for regulatory use in
1022 human safety assessment, the majority of end-users will be GLP compliant test facilities. If a
1023 safety study is conducted for regulatory purposes, compliance with the principles of GLP is
1024 required. GLP test facilities are covered by national and/or international GLP regulations and
1025 must adhere to the responsibilities as defined in these regulations.

1026 The OECD Mutual Acceptance of Data (MAD) is a multilateral agreement which states that
1027 test data generated in OECD countries and full adherent countries – (Argentina⁷, Brazil,
1028 India, Malaysia, South Africa and Singapore)⁸ in accordance with OECD Test Guidelines and
1029 the OECD Principles of GLP shall be accepted in other member countries by regulatory
1030 bodies for assessment purposes and other uses relating to the protection of human health and
1031 the environment⁹. However, end-users should also document their competency to perform a
1032 test in compliance with a specific TG, by running the proficiency chemicals and
1033 demonstrating the overall quality assurance system of the test facility (see Section 6.5/8.6).

1034 Non-GLP *in vitro* method users can profit from the use of the GIVIMP guidance. In these
1035 cases no regulations exist and no responsibilities are defined. It is highly recommended to
1036 apply all necessary good scientific, technical and quality practices that the guidance describes
1037 so as to reduce the uncertainties in the use of cell and tissue-based *in vitro* method: Examples
1038 of critical importance described in this document relate to guidance on test systems such as
1039 the documentation concerning the origins of cells and tissues (see Table 1), the evaluation of
1040 the performance of the *in vitro* method, adequate measures to ensure test item exposure and
1041 test item-test system compatibility, the guidance for suppliers of equipment, materials and
1042 reagents, evaluation of competence and training for executing a particular *in vitro* method,
1043 etc. Appropriate accreditation (ISO, 2015, 2005) may be requested or recommended in some
1044 other cases.

1045 In case no specific Test Guideline harmonised template is available, all generated test data
1046 should be submitted in an easily readable format to facilitate the decision making process of

⁷ Full adherence for Argentina only applies to industrial chemicals, pesticides and biocides

⁸ <http://www.oecd.org/env/ehs/non-memberadherentstotheoecdsystemformutualacceptanceofchemicalsafetydata.htm>

⁹ <http://www.oecd.org/env/ehs/mutualacceptanceofdatamad.htm>

Field Code Changed

Field Code Changed

1047 the risk assessors, preferentially according to the OECD guidance document for describing
1048 non-guideline *in vitro* methods Series on Testing and Assessment No. 211 (OECD, 2014).

1049 Data generated by non-validated *in vitro* methods, if submitted as supporting data, need to be
1050 accompanied by proof of relevance and scientific validity of the selected test system together
1051 with the description of the critical points of the test system used, including sensitivity, limits
1052 of detection etc.

1053 **1.6 Receiving authorities**

1054 Receiving authorities receive non-clinical safety data as part of regulatory submissions and
1055 they must ensure that the legal requirements are met. Receiving authorities in Europe include
1056 the European Chemicals Agency (ECHA), European Medicines Agency (EMA), European
1057 Food Safety Authority (EFSA), as well as various national agencies that are responsible for
1058 assessing safety data. Receiving authorities in the United States include the Environmental
1059 Protection Agency (EPA) and the Food and Drug Administration (FDA).

1060 The responsibility of the receiving authorities is to check that test data are obtained according
1061 to available OECD TG and guidance documents and they use the data accordingly in their
1062 evaluations and provisions in law. For GLP studies, they may verify whether the responsible
1063 test facility has been found in compliance by a national monitoring authority (see below) or
1064 request a study audit if the facility has not been subject to a GLP inspection by a national
1065 monitoring authority. The receiving authorities may additionally indicate to *in vitro* method
1066 developers where they see needs for new or better tests, and to validation bodies which tests
1067 deserve priority in the validation.

1068 The majority of European regulation requiring toxicological data, allow or even encourage
1069 the use of alternative methods (e.g. Regulation (EC) No 1107/2009 for plant protection
1070 products, Regulation (EU) No 528/2012 for biocides, Regulation (EC) No 1907/2006 for
1071 industrial chemicals, known as REACH) (Heringa *et al.*, 2014). Regulation (EC) No
1072 1223/2009 for cosmetics does not allow any *in vivo* animal testing.

1073 As a result of these developments European, but also other international regulatory bodies,
1074 tend to accept data generated by alternative methods, including validated *in vitro* methods
1075 and, in specific cases, also non-validated *in vitro* methods, especially as supportive
1076 information or when mechanistic data are needed. They have adopted the 3Rs principles and
1077 are now proactively supporting the use and implementation of alternative methods¹⁰.

1078 The European Medicines Agency (EMA) has recently drafted a number of documents
1079 expressing their vision and action plans towards the implementation of the 3R principles
1080 (EMA, 2014). According to the new EMA draft guideline on regulatory acceptance of 3Rs
1081 methods, multiple and flexible approaches are considered acceptable to demonstrate scientific
1082 validity of new testing approaches and their fitness for regulatory use, either as pivotal,
1083 supportive or as exploratory mechanistic studies. Besides established formal validation
1084 processes by recognised institutions such as the Centres for the Validation of Alternative
1085 Methods (CVAMs) and The European Directorate for the Quality of Medicines &
1086 HealthCare, the EMA draft guideline is offering a new approach for submission and
1087 evaluation of a proposal for regulatory acceptance of 3R testing approaches via an EMA in-
1088 house procedure. One interesting option of this process is the collection of real-life data by
1089 using a new 3R testing approach in parallel with the data generated using existing methods.

¹⁰ <https://echa.europa.eu/support/registration/how-to-avoid-unnecessary-testing-on-animals/in-vitro-methods>

Field Code Changed

1090 Data generated with the new method will however not be used as part of the regulatory
1091 decision making process (“safe harbour”) and will be used solely for the purpose of
1092 evaluation of the new method for possible future regulatory acceptance.

1093 **1.7 GLP monitoring authorities**

1094 For studies conducted for regulatory purposes, the responsibility for evaluating the results of
1095 the study lies with the regulatory reviewer at the receiving authority. However, this
1096 evaluation can only be effective if the study data can be relied upon. The principles of Good
1097 Laboratory Practice (GLP) ensure that the quality and scientific integrity of the data can be
1098 demonstrated and the conduct of the study reconstructed.

1099 *GLP was developed in the 1970s in response to fraudulent scientific safety studies that were submitted*
1100 *to regulatory authorities in support of applications for the regulatory registration/approval of drugs to*
1101 *FDA. Subsequently the principles of GLP were developed by the OECD to ensure data reliability and*
1102 *reconstructability of safety studies. The principles apply to all non-clinical health and environmental*
1103 *safety studies required by regulations for the purpose of registering or licensing chemical products of*
1104 *various kinds. The principles have been published in 1981 as an annex to the OECD Council Decision*
1105 *on MAD¹¹. The decision states that ‘data generated in the testing of chemicals in an OECD Member*
1106 *Country in accordance with the OECD Test Guidelines (Annex I of this decision) and OECD*
1107 *Principles of Good Laboratory Practice (Annex II of this decision) shall be accepted in other member*
1108 *countries for purposes of assessment and other uses relating to the protection of man and the*
1109 *environment’. Since 1981 a number of additional guidance, consensus and advisory documents have*
1110 *been published in the OECD Series on Principles of GLP¹².*

1111 *The Advisory Document of the OECD Working Group on GLP n° 14 specifically addresses in vitro*
1112 *Studies (OECD 2004b): ‘The purpose of this document is to facilitate the proper application and*
1113 *interpretation of the GLP Principles for the organisation and management of in vitro studies, and to*
1114 *provide guidance for the appropriate application of the GLP Principles to in vitro studies, both for test*
1115 *facilities (management, QA, archivist, study director and personnel), and for national GLP compliance*
1116 *monitoring authorities.’*

1117 *In the European Union, the principles of GLP are included in Directive 2004/10/EC, while the*
1118 *compliance monitoring procedures are included in Directive 2004/9/EC. GLP provisions are included*
1119 *in legislation for chemicals, human medicinal products, veterinary products, detergents, feed additives,*
1120 *food additives, genetically modified food or feed, pesticides, biocides and cosmetics (Coecke et al.*
1121 *2016). The current European medical device directives do not require GLP, but the harmonised*
1122 *standard series ISO 10993 require that all in vitro tests ‘shall be conducted according to recognised*
1123 *current/valid best laboratory/quality practices, for example Good Laboratory Practice (GLP) or*
1124 *ISO/IEC 17025’. However, the current proposal for a new Medical Device Regulation does require*
1125 *GLP where applicable.*

1126 The role of GLP Monitoring Authorities (MAs) may vary considerably from region to region,
1127 however their responsibilities in general remains the same, i.e. to ensure the compliance of
1128 GLP studies. GLP Compliance Monitoring is intended to ascertain whether test facilities have
1129 implemented the GLP Principles and that studies are performed in compliance with the GLP
1130 principles.

1131 The OECD expects Member countries to establish national MAs, a body or bodies
1132 responsible for monitoring the GLP compliance of test facilities within its territories and
1133 according to national legal and administrative practices. In the European Union, the
1134 principles of GLP are included in Directive 2004/10/EC, while the compliance monitoring
1135 procedures are included in Directive 2004/9/EC, where each member state is required to
1136 implement a monitoring programme, where GLP registered facilities are inspected on a

¹¹ <http://acts.oecd.org/Instruments/ShowInstrumentView.aspx?InstrumentID=263&InstrumentPID=361>

¹² <http://www.oecd.org/chemicalsafety/testing/goodlaboratorypracticeglp.htm>

Field Code Changed

Field Code Changed

1137 regular basis, approximately every two to three years. Routine monitoring inspections will
1138 include study audits. In addition, MAs can be requested by a receiving authority to conduct
1139 specific study audits as a result of concerns raised following the review of a regulatory
1140 submission. The MA has ultimate responsibility for determining the GLP compliance status
1141 of test facilities and the acceptability of a study audit. The MA also has responsibility for
1142 taking any action based on the results of test facility inspections or study audits which are
1143 deemed necessary.

1144 The respective national compliance MAs are also responsible for the exchange of information
1145 on the compliance of test facilities inspected, and also should provide relevant information
1146 concerning the countries' procedures for monitoring compliance. They have the
1147 responsibility to facilitate the mutual acceptance of test data (MAD) generated for submission
1148 to regulatory authorities of OECD Member countries and other countries that are full
1149 adherents to MAD.

1150 **1.8 Accreditation bodies**

1151 The International Organisation for Standardisation (ISO) is an independent, non-
1152 governmental membership organisation and the world's largest developer of International
1153 Standards with a central secretariat based in Geneva, Switzerland. The ISO story dates back
1154 to 1946 when delegates from 25 countries met at the Institute of Civil Engineers in London
1155 and decided to create a new international organization 'to facilitate the international
1156 coordination and unification of industrial standards'¹³.

1157 In this organisation, different industries define their specific technical standards and quality
1158 management requirements and issue ISO standards to guide conformity. Companies and
1159 organisations working according to ISO guidelines can ask for a conformity check and
1160 certification by independent accreditation bodies. ISO itself is not a controlling body, but has
1161 established a committee on conformity assessment (CASCO) guiding certification
1162 organisations.

1163 While the OECD Principles of GLP and ISO/IEC 17025 (OECD, 2016a) both set out
1164 requirements for quality management systems under which testing is conducted, they have, as
1165 a result of their evolution and history, different purposes (OECD, 2016a).

1166 The OECD Principles of GLP are used as a regulatory control mechanism to assure the
1167 quality and integrity of non-clinical health and environmental safety studies regulated under
1168 law. Such testing, for the most part, is complex and variable. The Principles are specifically
1169 designed to be applied to individual studies to accommodate the complexity and wide variety
1170 of such studies due the different scientific disciplines involved, for the different chemicals
1171 under test. The OECD Principles of GLP are therefore, out of necessity, quite general in their
1172 requirements and take the form of a set of principles.

1173 ISO/IEC 17025 is an international standard intended to be applied by laboratory facilities
1174 conducting testing, according to established or specifically developed methodologies. The
1175 focus of the standard is on the on-going operation, monitoring and management of the
1176 laboratory itself, and on the capacity of the laboratory to produce consistent and reliable
1177 results that are scientifically valid. ISO/IEC 17025 can, in theory, be applied to any testing
1178 laboratory in any scientific discipline including those performing non-clinical testing. It is a

¹³ <http://www.iso.org>

1179 reliable indicator of technical competence, and many industries routinely specify laboratory
1180 accreditation for suppliers of testing services.

1181 **1.9 OECD**

1182 The OECD strives for international harmonisation of test methods for toxicity and risk
1183 assessment of new products, and issues globally accepted Test Guidelines (TGs) that are
1184 accepted by regulatory bodies of all OECD member states and MAD-adherent economies.

1185 The availability of a guidance document on Good *In vitro* Method Practices (GIVIMP) for
1186 "The development and implementation of *in vitro* methods for regulatory use in human safety
1187 assessment" was identified by OECD as a high priority requirement by both the Working
1188 Group on GLP (WG GLP), and nominated experts from the Working Group of the National
1189 Coordinators of the Test Guidelines Programme (WNT) to reduce the uncertainties in cell
1190 and tissue-based *in vitro* method derived predictions with the OECD *in vitro* test guideline
1191 methods.

1192 Other national and international organisations are developing policies and standards for
1193 scientific practice to assure quality in implementation of *in vitro* methods. OECD has by
1194 supporting the GIVIMP endeavour, taken the responsibility to seek consensus on the current
1195 best *in vitro* method practices and will help to ensure that *in vitro* methods used for
1196 regulatory purposes are scientifically sound, of high quality, reproducible, credible, and
1197 acceptable. Generic *in vitro* method OECD TGs incorporating performance standards are
1198 being written to allow acceptance of both non-proprietary and proprietary *in vitro* methods by
1199 regulatory agencies and to provide assurance that any *in vitro* cell and tissue culture system
1200 performs over time in a manner that is consistent with the test system as it was originally
1201 validated (Gupta *et al.* 2005; Rispin *et al.* 2006).

1202

1203 2 Quality Considerations

1204 2.1 Quality assurance versus quality control

1205 In order to adopt GIVIMP during the test development phase or the routine phase of an *in*
1206 *vitro* method, it is necessary to have a good understanding of the relevant QA frameworks
1207 and definitions. In particular, users should be aware that a quality assurance programme in an
1208 OECD guidance has a particular meaning (OECD, 2005a, 2004c, 1998a), that is a defined
1209 system, including personnel, which is independent of study conduct and is designed to assure
1210 test facility management of compliance with the Principles of GLP.

1211 There are numerous definitions of quality control, but for the purposes of GIVIMP, quality
1212 control is the documented activity which seeks to confirm that starting materials for assays
1213 (including cell lines) and key stages of individual assays and the final results meet prescribed
1214 specifications. It should be clear that this is not the same as Quality Assurance, which is the
1215 overall quality system designed to assure the quality of results (see Table 2). Specific
1216 requirements may be reviewed with respect to the criticality of any change, but it must be
1217 borne in mind that a good quality system should be under ongoing review to ensure current
1218 best practice is sustained and to enable the requirement for continuous improvement.

1219
1220

Table 2: Examples of differences between quality assurance and quality control

Feature	Quality Assurance (QA)	Quality Control (QC)
Nature	Proactive: prevents occurrence of errors, is process oriented and is a managerial tool.	Reactive: detects an error, is product oriented and is a corrective tool.
Format	Set of SOPs to assure that the required standard is met at every stage in the process.	Review and test quality in a product or service, against a set of given requirements.
Aims	Doing things right the first and every time to prevent errors.	Detect deviations/defects that need corrective action.
Scope	Focused on continuous improvement of processes.	Typically follows the process established as part of QA function.
Responsibilities	Management but also requires commitment from all staff.	May be performed independently but overall responsibility lies with organisation delivering the product or service.
Examples of relevant activities	Creating a process. Designing templates/check lists. Identifying defects in a process (i.e. QC). Root cause analysis. Making recommendations to improve a process or service. Developing a measurement system to assess process effectiveness. Checking a process.	May include: Reviews. “Walk-throughs”. Inspection. Testing (verification and validation). Monitoring control charts.

1222 2.2 Quality Control of test system

1223 Quality control of the *in vitro* cell and tissue test systems must be a shared responsibility of
1224 the manufacturer of for instance proprietary test systems or suppliers of cells, tissues or tissue
1225 constructs, the test facility which uses the *in vitro* test system, and the entity that submits the

1226 *in vitro* method and the related test system for regulatory acceptance. Tissue constructs or cell
1227 cultures may have a short shelf life. Proprietary *in vitro* methods and the related *in vitro*
1228 systems may be relatively expensive; therefore the number of replicate systems available for
1229 quality control efforts by *in vitro* test facilities may be limited by practical considerations. In
1230 light of these considerations, the user may sometimes be dependent on the manufacturer for
1231 many of the basic elements of quality control, including cell or tissue characterisation and
1232 functional performance of the *in vitro* test system. The manufacturer should be expected to
1233 provide adequate documentation of quality control testing of representative test system from
1234 each batch manufactured.

1235 In addition, the user must provide a quality control check in the test facility on a regular basis
1236 appropriate to the test system so that the *in vitro* method performs as expected after transport
1237 and handling of the test systems in that particular method. For cell cultures, records should be
1238 kept of cell passage, of preservation, and of quality controls.

1239 Moreover, records recommended by GCCP or other relevant guidance documents (e.g. ISO
1240 standards, GLP) should be kept. The developer/user should confirm the key quality attributes
1241 (essential characteristics) on receipt, and periodically confirm the genetic/phenotypic
1242 stability, identity and freedom from contamination. Guidance on cell and tissue culture work
1243 is available for either general (Coecke *et al.*, 2005) or specific (Andrews *et al.*, 2015;
1244 Geraghty *et al.*, 2014; ISCBI, 2009) applications.

1245 **2.3 Quality control of consumables and reagents**

1246 Consumables (plastic ware such as flasks, cryovials, culture dishes, culture slides, tubes, cell
1247 scrapers) and *in vitro* method reagents (test and reference items, media additives, compounds
1248 added to a system to induce a chemical reaction, etc.) may have passed through a number of
1249 entities (producers, suppliers, shippers etc.) thus, their quality must be assured by good
1250 documentation, traceability, confirmation of suitable storage/shipment, identification, quality
1251 control and preparation for use in an *in vitro* method. Attention will also need to be given to
1252 the suitability of reagents and to the safety and ethical provenance of cells (Annex 2 and
1253 Annex 3).

1254 Consumables and reagents should be reliably available, consistent in their properties and for
1255 critical reagents, alternative sources should be identified. Certain materials which are critical
1256 to the performance of a method and which may be subject to significant variation, such as
1257 conventional serum products and growth promoting reagents, should be quality controlled
1258 before use (by functional tests) and this may include in-house pre-use testing. In addition,
1259 SOPs should mention that expiry dates should be documented.

1260 Laboratories can perform quality control checks of consumables, but the process how to do
1261 this is not always evident. Some laboratories have established procedures whereby a
1262 percentage of consumables from each lot number are evaluated prior to use in *in vitro* work.
1263 This may be especially useful for test facilities inquiring the cause of contamination. While
1264 this approach will not prevent contamination, it can provide data from any profile(s)
1265 developed during these checks, which can be used for future evaluation of potential
1266 contamination events. Ideally, sterile consumables with appropriate certificates should be
1267 used where possible. Alternatively, some consumables can be treated with ultraviolet (UV)
1268 light, gamma irradiation and/or autoclaved. These preventive measures may be useful in
1269 limiting contamination events. Other consumables, such as centrifugal filter units and filtered
1270 pipette tips, cannot be pre-treated. In the case no commercial sterile and documented
1271 centrifugal filter units and filtered pipette tips are available, establishing a method for
1272 detecting contamination from these items is very important.

1273 Again, as far as possible, reagents that are documented as sterile should be used without
1274 further quality control. Alternatively, laboratories should also run quality control checks on
1275 reagents prior to use in *in vitro* work. These checks assist in determining if a reagent is free of
1276 contamination at that time. Negative controls and reagent blanks provide a means to detect
1277 contamination from reagents. Negative controls can then be assessed on an ongoing basis to
1278 demonstrate that they remain contaminant free. Including and assessing negative controls and
1279 reagents blanks are critical quality control steps. These controls provide a means of detecting
1280 reagent contamination and, on occasion, sporadic contamination. Because many
1281 contamination events are sporadic, negative results in these controls do not necessarily mean
1282 that samples from the same batch are contaminant free. Additionally, the detection of
1283 contamination in these controls does not mean that all batch samples have been affected.

1284 Reagents may be selected on the basis of history of use or reference in relevant documents
1285 associated with regulatory accepted *in vitro* methods (e.g. validation reports, *in vitro* method
1286 SOPs). Established reagents should have predefined acceptance criteria, which must be met
1287 and accepted before a new batch is used. Same reagent obtained from different suppliers may
1288 each have specific and not necessarily similar acceptance criteria. Acceptance criteria must
1289 be established for different reagents considering the degree of risk they represent to the final
1290 results. This risk can be assessed based on the manufacturer/supplier's information on
1291 traceability and control on sourcing, processing, quality control (QC) and shipment. This can
1292 be addressed by:-

- 1293 1) Considering the potential impact of the perceived risk to prioritise certain reagents,
- 1294 2) Formally evaluate (the quality management system of) the supplier and
- 1295 3) Establishing suitable service level agreements (SLA) with the provider ensuring
1296 quality, availability and shipment of the reagent. Acceptance of individual batches of
1297 reagents can be addressed by review of key elements of the certificate of analysis,
1298 compliance with specific conditions of the SLA provided by the manufacturer/
1299 supplier or a combination of these and supplementary evaluation which may include
1300 pre-use testing to assure that individual batches are fit for purpose.
- 1301 4) Consistency of lot qualification tests on critical reagents. Serum can be a critical
1302 reagent for cell-based *in vitro* methods and lot qualification testing is common
1303 practice. Standardised serum is available. However, functional tests including
1304 acceptance criteria needs to be defined.

1305 Any reagent stocks prepared in the test laboratory for use in routine *in vitro* methods should
1306 be documented, assuring use of acceptable reagents (as above) and documenting the
1307 preparation of these stocks in batch preparation records. A batch number, storage conditions
1308 and an expiry date should be assigned (and printed on the label), which can be documented in
1309 records of performance of routine assays. It is also wise to have procedures for maintaining
1310 and controlling laboratory stocks of reagents such as minimum stock levels and identification
1311 and qualification of alternate suppliers to support continuity of testing studies etc. The
1312 stability of test items and reference items under storage and test conditions should be verified,
1313 if appropriate, and expiry dates allocated as appropriate (see GLP principles II.6.2.4). If a test
1314 item is administered in a vehicle, the homogeneity, concentration and stability of the test item
1315 in that vehicle should be determined (see GLP principles II.6.2.5).

1316 **2.4 Data management**

1317 Before beginning to collect data from *in vitro* test procedures, it is important to assess the
1318 format of collection, the complexity involved and requirements for traceability, storage,

1319 verification and transmission of data. Data should be recorded concurrently with the
1320 performance of the procedures. Specific standards may apply for data from regulatory testing
1321 and manufacturing (Coecke *et al.*, 2005; FDA, 2003; OECD, 1999). Data from material
1322 provided by tissue donors may also be subject to the requirements of data management and
1323 control under local, regional, national or international rules and regulations such as the EU
1324 Directive on Data Protection (NB. national and regional rules should be consulted as these
1325 may vary). It should be ensured that data reported accurately reflect the results obtained
1326 during experimental work, by performing adequate quality control of the data.


1327 **2.5 Types of documentation**

1328 Documentation in a quality system typically involves documents and records at several levels
1329 (Figure 1). The overarching document is a high-level, accurate description of the types of
1330 work performed by the organisation or group, key policies and standards adopted for
1331 delivering the work and the structure of the quality system used to deliver all aspects of test
1332 performance and provision of results. In some systems, this may be called a “quality
1333 manual”. Another level may include overviews of procedures referring to the various specific
1334 testing methods involved at the next level. Finally, supporting the SOPs, there will be formal
1335 record sheets for test and control data and templates for reporting results. Another example
1336 for descriptive and prescriptive documents is provided in the WHO handbook on quality
1337 practices in Biomedical Research (WHO, 2013).

1338 In many formal quality systems, including GLP, there will need to be a document control
1339 convention which assures that all documents are developed and approved in a formal process,
1340 that active versions are accurately dated, authored and approved with specific version
1341 numbers that will avoid inadvertent use of obsolete testing methods.

1342

1343 **Figure 1: Hierarchy of governance and management documents**



1344

1345 NB: this schematic overview is intended to reflect some general aspects of documentation
1346 hierarchy and is not a prescriptive requirement of GIVIMP or GLP. Another example for
1347 descriptive and prescriptive documents is provided in the WHO handbook on quality
1348 practices in Biomedical Research (WHO, 2013). QA documentation may be structured in
1349 different ways as required in different organisations and regulatory jurisdictions.

1350 **2.6 Staff training and development**

1351 Training files are helpful to capture all procedures in which staff should be trained before
1352 they are considered competent. These may include regulatory requirements of the work,
1353 specific controlled procedures (e.g., SOPs, testing methods), use and storage of
1354 documentation, as well as general training in best practice such as indicated in GCCP (2005).
1355 When new staff is recruited to work in the laboratory, it is important to guide the staff and
1356 review and document any training requirements before assignment to carry out any tasks. A
1357 written training plan might serve as a basis for a formal training record, which should be
1358 reviewed periodically by their line manager. It may be helpful to demonstrate competence in
1359 each procedure by recording individual elements of training in a particular SOP, until the
1360 trainee is considered competent and the training records are signed to authorise that they can
1361 perform the procedure independently. Training may include:

1362 1) self-study of testing methods and SOPs and recording that the document has been
1363 read and understood,

1364 2) correct performance of the procedure, witnessed by a qualified supervisor a
 1365 minimum number of times and notification from the trainer or supervisor that the staff
 1366 member is competent to perform the task independently.

1367 A list of core training for *in vitro* cell culture laboratory staff is given in GCCP (Coecke *et*
 1368 *al.*, 2005) and special aspects of training are also referred to in other sections of this
 1369 document where relevant.

1370 It is good practice to record all training in staff training records or periodic competency
 1371 reviews, which should be regularly updated. In some circumstances, the host organisation or
 1372 relevant professional bodies may run or support attendance at supplementary training and
 1373 education which should also be documented to demonstrate maintenance of ongoing
 1374 professional development and support assurance that current best practice is maintained in
 1375 testing procedures. Annual review of staff performance is also a useful tool for considering
 1376 ongoing training needs.

1377 **2.7 Assurance of data integrity**

1378 Data integrity is fundamental for any quality system and starts with good quality
 1379 documentation (see Table 3).

1380 **Table 3: Key criteria to be addressed for data generated under a GLP environment (ISPE,**
 1381 **2012)**

Criteria		GLP
Attributable	Who acquired the data or performed the action	Initials/signature/UserID
Legible	Data must be recorded permanently in a durable medium and be readable	
Contemporaneous	Documented at the time of activity	Date and time (stamp)
Original	Original or true copy	
Accurate	No errors, no undocumented changes	QC, QA audits
Complete	All data available	QC, QA audits
Consistent	Traceability, dataflow, date timestamps	
Enduring	Recorded on paper or electronically	
Available	Accessible for the lifetime of the record	Archiving

1382
 1383 For validation of certain documents, such as spread sheets, there may be specific procedures
 1384 to adopt good practices e.g. (Esch *et al.*, 2010; OMCL, 2009). It may be necessary to have an
 1385 audit trail of modifications.

1386 If data is translated between different recording methods, systems and/or data bases and in
1387 particular critical phases like manual or semi-automatic transfer (e.g. Excel™ files to
1388 database for kinetic assessment, combination of information obtained from 2 or 3 databases
1389 to one database), correct resolution of pre- and post-translation data should be reviewed and
1390 confirmed by a qualified person. For handwritten data, translated into an excel sheet, it is also
1391 advisable for the changes to be verified by the same person who has made the observations.
1392 These issues are of special concern where data are exchanged between countries (e.g. for
1393 pharmacokinetic assessment). When data translation occurs between different software or
1394 data base systems their compatibility and inability to be altered in translation should be tested
1395 and will need to involve appropriate validation procedures (see i) below)

1396 Data Integrity (DI) closely correlates with good scientific practices. The same principles
1397 should be applied when using either paper based and electronic systems, or both. It should be
1398 assured that the obtained information is unchanged from the source, and has not been
1399 intentionally modified, altered or destroyed. To ensure DI for both systems, the following
1400 components of this process should be taken into consideration.

- 1401 a) Sample integrity: Information and decisions based on data from samples used must be
1402 robust, so sample control should be ensured from point of entry into the laboratory
1403 through to data review and data release.
- 1404 b) Integrity of materials/solutions: The right standard of materials/solutions must exist to
1405 ensure a valid reference, so the process for standardisation and requirements for re-
1406 standardisation should be established e.g. documentation of receipt, calculation and
1407 preparation, traceability, assignment of expiry or retest dates etc. Reference data
1408 needs to be complete (e.g. lot number, reference standard data). For further details see
1409 Section 4.2.
- 1410 c) Instrument and equipment: The data generated by instruments and equipment must be
1411 robust. The equipment unit(s) must be suitable for their “intended use”, therefore the
1412 equipment must be calibrated or qualified by appropriately trained personnel.
1413 Calibration should, where appropriate, be traceable to national or international
1414 standards of measurement. In order to operate within specification, scheduled
1415 maintenance of the equipment and instruments once again must be performed by
1416 qualified personnel. For examples and more information see Section 4.1.
- 1417 d) Documentation and result reporting: Records must be clear and accurate. All activities
1418 should be recorded at the time they are performed (see Annex 2). Records should also
1419 be chronological, traceable, and readily retrievable. Original documents must be
1420 clearly identifiable (e.g. time stamps, watermarks) and standardised, predefined,
1421 authorised forms and templates should be used wherever possible and applicable.
1422 Records should be signed and dated allowing for clear identification, no pencil, no
1423 recording on loose ‘post-it’ notes, “white-out” paste or scrap sheets of paper should
1424 occur. Any corrections written on documents should be signed and dated by a trained
1425 staff member. Transcriptions if performed need to be attached to the original results
1426 (full traceability) and reviewed. Chronology of recorded data and correct page
1427 numbering must be ensured.
- 1428 e) Proficiency of the analyst/operator: The person producing the data using samples,
1429 materials/solutions, instruments and equipment, documenting and reporting results
1430 must be relied upon to produce robust data. Therefore, the ability of the person to
1431 perform a task accurately and proficiently should be ensured, e.g. through appropriate
1432 training (see section 2.6). The person must also understand that destroying, deleting,
1433 hiding or using selected data without scientific justification (or even falsifying data)
1434 are not acceptable. Facility management policies should assure a working

- 1435 environment in which that when a mistake has been made, the individual concerned
1436 should feel comfortable to report and explain it and not conceal it.
- 1437 f) Method validation and verification: The *in vitro* method has to be validated or verified
1438 “suitable for its intended use”.
- 1439 g) Effective review and verification: A clear definition and understanding of raw data
1440 should be ensured. There needs to be traceability to the testing method used, source
1441 data and verification of raw data. SOPs need to be in place for data handling, record
1442 retention and good documentation practices and deviation handling etc.
- 1443 h) Additional considerations for electronic data: In addition to comments made above, if
1444 a system is required to maintain electronic data, it should be managed by unique user
1445 identity and password combination. If the system does not permit this, a paper-based
1446 log must be in place to record who uses the generic user and password combination,
1447 or who uses the unprotected equipment. Paper records can be reviewed for any
1448 amendments or crossings out/deletions plus the signature/date of and the reason for
1449 doing so. This is to be replicated in an electronic system in the same way by use of an
1450 electronic log (audit trail). Where there are multiple users, an approach to periodically
1451 review who has system access must be in place. There must be a periodic user
1452 account review procedure. There should be procedures in place for assigning access
1453 rights to each user. The level of access should be in line with the tasks that have to be
1454 performed.
- 1455 i) Data storage, archive & retrieval: Data must be stored in a safe and secure place for
1456 paper-based systems and in protected folders for electronic systems. An approach
1457 must be in place to ensure that data are protected against loss, damage or overwriting.
1458 Manipulation of stored paper or electronic records must not be possible. Electronic
1459 records must be held in a format that is not readily corruptible and protected from
1460 deliberate or accidental alteration (e.g. CFR 21 part 11, GLP: see OECD GLP
1461 Guidance Document 17).


1462 A system must be in place for the archive of electronic and paper data and records (see 10.3).
1463 It should be considered and tested how the data and records can be retrieved and used to re-
1464 create a complete accurate picture on the rare occasions that this is required.

1465

1466 **3 Facilities**

1467 A very important aspect in the process of *in vitro* method development to *in vitro* method
1468 routine use is the need for a well-designed facility to ensure that good quality results are
1469 produced in a safe and efficient manner. The types of laboratory areas which might need to
1470 be separated are indicated in Figure 2.

1471
1472 **Figure 2: A typical "onion ring" structure used to separate different operational areas.**



1496 It may not be possible or acceptable to separate such laboratory functions but they are among the
1497 considerations for separation (physical or process/training) that someone establishing or running the
1498 facility should be aware of and there may be other issues which must be considered for different kinds of
1499 work. It is wise to avoid physical contact between materials transfers (blue arrows) and waste removal
1500 (red arrows) so that there is very low risk of contamination from waste affecting reagents, cultures and
1501 test materials.

1502 As stated in OECD No 14 (OECD, 2004c), facilities should be suitable to maintain the assay,
1503 resulting data and all archived items.

1504 Facilities must be fit and suitable for the purpose of the work; that is, size, construction, and
1505 location should be appropriate, and the building should allow for the separation of activities.

1506 Both test developer facilities and routine use facilities should ensure to be designed or
1507 adapted to have separate areas for similar but unrelated work and sample preparation. In
1508 addition, buildings should be validated for the required functionality (e.g. air handling) and
1509 properly monitored (e.g. air pressure differences, flow patterns, etc.), with easily accessible
1510 results. Validation of facilities should be formally documented. Environmental control
1511 systems should be regularly maintained and serviced, with full records of maintenance and
1512 any modifications to demonstrate appropriate upkeep and function. All the necessary permits
1513 should be in place before any activities are initiated. Finally, there should be dedicated areas
1514 for data storage and archiving.

1515 **3.1 Containment**

1516 Cell lines and primary tissues may carry a variety of different microorganisms or pathogens,
1517 which can potentially cause human disease, pose hazard to employees and distort the *in vitro*
1518 method results. Cell lines and primary tissues should be handled at biosafety (hazard) level 2,
1519 unless the cells are known to be specific pathogen-free. This level of containment is also
1520 appropriate for monoclonal antibody-containing supernatants and cell homogenates. Access
1521 to level 2 facilities should be restricted to authorised personnel only, and specific risk
1522 assessment and training activities should be followed according to the national legislation on
1523 Level 2 containment (Coecke *et al.*, 2005; Geraghty *et al.*, 2014).

1524 As a minimum, all cell and tissue work should be performed in a Class II biological safety
1525 cabinet as even screened tissues or cell cultures may carry infectious agents not covered by
1526 virological screening. Depending on the kind of test items in the *in vitro* method, several
1527 subtypes of Class II biological safety cabinets must be considered. Horizontal (Class I) flow
1528 cabinets, where the airflow is directed at the operator, are intended for sterile media
1529 preparation and not for use with tissues and cell cultures due to the inability to exclude all
1530 infectious agents in screening and to avoid the widespread laboratory contamination with
1531 mycoplasma or bacterial and fungal contaminants. For this reason they should only be used if
1532 the cells are known to be pathogen-free. Class I cabinets are also not appropriate when
1533 potentially toxic chemicals might be used (which is often the case for *in vitro* methods for
1534 regulatory use). In the case of using toxic chemicals it is recommended to let the outflowing
1535 air pass through a volatile organic compound-filter (active charcoal) after a High-Efficiency
1536 Particulate Arrestance (HEPA) filter.

1537 If microscopes or other equipment are to be installed in biological safety cabinets, the cell
1538 containment equipment should be checked for flow disruption as well as operator and cell
1539 culture protection.

1540 Splashes and aerosols carry contamination and infection risks, which can not only endanger
1541 the operator, but also compromise the integrity of the *in vitro* method (i.e., cross
1542 contamination of cell lines or introduction of adventitious agents). Therefore, all procedures
1543 should aim at minimising aerosol production. Any procedures likely to produce aerosols
1544 should be contained (biosafety cabinet) or the material should be rendered harmless.

1545 As contaminated working surfaces can lead to microbial contamination or cross-
1546 contamination between cell lines and pose a risk to the *in vitro* method quality, working
1547 surfaces should be easy to clean, resistant to acids, alkalis, solvents and disinfectants. There
1548 should be appropriate documented procedures for disinfection of work surfaces, safety
1549 cabinets and equipment.

1550 Higher containment levels may be required depending on the biosafety risk level of the
1551 biological agents handled. If *in vitro* work is to be performed with group 3 or 4 human
1552 hazards, which can cause severe human disease and may be a serious hazard to employees or
1553 spread to the community, then separated facilities, appropriate levels of biosafety, such as air
1554 filtration and negative pressure differences will need to be maintained. Groups 3 and 4 are
1555 more complex in complying with specific facility requirements and personnel skills.
1556 Therefore, *in vitro* methods for regulatory use in human safety assessment should be
1557 developed to require mainly level 2 biosafety.

1558 **3.2 Level of separation to avoid cross-contamination**

1559 Measures should be taken to ensure adequate separation of different biological agents and *in*
1560 *vitro* studies taking place in the same physical environment (OECD, 2004a). The integrity of
1561 each test system and study should be protected by spatial or temporal separation from other
1562 studies and test systems to avoid cross-contamination and mix-up. Air flow requirements, for
1563 example, are opposite for spaces where human cell lines are used (flow out of the space to
1564 avoid their contamination) and for spaces where pathogenic micro-organisms are used for
1565 tests (flow into the space to avoid spread of the pathogens).

1566 Tissues and cells from different studies can be kept in the same incubator provided that they
1567 are labelled appropriately, placed in a different position and none of the test items are volatile
1568 enough to cause contamination. Tissues and cells from different species or *in vitro* methods
1569 where yeast and bacteria are used would require a higher level of separation. The most
1570 important issue here is to separate the areas used for cell culture/tissue and microbiological
1571 culture (see Figure 2). Other degrees of separation may be achieved using the specific
1572 requirements described elsewhere for quarantine of untested material.

1573 Temporal separation of test systems should be used in Class II biological safety cabinets.
1574 This can be achieved by handling only one cell line at a time, followed by cleaning and
1575 decontamination of the cabinet, working surfaces and related equipment.

1576 Rooms and areas used for preparation and mixing of test and reference items with vehicles
1577 should allow aseptic working conditions in order to minimise the risk of microbial
1578 contamination of the test system.

1579 Appropriate training should be given to the authorised personnel regarding the necessary
1580 precautions to prevent contamination and cross-contamination.

1581 When performing molecular biology techniques and especially PCR-based assays, which are
1582 high sensitivity methods, extreme care should be taken in facility design and operation in
1583 order to avoid false-positive results. False-positive results can originate, for example, from
1584 sample-to-sample contamination, from carry-over of nucleic acid from previous amplification
1585 of the same or similar target. Cloned DNA or virus infected cell cultures may represent other
1586 source of contamination¹⁴.

1587 The greatest threat of contamination lies in laboratories that manipulate amplified or cloned
1588 DNA; laboratories exclusively performing real-time PCR and properly discarding all
1589 amplified products without opening the reaction tubes or sealed plates are less liable to
1590 contamination.

¹⁴ <https://www.gov.uk/government/publications/smi-q-4-good-laboratory-practice-when-performing-molecular-amplification-assays>

1591 It is vital that the correct workflow is followed in a molecular laboratory in order to minimise
1592 contamination and ensure good laboratory practises are followed. It is the responsibility of all
1593 laboratory staff to ensure that the workflow is followed. The equipment, consumables and
1594 laboratory coats should be dedicated to each area. If possible it is helpful to colour code
1595 racks, pipettes and laboratory coats in the different areas to be able to easily monitor
1596 movement between the different areas. Powder-free gloves should be used throughout the
1597 process in all the different areas as the powder on powdered gloves might result in assay
1598 inhibition. It is particularly important to always use powder-free gloves in the pre-PCR area,
1599 as this is prone to contamination by RNases.

1600 Perform physical separation of pre- and post-PCR assay stages. Between these two areas the
1601 work flow should be uni-directional and the relative air pressure and direction should differ.
1602 It is recommended for personnel working with post-PCR assay stages to not work with pre-
1603 PCR parts later the same day.

1604 PCR reactions should be set up in a separate room from that used for post-PCR manipulation.
1605 Bench areas should be wiped daily with hypochlorite solution following use and contaminated
1606 areas should be additionally decontaminated with ultra-violet radiation if fitted. Hypochlorite
1607 solution (20% or greater) should not be applied to stainless steel (types 304/347, 316 and 400
1608 series) as it leads to corrosion with repeated use. Reagents should be taken from clean storage
1609 into the pre-PCR area and never taken or shared with post-PCR areas.

1610 Where possible, the PCR facility should be organised in four discrete rooms/areas:
1611 (Requirements may vary with assay format e.g. real time PCR does not require post-PCR
1612 analysis).

- 1613 - **Reagent preparation clean room:** it is free from any biological material such as
1614 DNA/RNA, cloned material, etc.). Primers and reagents aliquoting is recommended to
1615 minimise contamination consequences. The air pressure should be positive and blow
1616 out of the room. The clean areas must be kept free of amplicon at all times, to ensure
1617 this occurs there should be no movement back from the dirty area to the clean area. If
1618 under extreme circumstances a consumable or reagent needs to be moved backwards
1619 it must be thoroughly decontaminated with bleach and ethanol. Returning racks
1620 should be soaked in 1% bleach overnight before soaking in distilled water and placing
1621 in the clean area. To ensure minimal movement between areas during the running of
1622 molecular assays, it is optimal to have dedicated storage (freezer, fridge and room
1623 temperature) for each area.
- 1624 - **Nucleic acid extraction room:** in this area samples are processed, reverse
1625 transcriptase step of RT-PCR is performed and DNA or cDNA and positive controls
1626 are added to the PCR reaction mix (prepared in the Reagent preparation clean room).
1627 The air pressure should be positive and blow out of the room. If chemicals are stored
1628 in this area appropriate facilities and storage requirements should be in place.
- 1629 - **Amplification room:** PCR machines are housed in this room. It may contain an area
1630 (cabinet with air pressure slightly positive) for the nested PCR.
- 1631 - **Product analysis room:** post-PCR manipulations such as agarose gel electrophoresis
1632 are performed in this area. It is thus a contaminated area and therefore no reagents,
1633 equipment, coats, etc. used in this room should be used in any other PCR areas. The
1634 air pressure should be negative and blow into the room.

1635 3.3 Air handling, water supply, environmental control, heating and cooling

1636 Air handling systems should be operated to ensure that the correct environment is maintained
1637 for the type of work conducted in the laboratory. These systems should be subject to regular

1638 maintenance by qualified personnel. In case *in vitro* work involves serious human pathogens,
1639 the laboratory should operate with specific trained personnel, using biosafety level 3 and 4
1640 and the room should be kept at negative pressure to guard against infection spread. In
1641 contrast, work with cell and tissue cultures, for example, may require positive pressure
1642 relative to other laboratory areas to minimise the risk of test system contamination from the
1643 outside. When HEPA filters are used in differential pressure isolation rooms, the filters and
1644 their fittings and seals need to be thoroughly examined and tested at regular intervals (e.g.
1645 annually). Decontamination should be carried out before servicing is carried out. The air
1646 handling system should also be designed to account for exhaust air from the Class II
1647 biological safety cabinets that are vented to the outside of the building.

1648 Cell culture work requires high-quality ultra-pure water, which is usually deionised via
1649 reverse osmosis, followed by passage through a series of carbon and micropore filters
1650 eliminating organic materials and pyrogens. Tissue culture grade water should be controlled
1651 for pH, conductance and total organic carbon, as well as absence of endotoxins. Note that
1652 pyrogens can be deleterious to cell cultures at concentrations below the level of detection for
1653 organic carbon. Where small quantities of purified water are required for cell culture this may
1654 be obtained by obtaining water for irrigation (WFI) or other medically approved pure-water
1655 preparations (Stacey and Davis, 2007).

1656 Heating, cooling and humidity should be adequate for the comfort of laboratory occupants
1657 and for operation of laboratory equipment, and should not adversely affect test system
1658 survival/behaviour and test item stability. For example, in some cases (e.g. preparation of
1659 microscopic slides) specific humidity might be required.

1660 Many tissue culture media components are sensitive to white light (especially sun light). The
1661 blue wavelengths are of particular concern. Filters can be used in the room and laminar flow
1662 cabinet light to reduce this exposure where necessary.

1663 Mid to long term storage of media is usually best at temperatures below ambient lab
1664 temperatures. Accordingly, an optimal solution may be to store all cell culture media at 4°C
1665 (refrigerator) or frozen (electric freezer) as recommended by the manufacturer. There may be
1666 exceptions to this general opinion but manufacturers' instructions should always be consulted.

1667 **3.4 Cell and tissue culture transportation and cryostorage**

1668 Because cells and tissue in culture are often transported across the world, it is very important
1669 to keep these test systems as healthy as possible during the long transport times. A quality
1670 cell culture incubator makes it possible to transport valuable and delicate cell cultures in a
1671 temperature-controlled environment, so that they are less likely to become damaged during
1672 the transport process. However, one should note that the incubator for transport has limited
1673 space (2-3 plates or flasks) and adequate sealing of plates is necessary to avoid leakage
1674 during transport. Because of their fragile nature, live cells and tissues cannot be shipped like
1675 a regular package in a cardboard box or shipping envelope. They must be transported in the
1676 special temperature-controlled environment of a mini-cell culture incubator, if they are
1677 expected to reach their destination in good condition. Technologies have provided such
1678 solutions during the last decennia. Alternatively, good low temperature transport practices
1679 can allow also for safe transport of cells and tissue systems. Cells are usually shipped on dry
1680 ice. Ideally, temperature should also be monitored (e.g. by using data-logger) during
1681 transportation, especially for long distance transport.

1682 Cryostorage systems should ensure the long term preservation of the stored test system. For
1683 cryopreserved cell cultures, the viability of mammalian cells is progressively lost within

1684 months at -80°C, thus, long term storage below the glass transition point of water (-136°C) is
1685 recommended. While true for mammalian cells, this is not the case for bacteria or yeast.

1686 Storage in the vapour phase of liquid nitrogen is generally advised for all cells and necessary
1687 for potentially infectious cells and tissues. This eliminates the chances of transfer of
1688 pathogenic material between vials which can occur in the liquid phase of nitrogen (Coecke *et*
1689 *al.*, 2005), Appendix 1). It is also considered safer as liquid nitrogen can enter storage vials if
1690 they are stored in the liquid phase and cause them to explode upon thawing. If vials need to
1691 be stored in the liquid phase, protection wrapping may be considered.

1692 Cryostorage requires temperature and liquid nitrogen level monitoring to ensure that the test
1693 system stocks are at optimal storage temperature. Cryostorage vessels can be fitted with
1694 alarms and data loggers and liquid nitrogen levels recorded at regular intervals (e.g. weekly).

1695 In cases when ultra-low electrical -150°C freezers are used, CO₂, liquid N₂ or electrical
1696 backup systems need to be in place to guard against loss of power supply.

1697 Storing valuable test system stocks in more than one cryostorage location is recommended
1698 for security/back up purposes, and off-site storage may also need to be considered in disaster
1699 recovery plans for the facility.

1700 **3.5 Quarantine for new test systems**

1701 All new cells and tissues should be quarantined in the laboratory and in storage until
1702 determined free of contaminating microorganisms (see Section 5). Early checks of cell
1703 authentication are also recommended to avoid wasted time and resources on unauthentic cell
1704 lines. If a separate quarantine laboratory is not available, a quarantine biological safety
1705 cabinet, a dedicated incubator and liquid nitrogen tank can be used. Alternatively, other steps
1706 can be taken to minimise contamination risks, such as handling the quarantine cells last on
1707 each day, rigorous post-manipulation disinfection of the work areas and placing cultures for
1708 incubation in a filter-sealed container into the general incubator (Geraghty *et al.*, 2014). Cells
1709 procured from a cell bank may carry a certificate of analysis for contamination tests
1710 performed. The certificate will list the tests performed and may give details of testing
1711 methods. As a minimum, a mycoplasma test should be performed upon receipt.

1712

1713 4 Apparatus, material and reagents

1714 4.1 Apparatus

1715 Apparatus, including validated computerised systems, used for the generation, storage and
1716 retrieval of data, and for controlling environmental factors relevant to the study should be
1717 suitably located and of appropriate design and adequate capacity. In general, all apparatus
1718 used should be operated and maintained correctly by trained staff to ensure correct use and
1719 knowledge on the procedures for equipment maintenance and calibration.

1720 Computerised systems should be formally validated prior to use in a GLP environment (see
1721 10.1). Apparatus should be periodically inspected, cleaned, maintained, and calibrated
1722 according to SOPs and records of these activities should be maintained (OECD, 1998b).

1723 As stated in OECD document *Number 14* (OECD, 2004b), the commonly observed, routine
1724 requirements for apparatus used in a GLP environment apply equally to apparatus used for *in*
1725 *vitro* development work, like laminar flow cabinets and incubators, there are some specific
1726 points and issues of particular importance for the latter. As an example, equipment such as
1727 microbalances, plate readers, centrifuges, micropipettes, laminar air flow cabinets, fridges
1728 and freezers, water baths, and incubators should be regularly maintained and calibrated
1729 (traceable to international standards) where possible, to ensure the integrity and reliability of
1730 the results. For each type of equipment, critical parameters (e.g. supply of gases for mass
1731 spectrometry; liquid nitrogen levels in storage containers, CO₂ levels in incubators; or
1732 functioning of filters to ensure sterility of the air in addition to the airflow in the laminar flow
1733 cabinets) should be identified as requiring continuous monitoring or the setting of limit
1734 values together with installation of alarms.

1735 Centrifuges which are routinely used in cell and tissue culture work (sub-culture work,
1736 cryopreservation etc.) produce aerosols and therefore it is important to consider models that
1737 have sealed buckets. Ideally, one should consider working with models where the condition
1738 of the load can be observed without opening the lid. Besides the containment issues for
1739 centrifuges, when documenting their operation, it is necessary to specify centrifugation
1740 speeds as x g rather than rpm (unless the rotor radius is stated), incubation conditions, time
1741 and volumes of centrifugation with tolerances when relevant, and any other information that
1742 enables the accurate reproducibility of procedures. In addition, procedures should be
1743 established on how to respond in case of an emergency (e.g. broken tube).

1744 Working with cell and tissue culture *in vitro* requires a strictly controlled environment for
1745 cell growth. This is achieved using specialised incubators which provide the correct and
1746 controlled growth conditions (temperature, humidity, CO₂ levels), which should be checked
1747 (and logged) on a regular basis. To provide the required humidity, incubators are frequently
1748 equipped with a water basin at the bottom. This combination of high humidity and
1749 temperature increases the risk of bacterial or fungal contaminations. To reduce the risk of
1750 microbial contamination, copper-coated incubators are now available. Antifungal or
1751 bactericidal agents can be added to incubator water trays to reduce the risk of bacterial and
1752 fungal growth (but any possible impact on the *in vitro* method to be carried out should be
1753 checked and documented). Incubators with self-sterilising cycles may also be used, although
1754 this does not replace regular cleaning and maintenance.

1755 Similarly, waterbaths used to thaw and/or to warm up stored solutions like medium and
1756 frozen stocks, or to defrost vials of cryopreserved cells and tissues, carry a high risk of
1757 introducing contamination. The use of bactericidal and fungicidal agents in waterbaths can

1758 aid in the control of contamination, but the impact on the test system should be checked and
1759 documented. Alternatively, cleaning procedures and regular change of water may be a better
1760 choice.

1761 A laminar air flow cabinet can be considered as a critical piece of equipment for cell and
1762 tissue culture work, since, when it is used correctly (see 3.1), it ensures a clean working
1763 environment protection for both the operator and for cells/tissues and other materials and
1764 reagents.

1765 For equipment such as refrigerators and freezers, temperatures should be checked regularly
1766 and preferably logged: simple data loggers are available to log the temperature at set
1767 intervals. In addition to the regular recording of temperatures, an alarm system to alert staff
1768 when acceptable operating limits are exceeded is desirable, and a backup system should be in
1769 place, such that materials may be transferred from one fridge/freezer to another, in case of
1770 malfunction.

1771 For all equipment used during an experiment, acceptable operating limits should be set,
1772 monitored and recorded. Equipment should be fit for purpose with respect to sensitivity and
1773 selectivity. Prior to use, it should be established that the equipment is functioning according
1774 to the (suppliers') specifications and it should be qualified and validated for its intended use,
1775 e.g. via a formal DQ/IQ/OQ/PQ process. Since DQ/IQ/OQ/PQ are not specifically defined in
1776 GLP, and its application is not always harmonised, several test facilities follow the OMCL
1777 guidelines for qualification of equipment (OMCL, 2011). The equipment needs to be
1778 maintained and calibrated at regular intervals depending on the type of equipment and the
1779 frequency of use. As an example, pipettes or micropipettors may need to be calibrated more
1780 frequently than centrifuges.

1781 It may be necessary to have separate procedures for regular checks (e.g. daily checks of
1782 pipettes used) and complete calibration (monthly or quarterly, depending on the frequency of
1783 use). The standard operating procedures need to describe how to deal with data produced
1784 with a potentially uncalibrated pipette. When equipment such as a pipette is out of
1785 specifications during a calibration, it is important to determine how to interpret data that have
1786 been generated since the most recent successful calibration and determine the impact of the
1787 potential deviation to the outcome of the study. Therefore, it is crucial to record every piece
1788 of equipment, uniquely identified, that has been used during the performance evaluation of an
1789 *in vitro* method. In general, facility practices should ensure that equipment is within
1790 specifications before the start of study and throughout the experimental phase to avoid the
1791 rejection of the *in vitro* study. Nevertheless, some test facilities use periodical calibrations for
1792 pipettes (e.g. every 3 months), given their low failure rate.

1793 To enable broader use of a new method, successful transfer to a range of equipment and
1794 different laboratories should be demonstrated. This increases the robustness of the method.
1795 To increase transferability, preference should be given to the use of generally widely
1796 available equipment. In addition, the impact of the use of a certain type or brand of
1797 equipment on the outcome of the individual assays or the overall *in vitro* method needs to be
1798 determined. An *in vitro* method should specify the acceptable limits that the equipment
1799 should meet to be used for a specific *in vitro* method¹⁵.

¹⁵ http://www.iso.org/iso/iso_catalogue/catalogue_tc/catalogue_tc_browse.htm?commid=48908

1800 4.2 Materials and reagents

1801 Critical materials and reagents can influence the outcome of *in vitro* studies and it is therefore
1802 important to identify and control them in each experimental phase. In using the term critical
1803 materials and reagents we refer to highly complex preparations used as a fundamental
1804 component in the *in vitro* method, which are not currently definable in terms of composition,
1805 stability or biological activity. Examples include processed tissue and cell culture
1806 preparations, which may contain viable cells. As with all reagents, the standards of
1807 preparation and testing of kits should be clear and traceable and they should be used before
1808 their expiry date. Further guidance on complex cell systems such as 3D culture, is given in
1809 Annex 2. For GLP studies, preparation of substrates and use of reagents and kits should be
1810 documented in sufficient detail to allow for complete reconstruction of these activities.

1811 Other examples of critical reagents are serum and growth factors. Each of the critical reagents
1812 should either be sourced from a reputable supplier, who accompanies the shipment with a
1813 certificate of analysis, or one should ensure that there are appropriate quality controls (see
1814 Annex 2). These controls may include in-house growth or functional characterisation and the
1815 controls should be done by trained personnel and according to procedures described at *in*
1816 *vitro* method development laboratories or in the standard operating procedures (SOPs) of
1817 routine or GLP facilities. It is recommended to perform a batch validation to reduce the
1818 introduction of unknown variables into a culture system that can interfere with assay or
1819 overall *in vitro* method performance. For this purpose a batch is tested first and when
1820 approved, a large quantity of the batch can be acquired to reduce variability during the
1821 performance of a certain number of assays.

1822 The use of good quality assurance practices applies to all laboratory reagents and
1823 consumables. For instance, preparation of reagents should follow a standardised and fully
1824 documented methodology that lists supplier information, lot numbers of component
1825 chemicals and reagents, dates of preparation, and the names of the staff involved in the
1826 preparation. Substances should be tested against reagents and standards of known reactivity
1827 before being released for use, and the results of these tests must be recorded (on, for example,
1828 special reagent preparation forms). All prepared solutions should have unique lot numbers (a
1829 laboratory-based system is acceptable as long as the provided numbers are actually unique).
1830 This information log prevents duplication and should comprise lot numbers, dates,
1831 descriptions, expiry dates, and signatures. For both reagents and reagents mixtures, the
1832 container should be inert to the stability of the substance or mixture and clearly labelled with
1833 the following details:

- 1834 • contents or identity
- 1835 • potency (titre, concentration, or activity, for example)
- 1836 • storage temperature
- 1837 • preparation date
- 1838 • unique lot number
- 1839 • container number (in case there are multiple containers of the same lot)
- 1840 • date first opened (as appropriate)
- 1841 • expiration date (determined by experimentation or reference to manufacturer's
1842 recommendation)
- 1843 • signature of the person who prepared the contents.

1844 Much of this information can be recorded separately with a unique identification number log.

1845 Even when reagents are sourced from a reputable supplier, it remains important to assure the
1846 stability of the reagents during shipment conditions, in addition to the storage. For example,

1847 reagents shipped frozen should arrive frozen and this should be documented on the receiving
1848 document. The presence of a data logger is the best practice in these cases.

1849 Storage should be done according to the manufacturer's specifications in the supplied
1850 certificate. Most solutions which come in a large quantity should be aliquoted, in order to
1851 minimize the number of times a bottle is opened and thus minimize risk of contamination
1852 (and spread of contamination). This is particularly important for solutions which require
1853 storage below 0°C, in order to avoid repeated freeze-thaw cycles. When reagents need to be
1854 thawed and possibly frozen again, it is recommended to determine the number of freeze-thaw
1855 cycles that the reagents can withstand (EMEA, 2011; FDA, 2001; Viswanathan *et al.*, 2007).
1856 Stability of aliquots should be verified in the lab performing the *in vitro* method and not rely
1857 solely on literature data.

1858 In these cases, quality controls need to be performed according to pre-defined procedures
1859 described in SOPs. Normally, stability of the analyte in the studied matrix is evaluated using
1860 at least triplicate samples of the low and high concentrations, which are analysed
1861 immediately after preparation and after the applied storage conditions that are to be
1862 evaluated. The thawed samples are analysed against a calibration curve, obtained from
1863 freshly prepared calibration standards, and the obtained concentrations are compared to the
1864 nominal concentrations. The deviation should be within previously established acceptance
1865 criteria (usually $\pm 20\%$ for large molecules). It is absolutely necessary that the number of
1866 cycles in the freeze/thaw stability evaluation should equal or exceed that of the freeze/thaw
1867 cycles of study samples.

1868 When pipetting problematic volatile/viscous liquids or suspensions, it is strongly
1869 recommended to use positive displacement pipettes. Certain chemicals may exhibit non-
1870 specific adsorption to the plastic tips of pipettes and the use of low-binding materials
1871 (including glass) or acoustic droplet ejection (Ekins *et al.*, 2013; Grant *et al.*, 2009) can be
1872 utilised to alleviate these issues.

1873 **4.3 Basal medium**

1874 Depending on the circumstances, the basal culture medium can be serum-supplemented (as in
1875 traditional cell culture methods) or serum-free, but supplemented with additives necessary for
1876 obtaining satisfactory cell proliferation and production, or for maintaining a desired
1877 differentiation status. Many slightly different formulations exist under the same general
1878 medium names, such as Minimum Essential Medium (MEM), and even subtle changes in the
1879 medium formulation can substantially alter the characteristics of certain cells and tissues. In
1880 many cases, these variations are deliberate for specific applications. Therefore, the medium to
1881 be used should be precisely specified, and it is important to check that new supplies of
1882 medium meet the required specifications (Coecke *et al.*, 2005).

1883 **4.3.1 The use of serum in cell culture**

1884 The use of serum has been discouraged in recent years due to the necessity of animals to
1885 produce it, while *in vitro* methods are mostly developed to replace animal use. Furthermore,
1886 serum contains many undefined factors and these could change every time a new batch of
1887 serum is ordered, even if it is from the same source. In this respect, serum starvation, use of
1888 serum-free media and serum replacements has become standard practice in many
1889 laboratories.

1890 Despite the availability of serum-free media and serum replacements, serum is still used in a
1891 lot of *in vitro* development work. Animal serum can be derived from adult, new born or

1892 foetal sources. Bovine sera are most commonly used, and during the last few decades, foetal
1893 bovine serum (FBS) has become the standard supplement for cell culture media. However,
1894 each new batch of serum may contain different concentration of growth factors and
1895 hormones, Therefore new batches should be tested on a relevant range of cell lines for cell
1896 attachment, spreading, cloning efficiency, growth rates and activity in functional assays
1897 (Geraghty *et al.*, 2014).

1898 Cell lines which have been derived or cultured in serum-containing media long-term may
1899 become dependent on the multitude of growth factors present in serum and may experience a
1900 phenotypic drift upon abrupt serum withdrawal. This may manifest as growth arrest or
1901 activation/inactivation of various signalling pathways. Serum can also be used to mimic the
1902 protein binding occurring in the blood *in vivo*, but this is a process for which specific
1903 attention and calculation should be taken into account. Specific test item aspects need
1904 consideration in choosing to work with serum: 1) If the test item is known to bind to protein,
1905 its effect might not be seen unless a very high concentration of test item is used (see section
1906 6.3 on biokinetic parameters); 2) If the test item antagonises an endogenous circulating
1907 hormone or factor, the serum might contain such hormone or factor and may thus affect the
1908 assay results.

1909 Test developers determine serum specifications that meet their particular needs (e.g. to allow
1910 the cells to function like in the *in vivo* situation, as much as possible, to enhance test result
1911 utility) and match the natural behaviour of the cells as much as possible, including the
1912 maximum acceptable levels of serum components, such as immunoglobulins (which may
1913 have inhibitory effects), endotoxins (indicative of bacterial contamination, but are also
1914 powerful cell mitogens), and haemoglobin (indicative of haemolysis during clotting).

1915 Animal sera are a potential source of microbiological contaminants, notably mycoplasma,
1916 bovine viruses, and possibly the agent which causes Bovine Spongiform Encephalopathy
1917 (BSE). Consideration should be given to impact of contamination with the most common
1918 viruses in bovine herds such as BVDV and BPV agents (WHO 2010). Suppliers use a variety
1919 of techniques, including filtration, irradiation and heat-inactivation, to reduce microbial
1920 contamination. Nevertheless, it is wise, and for some applications, obligatory, to specify
1921 sourcing of serum from countries where there is a low risk of infection, and, in the case of
1922 bovine sera, from not too old animals to reduce any potential risks (Festen 2007). The use of
1923 human serum is restricted to specialised applications as it carries additional risks, such as the
1924 potential presence of human pathogenic viruses. Its use must be subject to the strictest quality
1925 controls, including documentation to demonstrate origin and viral safety (Coecke *et al.*,
1926 2005).

1927 Variability between serum batches can lead to experimental variability and reduce inter-
1928 laboratory reproducibility, representing a major cost associated with cell culture (Usta *et al.*,
1929 2014). Notably, the use of serum can possibly lead to unexpected or undesired outcomes. In
1930 this regard, a study has shown that serum can inhibit transforming growth factor (TGF)- β 1-
1931 induced chondrogenesis in fibroblast-like type-B synoviocytes (Bilgen *et al.*, 2007). In
1932 another study, FBS, compared to autologous (human) serum, was found to induce a more
1933 differentiated and less stable transcriptional profile in human bone marrow mesenchymal
1934 stem cells, particularly at late passages, as shown by analysis of genome-wide microarray
1935 analysis (Shahdadfar *et al.*, 2005).

1936 Moreover, EU legislation, through the Directive 2010/63/EU on protection of animals used
1937 for scientific purposes, offers a certain amount of protection for foetal forms of mammals.
1938 Considering that blood harvesting from foetal calves can cause pain and distress in these
1939 animals, as previously reported (Jochems *et al.*, n.d.; van der Valk *et al.*, 2004), and

1940 considering Articles 1 and 3 of the Directive, harvesting of FBS from live bovine fetuses in
1941 the last third of their development for scientific purposes is a procedure under the Directive.

1942 In 2008 the ECVAM Scientific Advisory Committee (ESAC) stated that "*for methods*
1943 *forwarded to ECVAM for validation/prevalidation where [the use of non-animal alternatives*
1944 *to serum] is not fulfilled a justification for future use must be provided, including measures*
1945 *taken to seek non-animal alternatives to [FBS]*" (ESAC, 2008) The drawbacks of using FBS
1946 and recommendation to replace it with available chemically defined serum free media is also
1947 mentioned in the 2005 GCCP guidance issued by EURL ECVAM (Coecke *et al.*, 2005).

1948 Additionally, FBS is a by-product of the meat industry and for this reason is dependent on
1949 several external and generally unpredictable factors, such as weather conditions (e.g.,
1950 changes of climate with drought or flood) or variations in beef consumption (e.g., switch to
1951 poultry). All these factors may make FBS supply unpredictable (Brindley *et al.*, 2012).

1952 For all these reasons, serum-free media, (animal) protein-free media and chemically-defined
1953 media has gained interest over the years in order to avoid many of the disadvantages
1954 associated with the use of serum, as recently commented (Pamies *et al.*, 2016).

1955 **4.3.2 Serum-free media and serum replacements**

1956 Serum-free media are thought to circumvent the batch-to-batch variability issues associated
1957 with serum and offer better reproducibility and the potential for selective culture and
1958 differentiation of specific cell types (Geraghty *et al.*, 2014). Nevertheless, serum-free
1959 compositions may still need to be validated and monitored similarly to serum containing
1960 media as they are often not completely chemically defined. For instance, serum free media
1961 can include other poorly defined supplements such as pituitary extracts, chick embryo
1962 extracts, bovine milk fractions or bovine colostrum. Furthermore, some so-called 'defined'
1963 media contain complex serum replacement mixtures including chemically undefined agents.
1964 Notably B27 and its alternative NS21 used in the culture of neural cells contain bovine serum
1965 albumin and transferrin which can exhibit batch-to-batch variation in biological activity
1966 (Chen *et al.*, 2008). Therefore, it may be useful to carry out pre-use testing on new batches of
1967 reagents which could demonstrate variability that cannot be foreseen from manufacturers'
1968 information. Another example of an essential component prone to batch-to-batch variability
1969 is the so-called 'basement membrane extract', purified from Engelbreth-Holm-Swarm (EHS)
1970 mouse sarcoma cells and marketed under various trade names.

1971 Serum-free medium formulations play a critical role for the culture of stem cells, such as
1972 human embryonic stem cells (hESCs) and human induced pluripotent stem cells (hiPSCs),
1973 which currently hold great promise for toxicology and regulatory testing and biomedical
1974 research. Both hESCs and hiPSCs are often maintained on inactivated mouse or human
1975 embryonic fibroblasts or under feeder-free conditions (using extracellular matrices) in
1976 chemically defined, serum-free media, in order to avoid the presence of undefined or
1977 unknown serum components (which may compromise the differentiation towards desired cell
1978 lineages) and the risk of contaminations from pathogens (e.g., mycoplasma, viruses, and
1979 prions) (Pistollato *et al.*, 2012; Yamasaki *et al.*, 2014).

1980 **4.4 The use of antibiotics in cell culture**

1981 Routine culture of cell lines under GCCP (Coecke *et al.*, 2005; Geraghty *et al.*, 2014; Stacey
1982 and Davis, 2007) should not require the use of antibiotics and it can never be relied on as a
1983 substitute for effective aseptic techniques. However, its use is still widespread due to
1984 established routine procedures in many laboratories. Antibiotics are agents that may arrest or

1985 disrupt fundamental aspects of cell biology, and, while they are effective against prokaryotic
1986 cells (i.e. bacteria), they are also capable of causing toxic effects in animal cells. Not
1987 surprisingly, antifungal agents, being directed at higher order, eukaryotic micro-organisms,
1988 are likely to be more toxic to animal cell cultures. Given these obvious contra-indications, the
1989 use of antibiotics in cell and tissue culture should be focused in two areas: a) protection of
1990 materials at high risk of contamination such as tissues, organs and primary cultures in cases
1991 where sterility cannot be guaranteed; and b) the positive selection of recombinant cell clones
1992 based on the expression of antibiotic resistance genes (Coecke *et al.*, 2005).

1993 **4.5 Additional media components**

1994 Some media components are heat labile (e.g., L-glutamine), sensitive to light (e.g., retinoic
1995 acid) or have a limited half-life in diluted state or at high ionic strength, such as in prepared
1996 media (e.g., epidermal and fibroblast growth factors). These issues are best addressed by
1997 preparing a small volume of media necessary to cover the period of stability of the most
1998 sensitive component and discarding bottles after a set time period. Appropriate size aliquots
1999 of those labile components may be frozen by an appropriate method for long-term storage. In
2000 this respect, stock solutions with low concentrations of protein aqueous growth factors may
2001 require the addition of albumin or other excipients, to prevent adsorption to plastic and to
2002 increase stability in the frozen state. Stabilised forms of glutamine and retinoic acid are also
2003 available to avoid these issues.

2004 In case culture media or other reagents have to be sterilised via heat or filtration, the impact
2005 of the procedure should be assessed and recorded. For example, heat sterilisation may result
2006 in degradation (or denaturation) of one or more of the components and filtration can remove
2007 individual and/or essential components (e.g. Fe²⁺ or Fe³⁺ iron products that enhance growth
2008 of mammalian eucaryotic cells in serum-free cultures).

2009 **4.6 Dedicated media to particular cell lines**

2010 Different cell types or tissues need to be cultured in culture media containing various
2011 components at different concentrations to allow optimal growth. Although certain cell lines
2012 may be grown in media with the same composition, sharing media increases the risk of cross-
2013 contamination. Therefore, each cell line should be cultured with separate dedicated media,
2014 which must not be shared with other cell lines.

2015 5 Test Systems

2016 *In vitro* method cell and culture-based test systems are routinely used by all industries and
2017 regulatory bodies in toxicity testing, safety assessment, and risk evaluation. The greatest use
2018 of *in vitro* test systems, however, is for elucidating mechanisms of toxicity and/or
2019 demonstrating the biological process involved, when exposing test systems to toxicants of
2020 various kinds. With the advances in genetics and genetic screening approaches, routine *in*
2021 *vitro* methods include already the use of either genetically altered cells, stem cells, stem-cell-
2022 derived models and organ-on-chip models and other complex and sophisticated systems.

2023 5.1 GCCP

2024 Good Cell Culture Practice identifies a set of core principles of best practice for working with
2025 simple but also with more complex cell and tissue culture systems (Annex 1 and Annex 2).
2026 Failure to comply with GCCP can have serious consequences for individual researchers and
2027 their employers, which include invalidation of important data sets, ethically compromised
2028 research, loss of crucial cultures due to microbiological contamination, failed patent
2029 applications and injury or infection amongst laboratory workers. These could clearly have
2030 serious consequences for the institutions involved, including the threat to scientific
2031 reputation, legal and regulatory compliance and cost of wasted resources.

2032 The principles of the first version of GCCP published in 2005 (Coecke *et al.*, 2005) remain
2033 highly relevant to cell culture practice for *in vitro* methods today (see Annex 1), and aims at a
2034 far broader set of applications, including research, manufacture of medicines, and laboratory
2035 based GLP testing. In providing the detailed and specific principles of best practice for the
2036 handling and management of cell cultures systems, GCCP is a vital component of GIVIMP.

2037 A new and updated version of GCCP called “GCCP2.0” is under development through
2038 collaboration between European, Japanese and North American experts and incorporates
2039 more recent developments in cell culture, and in particular addresses the new technological
2040 developments of human pluripotent stem cell lines and complex 3D culture systems. During a
2041 workshop in 2015 (see Annex 2) new key elements for GCCP 2.0 were identified.

2042 5.2 Cell and tissue sourcing

2043 Complex biological substrates such as cell and tissue cultures may have passed through a
2044 number of stages of development before they are qualified for an *in vitro* method and it is
2045 important that this track is documented, thus assuring traceability to original source materials,
2046 source cell banks, absence of contamination by major classes of biological agents (e.g.
2047 mycoplasma, bacteria, fungi and cytopathic viruses), genetic identity/consistency and
2048 stability of desired functionality. See Annex 1, GCCP principle 3 and [Table 1](#) in
2049 section 1.2 for examples of document requirements concerning the origins of cells and tissues
2050 and see paragraph 5.5 below.

2051 Proposed sources of test systems should be qualified by the user to assure they are
2052 appropriately traceable, have been quality controlled for key features (see paragraph 5.7
2053 below) and in addition, the user should check that there is solid ethical provenance (e.g.
2054 hPSCreg registry¹⁶) and safety assessment for the cells. In addition, intellectual property

Formatted: Font: Times New Roman,
12 pt

¹⁶ <http://www.hpscereg.eu/>

Field Code Changed

2055 should be checked to ensure it does not impact on the use of the cell line and data. For more
2056 detailed information on these issues see (Stacey *et al.*, 2016).

2057 In the case of human tissues and primary cells, there will also be a requirement to assure
2058 donor consent and to manage sensitive personal data e.g. according to the EU Data Protection
2059 Directive. Human tissue is precious and difficult to secure. The broad range of issues in
2060 securing tissues for testing were addressed at the 32nd Workshop of the European Centre for
2061 Validation of Alternative Methods (Anderson *et al.*, 1998) and where tissues cannot be
2062 sourced via a qualified tissue bank, there should be an agreed testing method in place with
2063 clinical contacts regarding all aspects of harvesting, preparation, labelling, storage and
2064 transfer (for an example see (Stacey and Hartung, 2006)). It is also important to assess the
2065 risks of viral contamination of primary cells and tissues. Approaches to risk assessment are
2066 described in (Stacey and Hartung, 2006). Tissues should be obtained from tissue banks
2067 holding only materials from screened donors and this will significantly assist in managing
2068 viral safety issues.

2069 Availability of cell lines from a certified source (established cell banks with a high quality
2070 standard; reputable culture collections, commercial provider), that usually provide extensive
2071 documentation on the origins and characterisation of the test system, should be assured¹⁷.
2072 Moreover (or alternatively), master and working cell banks should be established to
2073 guarantee a cell supply of constant quality and records kept of the original source.

2074 If test systems used in *in vitro* studies are genetically modified the Directive 2009/41/EC
2075 (EU, 2009) is applicable. This Directive lays down common measures for the contained use
2076 of genetically modified micro-organisms (GMMs), aimed at protecting human health and the
2077 environment. A notification has to be sent to the competent authorities before any use
2078 commences in the premises. A risk assessment of the GMMs used has to be performed. The
2079 Annexes to the Directive detail the criteria for assessing the risks of GMMs to health and the
2080 environment, as well as the protective measures for each of the four levels of containment.
2081 The Directive lays down the minimal standards applicable to the contained use of GMMs.
2082 Individual European member States are also permitted to take more stringent measures.

2083 **5.3 Handling and maintenance of the test system**

2084 During routine handling and maintenance, growth and survival characteristics of the cell
2085 system (such as cell viability, doubling time, etc.) and subculturing details (e.g. date of
2086 subculture, subculture intervals, seeding density, passage number, etc.) should be recorded
2087 and documented in the study report, since they are required for the complete traceability of
2088 results. The documentation provided by the test system supplier should be taken into account
2089 together with the historical data, when available, and used to establish acceptance criteria.
2090 See [Table 1](#) for examples of document requirements concerning the handling,
2091 maintenance and storage of cells and tissues.

2092 Different cell lines have different growth rates which may depend on several environmental
2093 factors. Whether cells grow and divide in a monolayer or in suspension, they usually follow
2094 the same characteristic growth pattern in which four different phases can be recognized: lag,
2095 log (or exponential), stationary (or plateau) and decline. Growth during exponential growth
2096 or log phase is usually fairly constant and reproducible for a given set of growth conditions
2097 (ATCC, 2014). Each cell line will show different cell proliferation kinetics during the log

Formatted: Font: Times New Roman,
12 pt

¹⁷ http://wiki.toxbank.net/w/images/1/18/ToxBank_D4_6_final_10_04_13.pdf

Field Code Changed


2098 phase and it is therefore the optimal phase for determining the population doubling time
2099 (ECACC, 2010).

2100 Many dividing primary human cell cultures have a split ratio of one in two (1:2), while
2101 continuous cell lines have much higher splitting rates. In order to ensure viability, genetic
2102 stability, and phenotypic stability, cell lines need to be maintained in the exponential phase,
2103 i.e. they need to be subcultured before a monolayer becomes 100% confluent or before a
2104 suspension reaches its maximum recommended cell density.

2105 Many cell lines can be subcultured based on a rough estimate of cell density, and this is the
2106 usual practice unless stated otherwise in the cell maintenance protocol. Some cell lines
2107 require a fixed seeding density and subculturing scheme and counting the number of cells is
2108 required (Wilson *et al.*, 2015). Most commonly cell counting is performed using the Bürker
2109 Türk or Neubauer counting chambers. When automated cell counters are used, their correct
2110 functioning would need to be demonstrated (Cadena-Herrera *et al.*, 2015; Gunetti *et al.*,
2111 2012; Phelan and Lawler, 2001).

2112

2113 **Figure 3: Growth curve for cells grown in culture. Cells should be subcultured while still**
2114 **in the exponential phase (ATCC, 2014)**



2115

2116

2117 Each test facility should develop SOPs, where details are provided about how to thaw,
2118 handle, count, maintain and store cell lines. For example, the procedure to univocally assign
2119 progressive passage numbers as well as the selected assay(s) to determine the cell stock
2120 viability should be established.

2121 5.4 Cryopreservation

2122 Improved technologies that allow cryopreservation of *in vitro* cell and tissue cultures at
2123 different stages of differentiation, and their long-term storage has introduced new or more
2124 standardised *in vitro* test system into the pipeline of potential *in vitro* methods to be used in
2125 human safety assessment. Controlled-rate and slow freezing, also known as slow
2126 programmable freezing have been used all over the world for freezing down cell and tissue
2127 cultures to better preserve it for eventual thawing, before it is frozen, or cryopreserved, in
2128 liquid nitrogen. New methods are constantly being investigated due to the inherent toxicity of
2129 many cryoprotectants.

2130 As described in GCCP Principle 1 'Establishment and maintenance of a sufficient
2131 understanding of the *in vitro* system and of the relevant factors which could affect it'
2132 (Coecke *et al.*, 2005) it is important to prepare preserved banks of cells intended for use, to
2133 assure that reliable stocks can be obtained for testing, which are at a consistent passage level
2134 from the original 'seed stock'. This is in order to avoid the effects of changes or cross-
2135 contamination which may occur if cell lines are maintained indefinitely. Standard
2136 cryopreservation methods using DMSO (10%) and serum (20%) as cryoprotectants,
2137 combined with a slow cooling rate (e.g. -1°C/min), will usually be successful for most cell
2138 lines. However, it is important to check the viability of preserved stocks in case of freezing
2139 failure and also to try to assure consistency between individual vials in a cell bank regarding
2140 cell number, viability and desired functionality. It is important to note that viability
2141 measurements made immediately post-thaw can give misleadingly high values as many cells
2142 can be lost during the 24 h recovery phase post thawing. Therefore, it may be important to
2143 understand what losses occur at this stage.

2144 5.5 Cell line identity and genetic aberrations

2145 Ideally, there should be traceability to the original provider of the cell line and the related
2146 documentation. However, a frequent problem in the use of cell culture is the use of cell lines
2147 which have become cross-contaminated, misidentified (see ICCLAC database of cross
2148 contaminated or misidentified cell lines¹⁸), mixed-up, or underwent genomic instability
2149 (Allen *et al.*, 2016; Fratini *et al.*, 2015; Kleensang *et al.*, 2016; Vogel, 2010). This is not
2150 always detectable by cell morphology only. Establishing an early stock (or retention of a
2151 sample of original tissue) which is DNA fingerprinted will provide an important reference for
2152 future cell banks and for other centres. Short Tandem Repeat (STR) profiling is typically
2153 applied and has considerable background qualification for use in human samples (ISCBI,
2154 2009). STR analysis can be performed in most laboratories that have the capabilities to
2155 execute molecular techniques. It is an easy, low cost and reliable method for the
2156 authentication of human cell lines. An ANSI US standard exists for this technique with cell
2157 lines. For non-human samples, STR methods might be available but may need to be qualified
2158 for specificity to individual cell lines (ASN-0002: Authentication of Human Cell Lines:
2159 Standardization of STR Profiling can be found at <http://webstore.ansi.org>). Other techniques
2160 have also been used including isoenzyme analysis but probably the most common and well-
2161 qualified technique is cytochrome oxidase gene sequencing (Ono *et al.*, 2007). However, the
2162 field of genetic analysis is progressing rapidly and interested parties should maintain
2163 knowledge of current best scientific practice in this area as next generation sequencing begins
2164 to become a routine tool.

Field Code Changed

Field Code Changed

¹⁸ <http://icclac.org/databases/cross-contaminations/>

Field Code Changed

2165 Genetic instability is inherent in cell cultures and it is wise to minimise the number of
2166 passages over which cells are maintained (typically p15-20). Although passage number alone
2167 is not a reliable parameter to ensure good cell functioning, it is good practice to define a limit
2168 for the maximum number of passages, possibly in combination with defined performance
2169 characteristics. At that limit, new cultures should be restarted from a working cell bank. The
2170 use of cells at higher passage numbers must be justified and their integrity and functionality
2171 demonstrated. Cultures at passage numbers beyond which it is known that the cell line
2172 functionality is maintained should not be permitted for use. Where cells are known to be
2173 extremely unstable, some form of genetic testing, such as karyology or molecular analysis
2174 like single nucleotide polymorphism arrays (aSNP) or comparative genomic hybridisation
2175 (aCGH) may need to be performed. In particular, this applies to recombinant cell lines
2176 including those maintained with antibiotic selection.

2177 There are special issues for stem cells. Stem cell lines may contain a mixture of diploid and
2178 aneuploid cells, which may be unavoidable, but genetic testing (see above) can be used to
2179 screen for progressive change (e.g. between master and working cell banks) which could
2180 impact on the suitability of the cell culture. Human iPSC lines should also be tested for
2181 absence of ectopic expression of reprogramming genes and where produced by non-
2182 integrating vectors, for elimination of the vector.

2183 **5.6 Contaminants screening: sterility, mycoplasma, virus**

2184 Standard sterility tests are published and may be used for cell stocks and cultures. However,
2185 it is important to bear in mind that these are usually based on inoculation of broth cultures
2186 which may not support the growth of all contaminating micro-organisms. Alternative
2187 molecular methods such as identification by PCR and DNA sequencing of ribosomal RNA
2188 may be used.

2189 It is absolutely mandatory to quarantine (see 3.5) and to confirm to be negative for
2190 mycoplasma and free of bacteria, yeast and fungi of all cultures passaged in a particular
2191 laboratory environment.

2192 Viruses may arise as contaminants of cell cultures via the original donor used to produce the
2193 cell line or feeder cells and other biological reagents used in cell culture. They may cause
2194 cytopathic effects, in which case the culture should be discarded, or they may have no effect
2195 and become diluted out when fresh uncontaminated reagents are used. In certain cases they
2196 may establish persistent infections, although this is believed to be rare. Whatever the
2197 outcome, their presence and influence on cell biology may be significant as amongst other
2198 effects they may modify transcription factor networks and alter the cells' biology. To assure
2199 lab worker safety, some organisations require testing of all human cell lines for serious
2200 human pathogens such as human immunodeficiency virus (HIV) and Hepatitis B&C or
2201 evidence that the donors did not have these pathogens. However, such testing clearly does not
2202 cover more common human infections, and human pathogens may also be carried by cells
2203 from other species. Cell cultures should therefore always be handled and disposed of as if
2204 they are potentially infectious (Stacey, (2007) chapter 31: Risk assessment of animal cell
2205 cultures, In : Medicines from Animal Cells (Stacey and Davis, 2007). Cell line testing may
2206 be initiated if there are special hazards associated with the work or with the cells. Workers
2207 should always follow local rules for performing cell culture work, maintain their competence
2208 in aseptic processing, as well as carry out regular and careful inspection of cells for any
2209 unusual effects or morphologies that might indicate infection. It is an important part of a
2210 robust testing regime for contamination to have a procedure for managing positive results,

2211 whether this means immediate discard or quarantine until a means of action can be decided
 2212 along with detection of the root cause supplementary testing (Stacey, 2011).

2213 **5.7 Quality Control**

2214 It is important that certain key go/no-go points are established during the preparation and use
 2215 of the test system for an *in vitro* method.

2216 Typically, for use of cell lines, appropriate integrity checks should be applied at the sourcing
 2217 of new cell lines or cell bank preparation, to assure stable passage and expansion for use and
 2218 reproducible starting cells used to provide the test substrate, and, finally, to ensure acceptable
 2219 and reproducible performance against controls in formal tests. [Table 4](#) shows an
 2220 example of measures at different stages.

Formatted: Font: Times New Roman, 12 pt

2221

2222 **Table 4: Applicability of Integrity Checks on Cultures**

Characteristic	Source *	Early stocks of preserved vials**	Cell banks***	Routine in process testing of stock cultures****
Morphology	X	X	X	X
Viability	X	X	X	X*****
Identity	X	X	X	
Doubling time	X	X	X	X*****
Mycoplasma	X	X	X	X
Viruses	(X donor)		(X master bank)	
Bacteria and Fungi			X	X*****
Function/phenotype		X	X	
Genetic stability			X	(X)
Absence of reprogramming vectors (iPSC lines)		X	X	

2223 *This may be provided in writing but ideally with evidence of test results or a qualified service provider test
 2224 result. These tests should also be performed on the cells arriving in the lab as soon as samples can be obtained.

2225 ** A small number of vials frozen as soon after arrival of a new culture to avoid loss in case of a lab accident.

2226 ***Ideally master and working banks should be established (Coecke *et al.*, 2005) but testing may be focused on
 2227 the master stock with more routine checks applied to working cell banks e.g. mycoplasma and viability.

2228 ****Stock cultures of cell lines maintained for use in *in vitro* methods should be tested routinely e.g. all lines in
 2229 culture tested once per month.

2230 *****Viability testing at passage will also be helpful to ensure consistent seeding of fresh cultures and assays
 2231 for more reliable maintenance of stock cultures and reproducibility of cell-based *in vitro* methods. For this, the
 2232 assays described under 6.2.1.1 can be applied.

2233 *****For diploid cultures, passage number is roughly equal to the number of population doublings (or
 2234 population doubling level) since the culture was started (ATCC, 2014).

2235 *****To avoid development of low grade contamination, sterility testing may be desirable for long term
2236 cultures. These may also be sustained as separate replicate sets of flasks to provide backup cultures in case of
2237 contamination.

2238 Where primary cell cultures and tissues are used, variation in properties between individual
2239 donors must be considered, and each new batch should be qualified or controlled for key
2240 functionality. Special care should be taken to note any unusual observations in case of
2241 contamination or viral cytopathic effects or transformation, and all primary cell cultures
2242 should ideally be cryopreserved and screened for mycoplasma. Human and animal tissues
2243 and primary cells used for testing will also need to be appropriately documented, As part of
2244 QC for tissues, their differentiated state should also be documented, which may require a
2245 range of assays, which may include, morphology, histochemistry, cell markers, specific tissue
2246 function and cell-cell/matrix interactions (Stacey and Hartung, 2006). For primary cells
2247 prepared from tissues stored as banks of cryopreserved vials of cells, similar QC approaches
2248 can be used as adopted for banks of continuous cell lines (see [Table 4](#)).

Formatted: Font: Times New Roman,
12 pt

2249 **5.8 Biomarkers and functional tests to confirm the required cell function state**

2250 It is important to recognise that cell quality can vary during passaging, and in particular the
2251 time point in the growth curve at which cells are harvested may affect performance. In
2252 general, cells are best harvested in the logarithmic phase of growth. Accordingly, each culture
2253 used to set up an *in vitro* method should be subject to a key control regime measuring or
2254 indicating functionality. Acceptability criteria should be defined for functional tests and
2255 biomarkers that indicate the correct cell state. These may for example include: neuronal
2256 activity, competency of biochemical transformation, response to reference bioactive
2257 compounds, response to reference items in the particular *in vitro* method the cells are to be
2258 used for etc. In this way, each culture can be controlled, and consistency in *in vitro* methods
2259 is supported. Additionally, key markers which are associated with poor performance may be
2260 identified for future improvement. For example, expression of self-renewal genes (e.g. Oct4,
2261 Nanog, Sox2) in stem cell cultures is crucial to the functionality of the cell population
2262 (further examples for stem cells are laid out in Pistollato *et al.*, 2014 and Stacey *et al.*, in
2263 press).

2264 **5.9 Special issues for microbial strains**

2265 Microbial strains are used in many *in vitro* methods, and all the generic aspects for *in vitro*
2266 methods given here are also applicable. However, there may be special requirements needed
2267 to be applied to assure selection of certain mutant or recombinant strains (as for some
2268 recombinant cell lines). Moreover, functional tests (biochemical, cell biological, etc.) may
2269 need to be necessary to assure correct identity of each species and/or strain. Preservation of
2270 cultures is usually achieved by freeze-drying, for which storage conditions and their stability
2271 may vary for different organisms. Some strains may require cryopreservation and/or more
2272 careful testing for genetic integrity (including maintenance of plasmids) and relevant function
2273 (e.g. sensitivity to UV-light or antibiotic resistance) on recovery.

2274 **5.10 Qualification of reference strains**

2275 Where a common cell line or organism is recommended for a particular *in vitro* method it
2276 may be available from numerous sources. The selected cell line should be qualified for use in
2277 a way which can assure consistent function in the method. The idea of reference cell banks to
2278 act as central sources for all users of cell substrates has been developed by WHO (WHO,
2279 2010) and formal or regulatory *in vitro* methods may quote catalogue references from culture

2280 collections for both microorganisms and cell lines. However, culture collections may not
2281 necessarily check the performance of such strains using reference testing methods and as
2282 such, this still will need to be conducted by the *in vitro* method developer or user. If an
2283 original and quality controlled source of a cell line is not available, it is wise to obtain cells
2284 from more than one source in order to compare their performance and authenticity, as cross
2285 contamination of cell lines is very common. Moreover, cell lines with different histories of
2286 use in test facilities have shown to result in different phenotypic characteristics.

DRAFT

2287 **6 Test and reference items**

2288 This chapter describes the characterisation and preparation of the test item (e.g. a chemical to
2289 be characterized for potential hazard) and of relevant reference items for the *in vitro* test.
2290 Furthermore, details will be given as to how test item stability must be monitored, how the *in*
2291 *vitro* test environment can affect the test item and how the biokinetics of the test item in the
2292 *in vitro* method should thus be assessed. These aspects are important to ensure reproducibility
2293 among laboratories and certainty that the outcome of the *in vitro* test is indeed related to the
2294 test item. Additionally, it is described which control and additional items should be applied in
2295 general to verify correct function of the *in vitro* test, such as negative and positive controls
2296 (OECD, 2004c).

2297 A distinction is made between aspects that are important in the development phase of an *in*
2298 *vitro* method, and those that matter when the *in vitro* method is used routinely for regulatory
2299 purposes. For instance, the use of test and control items may differ between test development
2300 and routine use.

2301 The stability of the test item (or its representative preparation) and whether it is in contact
2302 with the reactive part of the test system in a known or assessable concentration should be
2303 determined when applying the *in vitro* method, in order to collect robust data from the *in*
2304 *vitro* method.

2305 **6.1 Test item**

2306 **6.1.1 Considerations during the development of the method**

2307 Test items are not used in the phase of *in vitro* method development, i.e. before the method
2308 performance has been validated. During the phase of test development, chemicals or products
2309 with well-known characteristics are evaluated to assess the relevance of the method and to
2310 amass results which will be used to set the acceptance criteria (see 6.4). Nevertheless, during
2311 the development phase, it is important to determine:

- 2312 • The different natures of test items for which the method is intended (define if the *in*
2313 *vitro* method can be used for liquids, solids, certain powders, mixtures, multi-
2314 component test chemicals, certain preparations, suspensions, nanoparticles,
2315 emulsions, etc. (OECD, 2000)). Moreover, the process of preparation of these test
2316 items should be explored.
- 2317 • The process of preparation of the test item, if needed before the test.
- 2318 • Which chemicals are suitable as reference and control items.

2319 **6.1.1.1 Nature of test items for which the method is suitable**

2320 The nature of the substances for which the developed *in vitro* method is suitable for may be
2321 characterised, but not necessarily limited to, by describing it by using a number of different
2322 properties (see below). The list given below is not exhaustive and may need to be developed
2323 depending on the nature of the test. It is important to note that solubility is a highly important
2324 yet often neglected characteristic and is therefore described in more detail in a separate
2325 section (6.1.1.2):

- 2326 • State: solid, liquid, gas and all in between-states such as aerosol, dust, or viscous liquid
2327 (see OECD TG 114 for determination of the viscosity of liquids); depending on its state,
2328 the substance could require preparation steps before the test (see also [6.1.1.36-1.1.3](#)) or a
2329 specific administration mode in the method, such as dry dispersion with pressurised air,
2330 nebulisation of a liquid formulation, or spark generation.

Formatted: Font: Times New Roman,
12 pt

- 2331 • Appearance: nominal size, morphology, size distribution, aggregation and agglomeration
 2332 phenomena and surface characteristics (surface area, surface charge, surface chemistry)
 2333 are essential characteristics to know the nature of a certain nanomaterial (OECD series
 2334 on the Safety of Manufactured Nanomaterials, n°36).
- 2335 • Colour: some test items interfere due to this characteristic with the endpoint detection
 2336 method.
- 2337 • Physicochemical characteristics:
- 2338 - pH for test item in solution (see OECD TG 122 for pH determination) and pK_a (pKa
 2339 indicates to what extent the test item may become ionised at the pH of the test
 2340 system). Changes of pH can also affect the test item in other ways than its ionisation
 2341 (see OECD TG 111 for sensitivity of hydrolysis to the pH, for example).
 - 2342 - Osmolality
 - 2343 - Volatility
 - 2344 - Solubility: see 6.1.1.2
 - 2345 - Dissociation constants in water (see OECD TG 112): dissociation is the reversible
 2346 splitting into two or more species which may be ionic. The dissociation governs the
 2347 form of the substance in the test system, which in turn determines its behaviour and
 2348 transport and which thus may affect the adsorption of the substance to culture dishes
 2349 or the penetration into cells or adsorption onto proteins in solution or resulting in a
 2350 aspecific aggregation behaviour.
 - 2351 - Lipophilicity (see OECD TG 123 and 107 for determination of the partition
 2352 coefficient = K_{ow}).
 - 2353 - Homogeneity and conditions of stable homogeneity
 - 2354 - Sensitivity to photolysis (OECD TG 316 is meant for environmental fate, thus less
 2355 suitable for this purpose, but can be used as a basis)
 - 2356 • Composition and purity: chemical purity/contaminants, microbiological
 2357 contaminants (including e.g. cell walls of decomposed microorganisms),
 2358 biological purity (e.g. of cells lines or test microorganisms, or complex
 2359 protein mixtures (vaccines)), composition of complexes (vegetal extracts,
 2360 products of fermentation, etc.). In case of a mixed solution, the list of
 2361 ingredients with percentages of each component can be relevant to
 2362 describe the composition. For each component, information like molecular
 2363 weight, chemical formula, CAS number, etc., is useful. Complex
 2364 substances could require different information. For example, substances of
 2365 Unknown or Variable composition, Complex reaction products or
 2366 Biological materials (UVCBs) cannot be sufficiently identified by their
 2367 chemical composition, because the number of constituents is relatively
 2368 large and/or the composition is, to a significant part, unknown and/or the
 2369 variability of composition is relatively large or poorly predictable. The
 2370 composition could then be defined by the manufacturing process
 2371 description¹⁹.
 - 2372 • Conditions of stability: the limits of temperature, pressure, and humidity to maintain
 2373 stability of the test item (to be compared with the *in vitro* method conditions).
 - 2374 • Microbiological status: requiring aseptic conditions or not.
- 2375
 2376 In general, the limits on test item suitability are determined so that reproducible and definable
 2377 interactions between test item and test system can be guaranteed (see chapter 6.2).

¹⁹ https://echa.europa.eu/documents/10162/13587/10_sb_siduvcb_d1_lrws_20120203_en.pdf

2378 6.1.1.2 Solubility

2379 The solubility defines how much of the test item (molecular and ionized forms) that can be
2380 maximally dispersed in the solvent to be used for the test (=thermodynamic solubility). The
2381 rate of dissolution (in practical terms: fraction of dissolved test item after a fixed time) is
2382 called kinetic solubility (Jouyban and Fakhree, 2012). This section is focused on
2383 thermodynamic solubility.

2384 The reason why solubility is an important property to know of a test item is that, above this
2385 limit, the substance precipitates and the effective concentration in the test medium is lower
2386 than the nominal concentration. Precipitates may also affect read-outs of the *in vitro* method
2387 and lead to impaired reproducibility within and between labs. It is thus important to ensure
2388 that all test items are properly dissolved for use in the *in vitro* method. The highest test
2389 concentration has to be below the solubility limit.

2390 Different ways to determine the solubility of a substance in a defined solvent are available.
2391 While computational methods for solubility predictions in different types of solvent or
2392 matrices are available (Bergström et al., 2002; Persson et al., 2013), it remains difficult to
2393 predict the solubility in the medium used for the *in vitro* method. Most frequently, solubility
2394 of a compound is determined visually, by determining whether a clear solution is formed
2395 after adding the substance of interest. However, while the visual test is very simple to apply,
2396 it gives only a rough impression of the solubility. In most cases however, it is sufficient for
2397 simple checking of solubility where reliability can be enhanced by use of microscopy to
2398 detect solid particulates or liquid droplet suspension (indicative of insolubility). However, it
2399 is a rather subjective operator-dependent judgement. Reliability can also be improved in a
2400 relatively easy way by centrifugation, particularly for detection of precipitation in medium
2401 dilutions, where foaming may obscure visual observation. Note that solubility is affected by
2402 the composition of the substance (presence of impurities) and by the test conditions
2403 (temperature, incubation time, possible adsorption to the test vessel or to medium
2404 constituents (e.g. albumin)). OECD TG 105 can be used for the determination of the aqueous
2405 solubility of pure substances which are stable in water and not volatile, and OECD TG 116
2406 can be used for fat solubility determination (fat solubility is the mass fraction of substance
2407 which forms a homogeneous phase with a liquid fat (oil) without giving rise to chemical
2408 reactions). Nephelometric, UV-spectroscopic and HPLC methods can also be used to
2409 determine solubility (Hoelke *et al.*, 2009).

2410 Nephelometry facilitates solubility determination, particularly suited to serial measurement
2411 (e.g., ranges of chemicals and/or concentrations) allowing systematic and precise evaluation
2412 of turbidity due to dispersed precipitation, independent of matrix composition. However, the
2413 measurement is relative, requiring a definition of threshold turbidity for insolubility based on
2414 expedient practice with detection limit dependent on instrument sensitivity. Moreover, even
2415 nephelometry may not detect chemicals such as transparent immiscible liquids for which
2416 visual inspection, enhanced by experienced microscope observation, remains a reliable
2417 approach.

2418 HPLC and UV spectrophotometry provides a quantitative determination of concentration
2419 with the use standard curves. While both methods are valid for solutions prepared in solvent,
2420 they may not be valid for preparations in biological media, which contain many components
2421 that often interfere with solubility detection. Cell culture media cannot be injected into HPLC
2422 columns and their multiple components will likely obscure the compound of interest through
2423 their inherent UV absorbance over a range of wavelengths. This necessitates pre-purification
2424 and extraction steps for quantifying test item concentration in media via HPLC.

2425 Test item chemicals are generally dissolved in solvent (e.g., DMSO, Ethanol) to create a
2426 stock solution at a predetermined target concentration (e.g., 50mg/mL or 100mM). The test
2427 item should have a relatively high solubility in the solvent of choice and the solvent should
2428 not interfere with the test item (e.g. inactivate the compound). For example, the commonly
2429 used solvent DMSO can reduce the effects of the platinum complexes (Hall *et al.*, 2014). In
2430 addition, the solvent should not affect cell health or the phenotype of the cells in the assay
2431 when diluted in media. Furthermore if the stock solution is diluted in media up to a
2432 concentration exceeding the solubility of the substance, the test item may precipitate.

2433 Regarding the sample preparation procedure, the following issues are key for reproducible
2434 results:

2435 - Optimal time for dissolution in solvent: Does the drug dissolve immediately (1 min
2436 vortex) in the solvent or does it require additional treatment (sonication and warming)
2437 and time (e.g., 2 hours or overnight)?

2438 - Solubility in media upon dilution: Including incubation to mimic assay conditions may
2439 be relevant: i.e., monitoring stability on incubation (e.g. at assay temperature and in the
2440 presence of CO₂) over a time period (e.g. 24 hours).

2441 - Visual inspection of solubility: Unless sample material is expensive or available
2442 quantities are limited, stock solutions for visual inspection of solubility should be
2443 prepared with a minimum weight of 25mg, and in a minimum volume of 0.5 ml.

2444

2445 The solubility in assay medium may be higher than in pure water, due to adsorption to
2446 medium contents such as proteins. In addition, the higher temperature applied in the *in vitro*
2447 method (37°C) than in OECD TG 105 (20°C) will give rise to a higher solubility in the assay.
2448 However, assay media typically have a rather high ionic strength and an inherently a complex
2449 composition. This makes it difficult to predict test item solubility upon dilution from stock
2450 solution in various media solvents. Therefore, it is also necessary to determine the solubility
2451 of the final test concentrations in the assay medium under assay conditions. In case of
2452 inorganic substances, the anion and cation part of the test item may precipitate with other
2453 cations and anions present in the culture medium, if the solubility of these newly combined
2454 salts is exceeded. It is therefore recommended to visually monitor the test system for
2455 precipitation, and to verify whether the anions and cations present in the medium can form
2456 low-solubility salts with the test item.

2457 Test items insoluble within an acceptable concentration range should be considered
2458 incompatible with the *in vitro* method. A remedy may be to increase the solvent
2459 concentration in the *in vitro* method. However, the tolerable solvent concentrations will
2460 depend on the solvent and the test system used: As a general rule, the final solvent
2461 concentrations should be as low as possible to avoid any potential interference with the *in*
2462 *vitro* method. The concentration of solvent present in the final test concentration should be
2463 considered in the test controls.

2464 As for nanomaterials, special issues on measuring solubility and dispersion characteristics
2465 may arise. For these materials, the specific guidance documents are best followed, which are
2466 continuously being developed (Scenhir, 2015). Any toxicity testing using *in vitro* methods
2467 should pay special attention to the agglomeration/aggregation behaviour, and the insoluble/
2468 partially-soluble nature of nanomaterials (Scenhir, 2015). Possibilities for dis-agglomeration
2469 and re-aggregation of nanomaterials should also be considered: some properties of
2470 nanomaterials may change due to interaction with the surrounding media.

2471 **6.1.1.3 Test item preparation**

2472 Test items may have to go through various steps of preparation, such as dissolution, dilution,
2473 extraction by wetting or centrifugation before being suitable for use in the *in vitro* method.
2474 These steps must preserve the characteristics of the test item. The purpose of each step of the
2475 preparation has to be explained, and the critical limits step/procedure should be determined.
2476 The impact on the test item stability, homogeneity and integrity should be assessed.

2477 Once prepared as a dilution in a solvent, chemical stability in solution is crucial, since it
2478 determines the concentration of the work solution to be applied to the test system. The
2479 concentration in stock solutions and exposure medium has to be measured over time, to check
2480 for stability and avoid aberration in cell treatment. Relevant examples of compounds
2481 hydrolysed in aqueous solutions can be found in (Crean *et al.*, 2015; Pomponio *et al.*, 2015).

2482 For more complex test items there are existing guidelines to aid this process: e.g. ISO 10993-
2483 12 gives conditions of the extraction to obtain a representative extract of medical devices
2484 depending on their composition, and The OECD series on the Safety of Manufactured
2485 Nanomaterials, n°36, gives advice on how to prepare and characterise a nanomaterial
2486 solution.

2487 **6.1.2 Considerations for the final user of the validated method**

2488 The type of test item for which the developed *in vitro* method is suitable has to be defined in
2489 the description of the validated method. Test items at the user stage should be characterised
2490 as required by the GLP principles (see OECD n°1, chapter 6.2): “each test and reference item
2491 should be appropriately identified (...); the identity, including batch number, purity,
2492 composition, concentrations, or other characteristics to appropriately define each batch of the
2493 test or reference items should be known (...).The stability of test and reference items under
2494 storage and test conditions should be known (...). If the test item is administered or applied in
2495 a vehicle, the homogeneity, concentration and stability of the test item in that vehicle should
2496 be determined (...).”

2497 Interactions with the test system should be predicted and a way of assessment of such
2498 interactions should be implemented (cf section 6.2). Important examples are adsorption to
2499 protein, precipitation with carbonate in the medium and hydrolysis by aqueous environment.

2500 **6.2 Interaction between test item and test system**

2501 A two-way interaction has to be assumed between the test item and the test system. In one
2502 respect, the test system can affect the test item (in analogy to pharmacokinetics in *in vivo*
2503 models; detailed in 6.4). In another respect, the test item can affect the test system in specific
2504 ways (alteration of a readout, in accordance with the design and intended application of the
2505 test system; (see Sections 5 and 8) or in non-intended ways, by interfering with the overall
2506 performance of the biological model on which the test system is based, or by disturbing a
2507 correct readout of the test endpoint. There are endless possibilities for artefacts to be created
2508 in this way. As not all of these can be controlled automatically, experienced operators and
2509 personnel interpreting the test data are required to detect potential problems. Problem
2510 detection is also facilitated by regular inclusion of consistency controls and plausibility
2511 considerations (e.g. do compounds with similar structure or similar mode of action behave
2512 similarly?; can effects be reversed; does another test system for the same biological process
2513 give similar results?; are findings consistent with biological expectations concerning
2514 concentration and timing of effect?; etc.). For more detailed discussion, three elements of an
2515 *in vitro* test will be considered separately (Aschner *et al.*, 2016; Schmidt *et al.*, 2016): the test

2516 system (biological model), the test endpoint, and the analytical endpoint (Leist *et al.*, 2010).
2517 These issues are relevant both in the development stage of an *in vitro* method and in the
2518 application stage.

2519

2520 **6.2.1 Interference with the test system**

2521 Test items can disturb the test system, especially if it is based on living cells, as they are
2522 highly responsive to changes in their environment. The most frequent and serious disturbance
2523 is general cytotoxicity often leading to cell death.

2524

2525 **6.2.1.1 Cytotoxicity**

2526 If the test endpoint is not cytotoxicity, then cytotoxicity triggered by test items is a serious
2527 confounder and needs to be controlled for. Indeed, changing cell numbers *in vitro* is known to
2528 affect observed effect concentrations (Gülden *et al.*, 2015, 2001). This is particularly
2529 important in repeated treatments (Kramer *et al.*, 2015).

2530 Measurement of cytotoxicity should be done using the same conditions as used for the
2531 specific test endpoint (i.e. in identical samples, ideally during the same run, or even better on
2532 the same plate), so as to obtain reliable and relevant cytotoxicity data. Alternative approaches
2533 use measurements in parallel cultures. Viability controls in related, but not identical, culture
2534 conditions (different plate format, different cell preparation, etc.) should be avoided. The
2535 definition of cytotoxicity is not trivial. Often 100% cytotoxicity is referred to as 100% cell
2536 death, however, for a cytotoxicity assay with a metabolic activity endpoint, a 100%
2537 cytotoxicity just means that the metabolic activity is 0% compared to the healthy control
2538 cells. Therefore, the choice of method used for cytotoxicity determination, but also the
2539 interpretation of the results needs careful consideration (see [Table 5](#)~~Table 5~~). A single
2540 endpoint is usually not sufficient to be fully conclusive. A combination of cell counting and a
2541 population measure (e.g. resazurin reduction), or a combination of a viability measurement
2542 (calcein staining; dye exclusion, neutral red uptake) and a cell death measurement (propidium
2543 iodide uptake; LDH-release; annexin V staining) provides a greater level of certainty.
2544 Importantly, positive and negative controls for the viability assay have to be included and
2545 need to be considered for normalisation of viability data.

2546 An important point to consider is the timing of cell death. For instance, a test item may
2547 trigger definitive cytotoxicity (e.g. > 70% of the cell population dying within 16 h). If the
2548 specific test endpoint (and cytotoxicity) is measured after e.g. 4 h, then cytotoxicity will not
2549 be measurable by the method of choice in the test, although the cascade of cell death may
2550 have been triggered. There are no established rules on how to deal with this (relatively
2551 frequent) situation. One solution is to follow up on test results by alternative tests for the
2552 same endpoint, or by using the same test with a changed incubation scheme (e.g. prolonged
2553 incubation). This is particularly important, if data are used for risk assessment and far-
2554 reaching regulatory decisions.

2555

2556 **Table 5: Viability testing of cell cultures**

Endpoint	Assay	Mechanism and comments
1. Structural cell damage (non-invasive)	Evaluation of overall cell shape, cytoplasmic structure, flatness and outline	Screening assay that covers many forms of damage with high sensitivity, if observer is experienced. May be automated and rendered quantitative to some extent by high

Formatted: Font: Times New Roman, 12 pt, Font color: Text 1

	properties on a good phase contrast light microscope	content imaging. Advantages: high throughput, non-invasive, repeatable on same well over time. Disadvantages: No clear prediction model (only qualitative data, no exact cell death definition), no stand alone approach; requires extensive experience of operator.
	LDH-release test	Cells with intact membrane retain their content of LDH enzyme; LDH is released when cell membranes rupture (non-viable cells), and the enzyme can then be measured in the supernatant. To give fully quantitative data, the assay requires normalisation to the total LDH content of a culture well. It can to some extent be repeated for the same culture at different time points. Advantages: Measurement of a definite/unambiguous cell death endpoint; can be combined with cell function assays. Allows cells to be used for other purposes, if only supernatant is sampled. Disadvantages: Information only for cell populations. Normalisation necessary (extra wells for controls). Frequently high background LDH levels are observed (e.g. from serum components; signal/noise ratio can be bad in some culture media or with some cell types). Problems with long-term assays involving medium changes.
2. Structural cell damage (invasive)	Membrane penetration by dyes to detect 'cytotoxicity' (e.g. naphthalene black, trypan blue, propidium iodide, ethidium bromide, EH-1)	Dyes are selected so that they stain non-viable cells, but do not enter viable cells with an intact cell membrane. Some of the dyes stain the entire cell (e.g. trypan blue), others stain the nucleus/DNA (e.g. propidium iodide). Dyes that only stain dead cells usually need a combination with a method that stains/identifies all cells (such as phase contrast for trypan blue, or a nuclear counterstain (H-33342, acridine orange, SYTO-13) for fluorescent dyes. Advantages: Rapid and usually easy to interpret. Gives information on the single cell level. High throughput and absolute quantification are possible (high content imaging). Disadvantages: May overestimate viability since apoptotic cells continue to have intact membranes and may appear viable. Some dyes (e.g. trypan blue, H-33342) are cytotoxic, so that the evaluation has to be performed rapidly.
	Retention of dyes within intact cells to detect 'viability' (e.g. fluorescein diacetate or calcein-AM)	After dye exposure, viable cells fluoresce when observed under UV light. The lipid-soluble dyes are transformed by cellular enzymes (esterases) into lipid-insoluble fluorescent compounds that cannot escape from cells with intact membranes. Thus, cells can be observed under a microscope (single cell analysis) or with a fluorescent plate reader (population analysis). The dyes are often used in combination with a cytotoxicity stain (e.g. propidium iodide). Advantages: Rapid and usually easy to interpret. Gives

		<p>information on the single cell level (including morphological information on the cell shape). High throughput and absolute quantification are possible (high content imaging, fluorescent plate reader or FACS).</p> <p>Disadvantages: Some cells leak the dyes; some cells actively export the dyes through P-gp activity. Many fluorescent dyes are prone to photo-bleaching, and some may be sensitive to their local environment (pH etc).</p>
	Evaluation of programmed cell death/apoptosis markers	<p>As programmed cell death is a universal cell biological process based on defined biochemical pathways and organelle changes, the activation of cell-death-associated pathways is often used as surrogate marker for cell death. An example for such a pathway is the activation of caspases (detectable in populations by enzymatic analysis or in single cells by staining) or the activation of endonucleases (detectable on population level as DNA-fragmentation). Moreover, a typical type of chromatin condensation (detectable by DNA stains) and the display of phosphatidylserine on the outside of the plasma membrane (detectable by annexin staining) is highly correlated with apoptotic death.</p> <p>Advantages: Adds mechanistic information to cytotoxicity data. Several endpoints are easy to quantify and useful for high through-put measurements.</p> <p>Disadvantages: Not all types of cell death may be detected by a given endpoint. Needs to be combined with a general cytotoxicity test. Some endpoints are prone to artefacts (annexin staining) and some staining techniques (TUNEL, caspase-3) lead to an un-intentional selection of subpopulations. Caspase activity measurement does not easily yield a prediction model for the extent of cell death.</p>
3. Cell growth	Cell counting	<p>For some cell populations, continued growth is a defining feature, and thus impaired growth needs to be considered as a reduction of viability. Notably, impaired growth/proliferation is not necessarily correlated with cell death; it is thus rather a functional viability endpoint than a cytotoxicity measure. A special case for growth is the increase in cell size without proliferation. This feature is e.g. seen for the extension of neurites by neurons. The gold standard analytical endpoint for the growth/proliferation endpoint is counting (or morphometry). There are many ways of counting cells, either as single particles (e.g. by FACS or HCI) or by assessing a biochemical parameter correlated to cell number (e.g. DNA content).</p> <p>Advantages: growth can be a sensitive parameter of cell well-being.</p> <p>Disadvantages: growth is a functional endpoint, not necessarily linked to cytotoxicity; artefacts for growth endpoints may arise from inhomogeneous growth of subpopulations: moreover, growth may hide ongoing cell</p>

		death, and thus needs careful control in combination with cytotoxicity assays.
	BrdU or EdU incorporation	<p>Measures new DNA synthesis based on incorporation of the easily detectable nucleoside analogs BrdU (or EdU) into DNA. BrdU can be detected e.g. by fluorescent-labelled antibodies in permeabilised cells. Alternatively, radiolabelled thymidine can be used.</p> <p>Advantages: Measurement on single cell level. Easy to quantify and use at high throughput.</p> <p>Disadvantages: BrdU/EdU can be cytotoxic; no information available on how often one given cell has divided. High cost and effort compared to counting.</p>
	Staining of cellular components that are proportional to overall cell mass (proteins by e.g. sulforhodamine B or crystal violet; DNA by Hoechst H-33342)	<p>These assays evaluate a surrogate measure of overall cell mass and assume that it correlates with total cell number. In non-proliferating cells, or with continuous ongoing proliferation, the endpoints are also frequently used as indicators of cytotoxicity, as dead cells often detach from plates and reduce the overall cell mass.</p> <p>Advantages: Simple and cheap; lots of historical data</p> <p>Disadvantages: Mostly not a single cell measure but only population level. Protein staining is only a surrogate endpoint of real cell number. For DNA quantification with Hoechst 33342: fluorescent probe penetration, bleaching, and cytotoxicity are issues to be considered.</p>
4. Cellular metabolism	3-(4,5-dimethylthiazol-2-yl)-2,5-diphenyltetrazolium (MTT) assay, or similar tetrazolium dye reduction assays from multiple suppliers	<p>Biochemical activity (mostly mitochondrial metabolism; production of reducing equivalents like NAD(P)H) in viable cells causes reduction of the tetrazolium dye. The resultant formazan is extracted and measured spectrophotometrically. The rate of formation of formazan corresponds to the function of essential cellular processes like respiration.</p> <p>Advantages: High throughput, easy, robust, low cost. Used in several ISO standards and OECD test guidelines. High sensitivity. Can be used for tissue constructs.</p> <p>Disadvantages: Measures amount of viable cells (only indirect measure of cell death), and needs control for contribution of proliferation. Cells with reduced mitochondrial function may appear non-viable. Inhibition of cell metabolism by the test item causes low values in the assay which is not necessarily related to cell viability. Some test items interfere with the assay e.g. by reducing the dye why interference testing is recommended. Measurement usually not on single cell level. Some cell cultures need long time to reduce sufficient amount of dye (no sharp time point for viability definition). Assessment of kinetic of the reduction may be necessary to ensure proper selection of incubation time with a tetrazolium dye (to avoid reaching plateau of OD).</p>

	Resazurin reduction assay (sometimes called Alamar blue)	<p>Similar to tetrazolium reduction assays. Fluorescent/absorbent resorufin is formed from resazurin through mitochondrial metabolism of viable cells.</p> <p>Advantages: Many tests can be performed rapidly in multi-well dishes. Cells can be tested repeatedly (non-invasive measurement). High sensitivity.</p> <p>Disadvantages: Cells with reduced mitochondrial function may appear non-viable. Some test items interfere with the assay (e.g. superoxide also reduces the dye) why interference testing is recommended. Measurement only on population level. Some cell cultures need a long time to reduce sufficient resazurin (no sharp time point for viability definition).</p>
	Mitochondrial depolarisation assays (based on fluorescent indicator dyes)	<p>Many organelle functions are used as endpoints of general cell health. Most frequently, mitochondrial function is assessed (see MTT, resazurin). One mitochondrial test on the single cell level is the measurement of mitochondrial membrane potential by addition of potential sensing fluorescent dyes like JC-1, TMRE, MitoTracker, etc. Quantification is by HCl or FACS</p> <p>Advantages: fast, cheap, high throughput; single cell information.</p> <p>Disadvantages: as for MTT (measures cell function, not cytotoxicity). Artefacts by test items that affect mitochondria specifically. Artefacts by test items that affect plasma membrane potential. Artefacts due to bleaching, quenching and unquenching, and due to shape changes and clustering of mitochondria.</p>
	Neutral red assay (ISO 10993)	<p>A cell organelle function assay assessing lysosomal function. Active cells accumulate the red dye in lysosomes and the dye incorporation is measured by spectrophotometric analysis.</p> <p>Advantages: Low cost. Used in several ISO standards and OECD test guidelines. Historic data base.</p> <p>Disadvantages: Normalisation required for quantitative measurement, e.g. with protein content or number of cells. Gives usually information only at the population level. Not suited for tissue constructs and certain cell lines.</p>
	ATP assays	<p>Measurement of the total ATP content in a cell population. Dying cells fail to produce ATP, have an increased ATP consumption, and may lose ATP through perforations of the plasma membrane. For the test, cell lysates are prepared, and the ATP content is assessed by a luminometric assay.</p> <p>Advantages: fast, high throughput</p> <p>Disadvantages: no single cell data, expensive, requires a luminometer, as MTT: measurement of viable cell mass, not a direct measure of cytotoxicity. Artefacts as for other mitochondrial tests.</p>

2557

2558 **6.2.1.2 Functional disturbances**

2559 While the strict definition of cytotoxicity refers to cell death, a wider interpretation also
2560 includes adverse effects on cells that alter their functionality but do not lead to cell death
2561 (within the observation period). For instance, protein synthesis may be impaired, or
2562 mitochondrial function altered. Cytostasis, where dividing cells do not die but cease dividing,
2563 is another example of delayed cell death which can impact the endpoint measures. This can
2564 affect the specific endpoint of a test system (e.g. luciferase reporter assay), without being
2565 relevant for the intended test objective. Test items with such properties can lead to erroneous
2566 readouts.

2567 **6.2.1.3 Disturbed differentiation state or gene/protein expression**

2568 A special case of artefacts caused by test items is the change of biological properties of the
2569 test without overt cytotoxicity. The most common example is an altered differentiation of
2570 cells or an altered composition of cell sub-populations. For instance, a test item might alter
2571 cell differentiation state in a migration assay, and this alteration might lead to altered
2572 migration. The item would be wrongly classified as modifying cell migration. Another
2573 example would be measurement of monocyte function (e.g. cytokine release) in a whole
2574 blood assay. If a test item leads to platelet degranulation, it might influence the overall
2575 endpoint of the *in vitro* method without affecting the monocyte response as such.

2576 **6.2.1.4 Altered communication/adhesion properties**

2577 Another special case of artefacts can be generated by interference of the test item with cell
2578 adhesion or communication. This is listed here separately, as it would not normally be
2579 detected by cytotoxicity assays, but it would strongly alter the behaviour of the test system
2580 (biological model) in the test situation. An example is binding of test item to molecules used
2581 for the coating of culture dishes. This would then alter readouts such as migration or neurite
2582 growth, without really affecting such processes within the cells (and without necessarily
2583 being relevant *in vivo*).

2584 **6.2.2 Interference with *in vitro* method endpoint**

2585 For pure test items, most of the unwanted interactions with the *in vitro* method are covered by
2586 considering undesired interactions with either the test system or the *in vitro* method endpoint.
2587 The majority of interferences with the *in vitro* method endpoint will be related to cytotoxicity
2588 (immediate or delayed cell death or functional impairment), as covered in 6.2.1.1. For test
2589 items containing impurities or non-inert additional substances in their formulation (see 6.2.4)
2590 the situation can become more complex: this highlights the need to have clear specification
2591 for the test item, as stated above. For instance, impurities (e.g. detergents or solvents) may
2592 alter skin or blood-brain-barrier (BBB) permeability (without being cytotoxic) and thus result
2593 in incorrect data on the pure test item of interest, if the *in vitro* test assesses e.g. skin or BBB
2594 permeation capacity. In other cases the test item is a finished product. Then potential
2595 impurities and contaminations are part of the product and their effect on the response is
2596 important and has to be evaluated.

2597 **6.2.3 Interference with the analytical endpoint**

2598 Interference of the test item with the *in vitro* method endpoint means that the test item
2599 disturbs the normal measurement results. This can be controlled for by performing the *in*
2600 *vitro* method using adequate positive, negative, blank or vehicle controls. If the endpoints are

2601 of analytical nature, the controls can also be spiked with test item to verify that the test item
2602 do not in any way hinder the normal function of the test system or interfere with the readout.

2603
2604 Examples of such kind of interference include:

2605 **Fluorescence/absorbance-based methods:** disturbance by test items that fluoresce or absorb
2606 at the evaluation wavelength, or test items that quench fluorescence.

2607 **Enzymatic assays:** alteration of enzyme function, of co-factor, or of other limiting reagents
2608 by test item; display of enzymatic activity (or chemical reactivity) by test item itself.

2609 **Resazurin/ or MTT reduction:** strongly reducing agents directly reduce resazurin/ (or MTT)
2610 non-enzymatically. Compounds that trigger the release of superoxide can trigger reduction of
2611 resazurin by superoxide. This results in erroneous cytotoxicity data.

2612 Another relevant example of this kind of interference is provided by the interactions between
2613 test reagents and nanomaterials in colorimetric assays for cytotoxicity (such as
2614 sulforhodamine B dye, or MTT used in the viability assays) (Scenhir, 2015). Moreover, some
2615 nanomaterials may themselves disperse/absorb light and therefore interfere with the
2616 measurements in colorimetric assays. Some of these problems might be overcome by either
2617 adding appropriate controls or modifying existing protocols: as an example removal of
2618 nanomaterials via centrifugation before reading the assay can reduce the variations in data
2619 generated for the same nanomaterials (Scenhir, 2015).

2620 **6.2.4 Consideration of interferences not coming from the active ingredient**

2621 With test items that are not pure, the interferences with the test may come from impurities or
2622 from ingredients of the formulation. Particularly difficult cases arise when such additional
2623 chemicals are inactive alone, but synergize somehow with the effect of the test item.

2624 This can also occur for the solvent of the test item. Frequently, a solvent concentration that
2625 does not affect the standard endpoint of a test as such (e.g. 0.1% DMSO) may still alter the
2626 effect of a test item on the test system (e.g. in the case of DMSO: through the antioxidant
2627 properties of DMSO; or through its effect on cell membranes; or through other activities
2628 including cell differentiation).

2629 In case of a test item consisting of a natural mixture (e.g. essential oils), it should be
2630 considered to test the mixture as well as the known pure substances present, since the other
2631 ingredients of the mixture can change the overall effect of the test item. The different kinetics
2632 of the ingredients must then be considered; however, although ingredients that are not
2633 absorbed *in vivo* will not be able to have an effect on the test item systemic toxicity, they may
2634 affect test item toxicity *in vitro*.

2635 **6.3 Biokinetics / dose extrapolation**

2636 Just like the biokinetics *in vivo* are about what the body of the organism does to the test item,
2637 the biokinetics *in vitro* concerns what the *in vitro* test environment does to the test item. A
2638 central issue in biokinetics is that generally only the freely dissolved molecules of a chemical
2639 can pass membrane barriers and reach a target inside a cell. Thus, in an *in vitro* system, the
2640 freely dissolved concentration of the test item in the medium or in the cell (as close to the
2641 target as possible) is the central parameter. Different processes result in a freely dissolved
2642 concentration that is not the same as the nominal concentration, (i.e. the added
2643 concentration). These processes are described in e.g. (Heringa *et al.*, 2006) and (Groothuis *et*
2644 *al.*, 2015), and were one of the main topics investigated by the FP7 EU Project Predict-IV,
2645 aimed to improve the predictivity of *in vitro* assays for unwanted effects of drugs after
2646 repeated dosing integrating biokinetics and biodynamic data. As one of the project outputs, a


2647 step-wise strategy was applied to measure and model cell exposure levels over time of a
2648 selected number of drugs in the developed *in vitro* assays. The strategy and the major
2649 obtained results are described in (Kramer *et al.*, 2015).

2650 Figure 3 schematically depicts the different processes, which can affect a xenobiotic *in vitro*
2651 bioavailability. Also the processes that change the identity of the test item, affecting its
2652 stability, are described briefly in the following.

2653

2654 **Figure 4: Schematic representation of the processes that can cause the final target**
2655 **concentration to be different than the nominal concentration in an *in vitro* test (Kramer**
2656 ***et al.*, 2012)**

2657



2658

2659 6.3.1 Kinetic processes

2660 6.3.1.1 Evaporation / plastic binding / precipitation

2661 *In vitro* systems are often open, with a small gap between the well plate and the lid, to allow
2662 air circulation for provision of oxygen for the cells. This air circulation allows volatile
2663 substances to evaporate into the air of the incubator. This may decrease the concentration in
2664 the medium in the test system, but can also contaminate medium in e.g. blank wells, as the
2665 substance can dissolve from the air into the medium of other wells present in the well plate or
2666 even the incubator. An example showing the effect of evaporation on test results can be
2667 found in Tanneberger et al 2010 (Tanneberger *et al.*, 2010).

2668 Lipophilic substances tend to bind to the plastic the cell culture plates are made of, although
2669 differences exist among the types of plastic used. The adsorption to polycarbonate is limited,
2670 but in organ-on-a-chip devices made of Polydimethylsiloxane (PDMS), there will be
2671 partitioning between the PDMS and the medium. PDMS is even used as an extraction
2672 material for-solid phase microextraction (SPME) (Heringa and Hermens, 2003), it is
2673 therefore not suitable for *in vitro* test devices for testing of chemical substances. Glass is a
2674 better material to avoid binding to some extent, but very lipophilic substances are known also
2675 to bind to glass. Silanized glass can decrease this binding even further. Using glass has other

2676 practical downsides in *in vitro* tests. Examples where considerable binding to plastic was
2677 measured are the one of Kramer *et al.* (Kramer *et al.*, 2012). In this study, it is also shown
2678 how the addition of serum to medium decreases the binding to plastic. Other examples are
2679 reviewed in (Kramer *et al.*, 2015), reporting results of the Predict-IV project on cyclosporine
2680 A, amiodarone and chlorpromazine. The addition of serum to medium decreases the binding
2681 to plastic, but likely also the uptake into the cells (Pomponio *et al.*, 2015).

2682 Sorption of the test item to cell-attachment matrices (e.g. collagen or matrigel layer used with
2683 hepatocytes in culture) is a specific aspect of interaction with the test device, although the
2684 relationship between a test item's lipophilicity and binding to is not as clear cut as it is for
2685 binding to plastic labware. The possible physical sequestration of test items, can lead to
2686 overestimating intracellular concentrations (Kramer *et al.*, 2015).

2687 **6.3.1.2 Chemical degradation**

2688 The aqueous environment of the medium in an *in vitro* test enables spontaneous hydrolysis
2689 (i.e. without the aid of an enzyme) of substances with structures sensitive to this chemical
2690 reaction. During the time the test system, e.g. the well plate, is outside of the incubator, light
2691 will reach the medium and photolysis can take place for light sensitive substances. Therefore,
2692 information on hydrolysis and photolysis sensitivity is necessary before a substance is tested
2693 in an *in vitro* method (see 6.1.1.1). More generally, each test facility should have adequate
2694 test item characterisation procedures in place to identify if the test item characteristics are
2695 compatible with the *in vitro* method.

2696 **6.3.1.3 Metabolism/metabolic stability**

2697 Some cell types have metabolic capacity, meaning that they contain significant levels of
2698 enzymes that convert the test substance to another substance. Especially cells originating
2699 from liver, intestine and lung are known to possess metabolic capacity, in decreasing order.
2700 In test systems with such cells, especially from these tissues, the concentration of the test
2701 item may decrease because of this metabolism, and the concentration of metabolites will
2702 increase. When a positive hazard response is obtained in such a cell system, it may thus either
2703 be caused by the test item itself, or its metabolite(s). The time profile of the response can
2704 reveal which is the main causative agent: when there is a lag time for the response (compared
2705 to the positive control or other reference items), it could be that a metabolite is responsible
2706 for the response. A good example is described in (Pomponio *et al.*, 2015)

2707

2708 **6.3.1.4 Protein binding**

2709 Serum is often added to cell culture medium to supplement it with important factors required
2710 for cell proliferation and maintenance. Serum-free medium is available and used, but not all
2711 cell types thrive in such culture conditions. Serum contains proteins, including albumin,
2712 which has non-specific binding sites, to which most organic substances tend to bind. As
2713 proteins are large molecules that do not transfer across a membrane, the binding to a protein
2714 renders a test item unavailable for cellular uptake, thus unable to reach any target inside the
2715 cell. An example of the effect of serum protein binding can be found in Heringa *et al.*,
2716 (Heringa *et al.*, 2004) and in Pomponio *et al.*, (Pomponio *et al.*, 2015). On the other hand,
2717 serum proteins can also make some test items more accessible or more stable. E.g., for
2718 medical devices a medium with serum is preferred for extraction because of its ability to
2719 support cellular growth as well as to extract both polar and non-polar substances. In addition,
2720 protein binding also occurs *in vivo*. Gulden *et al.* 2005 therefore developed an *in vitro-in vivo*
2721 extrapolation method extrapolate nominal effective *in vitro* concentrations equivalent *in vivo*

2722 plasma concentrations by accounting for the differences in protein concentrations (Gülden
2723 and Seibert, 2003).

2724 **6.3.1.5 Cell membrane absorption**

2725 Cell membranes are composed of fatty acids, thus provide a lipid environment in which
2726 lipophilic substances will like to absorb. These absorbed molecules are then also not
2727 available for a target inside the cell. An example showing the effect of membrane sorption
2728 can be found in Gülden *et al.* (Gülden *et al.*, 2001) and in Bellwon *et al.*, (Bellwon *et al.*,
2729 2015).

2730 **6.3.2 Measurement of free concentration /passive dosing**

2731 Clearly, several processes can influence how much of the added test substance will actually
2732 reach the target and this relates to its saturation concentration. If test results are based on the
2733 added, or nominal, concentrations, considerable variation between laboratories may be
2734 obtained. Furthermore, *in vitro* processes are included in these results (e.g. an EC₅₀),
2735 rendering these unfit for extrapolation to *in vivo* (see IVIVE below and (Kramer *et al.*,
2736 2015)). For example, if there is considerable evaporation, the EC₅₀ *in vitro* will appear to be
2737 much higher than it will be in the same tissue *in vivo*, as there is no evaporation *in vivo*.
2738 Thus, in order to obtain pure EC₅₀ values, that relate target concentrations to responses, these
2739 target concentrations should be measured. As the precise concentration at the target site
2740 inside the cell is too difficult to measure, the best approximation should be measured: the free
2741 concentration in the cell or in the medium. The free concentration in the cell is often still
2742 difficult to measure, therefore the free concentration in the medium (similar to the cellular
2743 concentration for membrane permeable substances), or the total concentration in the cells
2744 (often for metals) are usually measured. Further information can then be added by
2745 calculations that take physicochemical and biochemical properties (e.g. transporter
2746 substrates) of the substances into account.


2747 Methods with which the free concentration can be measured have been reviewed in Heringa
2748 *et al.* 2003 (Heringa and Hermens, 2003). This review also describes how negligible
2749 depletion-solid phase extraction (negligible depletion SPME) should be applied to measure
2750 free concentrations. This method is very suitable for *in vitro* tests, as it is suitable for small
2751 volumes. Examples of its application in *in vitro* tests are (Heringa *et al.*, 2004), (Broeders *et*
2752 *al.*, 2011), and (Kramer *et al.*, 2012).

2753 Measuring the free concentrations does require extra effort and resources in the conduct of
2754 the *in vitro* test, as e.g. a chemical analysis method is necessary. This effort can be saved in
2755 some instances, depending on the properties of the test item: in case of very hydrophilic, non-
2756 volatile substances that hardly bind to serum proteins, there will hardly be any losses and the
2757 nominal concentration will be very similar to the free concentration. Groothuis *et al.*
2758 (Groothuis *et al.*, 2015) provide a decision scheme on which concentration should/can be
2759 used as dose metric.

2760

2761
2762

Figure 5: Flow chart to aid in choosing an appropriate dose metric for a specific *in vitro* toxicity test (Groothuis *et al.*, 2015)



2763

2764
2765
2766
2767
2768
2769
2770

First, a choice should be made for dose type based on the characteristics of the chemical and available knowledge. Then, the metric can be integrated or averaged in case of time-dependent exposure and irreversible mechanisms, or steady reduction over time. Peak concentration is defined here as the maximum concentration reached during the exposure period. BK/TD may be applied to model partitioning and assess concentration changes over time. The chart has been compiled by (Groothuis *et al.*, 2015) using literature data (Austin *et al.*, 2002; Glden *et al.*, 2010; Glden and Seibert, 2003; Knbel *et al.*, 2012; OECD, 2011, 2006a, 2006b; Reinert *et al.*, 2002; Riedl and Altenburger, 2007).

2771

2772
2773
2774
2775
2776
2777
2778
2779
2780
2781

To avoid the effort of measuring free concentrations in every sample, passive dosing can be applied. In this method, a disk or ring of absorbent material, which is loaded with the test substance, is added to the sample. After a time of equilibration, the free concentration will have become proportionate to the concentration in the disk or ring, governed by the partition coefficient between water and the disk or ring material. If this partition coefficient has been predetermined, and if the amount of substance in the ring or disk by far exceeds the amount to be dissolved in the medium, then the free concentration in each sample can be easily calculated, and does not need to be measured. A more detailed description of the method can be found in Smith *et al.* (Smith *et al.*, 2010). This reference, as well as Smith *et al.* (Smith *et al.*, 2013), give examples of how passive dosing can be applied to *in vitro* tests.

2782
2783
2784

In vitro to *in vivo* extrapolation (IVIVE) refers to the qualitative or quantitative transposition of experimental results or observations made *in vitro* to predict phenomena *in vivo*, on full living organisms. When the response of the *in vitro* test is plotted against the free

2785 concentration (or the nominal concentration only in case it can be demonstrated/estimated
2786 this approximates the free concentration), toxicity parameters such as the EC₅₀ or a
2787 benchmark concentration (BMC) can be derived from the obtained curve. This *in vitro*
2788 toxicity parameter can be used as point of departure (PoD) for *in vitro* test circumstances and
2789 directly applicable to *in vivo* extrapolations (Leist, 2014, Blaauboer 2012). The
2790 corresponding *in vitro* concentrations can be converted into relevant plasma concentrations
2791 by taking the protein and lipid concentrations in plasma and cell culture medium into account
2792 (Bosgra and Westerhout, 2015; Zimmer *et al.*, 2014). In a final step, this concentration can be
2793 used as input for physiologically based pharmacokinetic (PBPK) models to estimate the dose
2794 that would result in the respective plasma concentration in man. This way an external
2795 benchmark dose (BMD) can be obtained. PBPK models describe the kinetic processes *in*
2796 *vivo*, relating external doses to tissue concentrations in time. For these models, some
2797 physical-chemical properties of the test substance need to be known, as well as some kinetic
2798 parameters such as the fraction absorbed, rate of metabolism, tissue partition coefficients,
2799 protein binding coefficients and urinary excretion rate. An example of how this can be
2800 performed is described in Louisse *et al.* (Louisse *et al.*, 2010). Good modelling practices for
2801 PBPK models have been described by Loizou *et al.* (Loizou *et al.*, 2008). The
2802 recommendations from a joint EPAA - EURL ECVAM on how PBTK modelling platforms
2803 and parameter estimation tools could enable animal-free risk assessment are reported in
2804 Bessems *et al.*, (Bessems *et al.*, 2014).

2805 **6.4 Reference and control items**

2806 The inclusion of relevant reference and control items, and the setting of acceptance criteria on
2807 the basis of historical data, is essential for regulatory applicability of *in vitro* methods. By
2808 including the correct reference and control items, the data set obtained from the *in vitro*
2809 method will demonstrate the correct functioning of the test system and the method used for
2810 analysis and therefore the validity of the executed experiments.

2811 Reference items can be one or more items where a specific readout and well-known response
2812 is expected (OECD, 2004a). The reference item(s) is used to provide a basis for comparison
2813 with the test item or to validate the response of the test system to the test item i.e., provide a
2814 known measurable or observable response. Notable requirements for the reference chemicals
2815 include having a well-defined chemical structure and purity and availability from commercial
2816 sources without prohibitive costs, hazards or disposal considerations. Reference items should
2817 be specific to the endpoint being measured. Documentation should be available to justify the
2818 selection of the reference items. While these reference chemicals are meant to represent the
2819 minimum to evaluate the performance of *in vitro* methods they can be substituted with other
2820 chemicals from the same class or substances possessing similar activity.

2821 Reference item(s) should be tested for batch-to-batch variability and be appropriately
2822 characterised (e.g., purity, stability) and identified (e.g., CAS number) (GLP Principles,
2823 No.1). Solubility, stability, and purity need to be established for each reference item used,
2824 and acceptance criteria based on historical data need to be developed. The continuous
2825 monitoring of the reference items, e.g. in the format of a control chart, is important to prove
2826 that the *in vitro* method continues to perform within the limits, and is consistent over time.

2827 Control items are used to control the proper performance of the test system (OECD, 2004a).
2828 Monitoring and recording performance against negative and positive control items may
2829 constitute sufficient proof for the responsiveness of a given test system. Non-response of the
2830 test system to the negative control and response to the positive control, within the acceptance
2831 criteria, show that the test system is “reactive” and behaves as expected. For positive,

2832 negative, and vehicle control items (in those cases where the negative control is different
2833 from the solvent control), it may or may not be necessary to determine concentration and
2834 homogeneity, since it may be sufficient to provide evidence for the correct, expected
2835 response of the test system to them. Such evaluation may consist of documented evidence
2836 that the response of the respective test systems to these positive, negative, and/or vehicle
2837 control items does not deviate from the historical control values recorded in the test facility,
2838 which should furthermore be comparable to published reference values. Guidance on how to
2839 compile and use historical data can be found in literature. Hayashi (Hayashi *et al.*, 2011)
2840 describes the compilation and use of historical data specifically for genotoxicity data, but this
2841 approach can also be applied in a broader context. A more general approach is described by
2842 Yoshimura (Yoshimura and Matsumoto, 1994).

2843 The positive control item is one that has a consistent and predictive effect in the *in vitro* test
2844 system. It is often used to assess aspects such as the characteristics of the test system and
2845 whether the *in vitro* methods gives reproducible results for this positive control over time in
2846 accordance with historical data. Moreover, a positive control-induced endpoint value that
2847 meets the predefined acceptance criteria, assures that when a test item has no effect on the *in*
2848 *vitro* method, the negative result is not caused by an error during the testing phase. Therefore,
2849 the results from the control test items are of utmost importance to show that a valid run has
2850 been performed when test data for each unknown test item are submitted to regulatory
2851 authorities.

2852 The response to reference items may be used to characterise the test system. Therefore,
2853 particular attention should be given to the definition of appropriate acceptance criteria for the
2854 response to each selected reference item. Acceptance criteria should be available when the
2855 method can be considered as fully developed and subsequently enter the validation process
2856 As such, acceptance criteria, should be available prior to the start of the *in vitro* study when
2857 used for regulatory purposes.

2858 Acceptance criteria for reference items, positive, negative and vehicle control items (e.g.,
2859 purity, known response, cytotoxicity) should be developed and detailed in the SOP.

2860 Records of receipt, storage, preparation and use should be available to allow for a full
2861 reconstruction of the history and use of each reference item.

2862 **6.5 Use of proficiency chemicals**

2863 While reference chemicals are used to validate the *in vitro* method, proficiency chemicals are
2864 used to test the suitability of a laboratory to carry out the *in vitro* method and produce
2865 acceptable results.

2866 Transfer of a new method to another laboratory is part of the development and validation
2867 process before a new *in-vitro* test can be accepted by regulators. The developer of the *in vitro*
2868 method should provide a list of proficiency chemicals. Proficiency chemicals should be fully
2869 characterised in the new *in vitro* method. By comparing the results of testing the selected
2870 proficiency chemicals by both laboratories, the inter-laboratory transfer of the new *in vitro*
2871 method can be evaluated

2872 Proficiency chemicals selected for validation of the assay should cover a range of responses
2873 within the *in vitro* method's dynamic range. In addition, a good proficiency chemical set
2874 should include representatives of several chemical classes within the applicability domain of
2875 the *in vitro* method [chemical classes for which the test is intended]. The number of
2876 proficiency chemicals should be chosen in such a way that a new laboratory can be confident
2877 that their results will be acceptable and robust. Since this greatly depends on the properties of

2878 the method, some methods may require 5 proficiency chemicals, while for others up to 20
2879 compounds should be tested. On average, a list of proficiency chemicals usually contains
2880 around 10 materials when the method is transferred to an OECD test guideline. In this way,
2881 new laboratories undertaking the *in vitro* method can demonstrate their proficiency.
2882

DRAFT

2883 **7 Standard Operating Procedures (SOPs)**

2884 According to the Principles of GLP, documented procedures which describe how to perform
2885 tests or activities normally not specified in detail in study plans or test guidelines are defined
2886 as standard operating procedures (SOPs). SOPs are a set of written documents describing
2887 routine or repetitive activities (e.g. *in vitro* methods and complementary procedures) that
2888 facilitate consistency in the quality and integrity of a product or end-result. SOPs are required
2889 by GLP. Outside the GLP environment, terms such as testing methods, instructions,
2890 worksheets, and laboratory operating procedures are often used. In principle SOPs are unique
2891 to an organisation, as they are part of an integrated quality assurance process of a facility.

2892 The aim of a SOP is to ensure that procedures are carried out in a reproducible way by
2893 qualified personnel. Therefore SOPs need to be written with sufficiently detailed but clear
2894 work instructions and be kept simple, to minimise the risk for misinterpretation.

2895 An *in vitro* method will be supported and documented by a number of different SOPs, forms,
2896 templates and worksheets. Besides the description of the main test procedure, also SOPs for
2897 supporting procedures (e.g. the handling of cell cultures, waste handling, cleaning
2898 procedures, operating and calibration instructions for the equipment, record keeping,
2899 reporting, archival, quality assurance procedures, etc.) need to be available and referred to.
2900 To avoid lengthy documents, the instructions are preferably divided into a series of short
2901 SOPs. SOPs should be readily available to personnel in each working area.

2902 **7.1 *In vitro* method SOP(s) development**

2903 The development of an *in vitro* method for regulatory testing purposes by using complex
2904 animal or human cells and tissues as test systems is a difficult and time-consuming task. In
2905 the initial stages of the development, the procedure will undergo many changes and each step
2906 needs to be described in laboratory records, which will crystallise into a laboratory procedure
2907 or a SOP along the test development process. During this period, also the historical data set
2908 of the model compounds will be collected. This dataset will be needed to define the critical
2909 and relevant end-parameters, the control and reference items and associated acceptance
2910 criteria.

2911 Once the method is sufficiently developed and all parameters are defined, the *in vitro* method
2912 SOP(s) is/are sufficiently "routine" and standardised to be used for an in house validation
2913 process during which the *in vitro* method is checked for repeatability (accuracy & precision),
2914 selectivity, sensitivity, and stability assessment over time. Likewise, its robustness is assessed
2915 (i.e. the influence of critical (external) parameters on the outcome parameters), as it is
2916 important to secure the test performances in different laboratory environments, albeit within
2917 defined boundaries.

2918

2919 **Figure 6: Evolution of a Standard Operating Procedure**

SOP Evolution
No routine => No SOP => no reproducibility
Method becomes routine => Draft SOP <ul style="list-style-type: none"> - Describe the <i>in vitro</i> method procedure - Historical data of reference items are generated in a controlled way
SOP Version 01 <ul style="list-style-type: none"> - Advanced <i>in vitro</i> method description - Acceptance criteria for valid/invalid experiments - Lists of needed equipment, reagents, consumables and reference items - Calculation of results
SOP Version xx <ul style="list-style-type: none"> - Further optimised procedure - Acceptance criteria for valid/invalid experiments - Calculation of results for test items - Data recording Forms, Data Calculation forms - SOP is robust

2920

2921 During the in house validation process, weaknesses can come to light that demand for
 2922 adaptation (optimisation) and which might also trigger the re-initiation of a new validation
 2923 cycle. Optimisation of the SOP should be performed by following a formal procedure. It is
 2924 critical that any parameter(s) to be changed should be chosen prior to the optimisation
 2925 process, including the steps to be performed: it is recommended to define these steps in
 2926 amendment in the validation study plan, so that every passage is recorded in a proper way.
 2927 Also, the historical data should be annotated to allow tracking, comparison and measurement
 2928 of the acquired optimisation.

2929 Upon a satisfactory completion of the validation process, the method development can be
 2930 finalised and the final set of SOPs associated with the *in vitro* method will be available.

2931 In addition, the test developer should be aware that if the established *in vitro* method makes
 2932 use of complex instrumentation and software, including data analysis and computer models
 2933 and if developed in-house (e.g. excel data analysis templates), this software will also require
 2934 documentation and full validation with SOPs for correct use, prior to transferring the *in vitro*
 2935 method to the validation bodies.

2936 Once an *in vitro* method has been validated and published, e.g. in the format of an OECD
 2937 Test Guideline (TG), the end-users will, from the published method, need to develop their
 2938 own set of SOPs which are applicable and integrated into their organisation to assure the
 2939 correct execution of the *in vitro* method within their facility environment.

2940 **7.2 Preparing Standard Operating Procedures**

2941 As indicated above, the guidance how to correctly perform the *in vitro* method is given in a
 2942 set of SOPs, covering how to execute the *in vitro* test but also SOPs referring to general
 2943 supporting procedures e.g. test system handling, solubility assessment, cytotoxicity

2944 measurement, equipment maintenance, calibration and cleaning; handling of test and
2945 reference items; record keeping, reporting, storage, and retrieval, etc. The reason for not
2946 having all these steps and processes described in one single SOP, but a set of
2947 methods/procedures referring to other specific SOPs is to have available and easy-to-handle
2948 documents, to be consulted by the personnel involved. However, to avoid deviating
2949 procedures over time, it is advisable to include cross-references between these SOPs rather
2950 than duplicating information as it might otherwise not be updated in all documents where it
2951 appears.

2952 SOPs should be written in the active voice and concisely explain a procedure in a step-by-
2953 step, easy-to-read format. The information presented should be unambiguous and not overly
2954 complicated. The document should not be wordy, redundant, or overly lengthy but simple
2955 and short. The inclusion of a flow chart to illustrate the process can help to make it clear.

2956 A SOP is best written by the people actually performing the work on a daily basis. The
2957 finalised SOP needs to be reviewed and approved by laboratory management. SOPs are not
2958 static documents and need to be systematically reviewed on a periodic basis and adapted
2959 whenever something changes (products/ equipment / facility...), in order to be possible to
2960 fully adhere to. As soon as a new version is approved, all concerned personnel need to be
2961 informed before the new version comes into force. From that date onwards the obsolete
2962 version needs to be made unavailable and adequately archived. To allow and control this, all
2963 SOPs need to have a unique identifier (Title / version number / approval date). It is also
2964 recommended to keep the revision history of the document in the document. SOPs should be
2965 formally authorised by test facility management.

2966 SOPs may be written according to a pre-set template and may contain (adapted from “EPA
2967 (USA) “Guidance for Preparing Standard Operating Procedures (SOPs) – EPA QA/G-6
2968 /2007”): (EPA, 2007) ;

- 2969 • Title page or page header (including: a unique name which describes the essence of the
2970 SOP; the name of the author or authors and the signatures of reviewers and approvers, an
2971 identification code -if applicable-, the revision number, the date of approval and effective
2972 date, and number of total pages in case it is not included in the subsequent pages).
- 2973 • Change log (list summarising the changes from each previous edition of the SOP).
- 2974 • Table of contents in particular necessary for longer SOPs. In some GLP test facilities short
2975 SOPs (e.g. 5 pages) do not always have a table of contents.
- 2976 • Purpose/objective – (identifying the intended use of the procedure).
- 2977 • Field of Application/Scope (identifying when the procedure is to be followed, and
2978 explaining who the intended user of the document is and what the method does and does
2979 not cover).
- 2980 • Summary of main steps (can be illustrated by a flow chart, if desired).
- 2981 • Definitions (defining any words, phrases, or acronyms used in the SOP).
- 2982 • Personnel qualifications/responsibilities (identifying any special qualifications users
2983 should have such as certification or training experience and/or any individual or positions
2984 having responsibility for the activity being described)..
- 2985 • Reference to all related SOPs and other documents (forms, templates, worksheets etc.) to
2986 be used together with the SOP (not necessarily as a separate section)
- 2987 • Health & safety warnings (indicating the possible risks to address and the personal
2988 protection equipment and containment equipment to use).
- 2989 • Procedure (including: a list of the equipment, reagents, consumables and reference
2990 materials, with their quality requirements; pre-analytical preparation: sample collection,

2991 handling and preservation; data processing; evaluation of results against acceptance
2992 criteria; post-analytical activities: reporting of results; sample and chemicals discarding).
2993 • Criteria, checklists, or other standards that are to be applied during the procedure.
2994 • Records management (specifically, e.g., as forms to be used and locations of files).
2995 • *In vitro* method acceptance criteria section - describe any control steps and provisions for
2996 review or oversight prior to acceptance of the results.
2997 • Reference Section - cite all references that have been consulted during the authoring of the
2998 SOP.
2999 • In addition to the first page, all the subsequent pages of a SOP should include the title, the
3000 identification code (if applicable), the revision number, the page number (and the total
3001 number of pages).

3002

DRAFT

3003 **8 Performance of the method**

3004 To assess *in vitro* method performance one needs to define the selection of reference and
3005 control items as these are used to check the performance of the method. Of critical
3006 importance is also how acceptance criteria (often based on historical data) are developed and
3007 defined for the method regarding its performance.

3008 The number of replicates for each condition (e.g. chemical concentration) should be
3009 specified. During method development the number of replicates must be chosen using
3010 appropriate statistical methods. For example, a statistical power analysis (Crawley, 2015) can
3011 be used to calculate the desirable number of replicates to detect a defined difference between
3012 treatments with pre-set levels of confidence (Krzywinski and Altman, 2013). However, one
3013 should be aware that this number can be too high to be useful in practice. Alternatively the
3014 statistical power is provided for the chosen number of replicates.

3015 Additionally, when multiple concentrations of a test item are tested, the mathematical curve
3016 fitting model (e.g. dose-response curve) can be computed with increasing number of
3017 replicates. The lowest number of replicates that gives satisfactory variability of the parameter
3018 of interest (e.g. IC₅₀ within acceptable limits) can be used in future studies (Assay Guidance
3019 Manual, HTS Assay Validation, Section 5.2 assay guidance (Iversen *et al.*, 2004)). Apart
3020 from these statistical considerations, sometimes practicalities such as cost and availability of
3021 replicates may also play a role in the selection process. However the impact of reducing
3022 replicates should always be subjected to careful analysis and corresponding power should be
3023 given.

3024 Similarly, the number of independent experiments needs to be evaluated. For instance, *in*
3025 *vitro* methods with a high degree of inter-experimental variability, such as those using
3026 primary tissues, may need a higher number of independent experiments compared to *in vitro*
3027 methods employing continuous cell lines.

3028 Statistical methods can be very useful in the process of optimising new *in vitro* methods.
3029 Newly developed *in vitro* methods can include steps that may be interpreted differently and
3030 result in significantly different readouts for the different interpretations. To obtain an *in vitro*
3031 method that leads to accurate, reliable and robust readouts, the results of several
3032 combinations of any changes in the *in vitro* method would have to be assessed. Given the big
3033 number of variables that may or may not be changed, it is nearly impossible to approach all
3034 these combinations experimentally. In such cases, factorial design of experiments may be of
3035 added value. These are efficient at evaluating the effects and possible interactions of several
3036 factors (independent variables). A statistical approach predicting the effect of changes in
3037 testing method steps on the observed readout (known also as method robustness assessment)
3038 would allow for the development of an efficient *in vitro* method design, since the
3039 experimental robustness check can be based on a much smaller subset of combinations (Box
3040 *et al.*, 2005; Groten *et al.*, 1997).


3041

3042 **8.1 Plate layout**

3043 The plate layout should be such that cross-contamination (e.g. between test items) can be
3044 controlled for (replicates). It should also allow for cross-plate comparison by using
3045 appropriate reference and control items. An example of an experimental 96-well plate layout
3046 using reference and control items is shown in [Figure 1](#) (Coecke *et al.*, 2014).
3047

3048
3049

Figure 7: Example of plate layout. PC positive control, RI reference item, NC negative control, UC untreated control, VC vehicle control, TI test item



3050
3051

3052 The example plate layout (Figure 6) minimises potential edge effects (difference between
3053 outer and inner wells due to evaporation). For strategies to assess edge effects and drift see
3054 Section [Error! Reference source not found.8.7](#). Another way to assess plate drift is to
3055 include vehicle controls (VC) on both the left and right side of the plate. Left and right VCs
3056 should not differ more than a certain percentage for the plate to be accepted. For example, in
3057 the test in appendices, the following criteria is used: a test meets acceptance criteria if the left
3058 and the right mean of the VCs do not differ by more than 15% from the mean of all VCs
3059 (National Institutes of Health, 2001).

3060 Comparison across plates, evaluation that different plates provide comparable data, should
3061 also be taken into account for the plate design.

3062 The inclusion of relevant reference and control items, and setting of acceptance criteria on the
3063 basis of historical data, is essential for regulatory applicability of *in vitro* methods and should
3064 be considered when developers decide on their plate lay-out. By including the correct
3065 reference and control items, the data set obtained from the *in vitro* method will demonstrate
3066 the correct functioning of the test system and the method used for analysis and therefore the
3067 validity of the experiments executed.

3068 In addition, certain reference chemicals may be volatile (e.g. solvents) or may contaminate
3069 neighbouring wells by capillary action, the wicking effect (Sullivan, 2001) and this may need
3070 to be taken into account in designing plate layouts. For instance, the commonly used cell lysis
3071 surfactant Triton X can affect cell viability in neighbouring wells and should be used at low
3072 concentrations or separated from cell-containing wells by placing wells containing media or
3073 buffer in-between.

3074 8.2 Data analysis

3075 When data need to be transformed by formulas for normalisation, computer scripts and/or
3076 any model equations fit to the data, this should be documented in the SOP (OECD, 2014).
3077 Formulas for normalisation (checked for accuracy) should be documented, validated (when
3078 implemented in electronic format) and disclosed along with a description and justification of
3079 the controls used in the calculation. Computer scripts used to process raw data (e.g. Excel
3080 spreadsheets, scripts, macros etc.) should be validated and fully documented.

3081 When a model equation is fitted to the data (as in the case of dose-response curves, standard
3082 curves or other models), the equations and reasoning behind their choice need to be stated.
3083 For example, when fitting a dose-response curve, the type of the equation fitted to the data
3084 should be stated (e.g., four parameter logistic curve), along with any constraints (e.g. top
3085 constrained to 100% in normalised data) and weightings (e.g. by inverse data uncertainty)
3086 applied (Motulsky and Christopoulos, 2004). Furthermore, the software name and version
3087 used to fit the equations should be listed/stated, as well as the confidence interval of the
3088 measurements of interest (e.g. IC50 values) and the relevant goodness of fit parameters (R-
3089 square, sum of squares etc.) stated. Also see section [Error! Reference source not found.8.8](#)
3090 of the current document and section 4.3 of OECD 211 (OECD, 2014) (OECD, 2014) as the
3091 same data will need to be reported in the late stages of method development during the
3092 assessment of method performance.

3093 8.3 Data-intensive *in vitro* methods

3094 The last decennia brought a paradigm shift in toxicity testing of chemical compounds, relying
3095 more on less expensive and higher throughput high-content screening *in vitro* methods. They
3096 allow the processing of hundreds or thousands of compounds simultaneously enabling the
3097 identification of mechanisms of action, and ultimately facilitating the development of
3098 predictive models for adverse health effects in humans. Furthermore, image analysis and
3099 genomics-based *in vitro* method read-outs are getting more popular for *in vitro* method
3100 developers due to the data rich information obtained with such methods.

3101 The utility of "big data" for regulatory safety assessment has been discussed recently, for
3102 example omics (ECETOC 2013) or High Throughput screening (Judson 2013). These data
3103 could be used in various contexts, such as supporting evidence for read-across, defining
3104 categories or to allow the design of Integrated Testing Strategies (ITS). Still, most
3105 applications have focused on screening and prioritisation as in the EPA ToxCast program
3106 (Judson 2010).

3107 Although some technologies have been extensively used for more than a decade, e.g.
3108 microarrays, debate is still ongoing about the reproducibility of experiments and the
3109 comparability of results at different sites and on different platforms. Consensus is still to be
3110 achieved concerning best practices in many critical aspects such as the experimental design
3111 and protocols for sample preparation and handling, data processing, statistical analysis and
3112 interpretation. One of the important challenges to be addressed for regulatory acceptance is
3113 related to quality assurance (Bouhifd *et al.*, 2015). The maintenance of high standards is
3114 essential for ensuring the reproducibility, reliability, acceptance, and proper application of the
3115 results generated. A certain level of standardisation is also needed since "big data" are
3116 generated using diverse technological platforms and various biochemical, analytical and
3117 computational methods, producing different data types and formats.

3118 8.4 Acceptance criteria

3119 Because the *in vitro* tests for the prediction of human health effects are intended to predict a
3120 quantitative endpoint response predictive of the degree of human hazard, it is essential that
3121 the test performs consistently over time and across laboratories. The positive control response
3122 plays a major role in determining that consistency ([Error! Reference source not
3123 found.Figure 40](#)). To do that, the positive control needs to induce a known change in the
3124 endpoint measure(s) that falls at a point within the dynamic range of the test. That way,
3125 increased and decreased magnitudes of response can be measured.

3126 In many cases, the acceptable range for the positive control response is based on a statistical
3127 analysis of the historical range for the assay (firstly produced at the developer's laboratory,
3128 then supplemented by data from validation studies, and finally produced in the laboratory
3129 performing the *in vitro* method). For data, which follows an approximately normal
3130 distribution, the mean and standard deviation are calculated and the acceptance range is set at
3131 for instance plus or minus two standard deviations from the mean response. For example, the
3132 Bovine Corneal Opacity and Permeability (BCOP) *in vitro* method using 100% ethanol as the
3133 positive control (selected because it induces both increased opacity and fluorescein passage),
3134 the mean published *in vitro* score (opacity + 15×permeability) was 51.6 ± 6.2 (mean \pm
3135 standard deviation SD), which would set the acceptable range ($\pm 2SD$) to be 39.2 to 64.0
3136 (n=1171 trials) (Harbell *et al.*, 2014).


3137 Development of the positive control should begin as early in the *in vitro* method development
3138 process as is practical (Hartung *et al.*, 2004). Positive control performance can help identify
3139 dependent variables that impact assay consistency. It is essential that the positive control
3140 should be incorporated into the assay every time it is performed. The positive control is
3141 particularly important during the *in vitro* method development stage where data are being
3142 generated to define the prediction model for translating the *in vitro* endpoint data into a
3143 prediction of human (animal) health effects. A robust positive control historical record will
3144 facilitate technology transfer to another test facility as part of the *in vitro* method transfer
3145 process.

3146 The challenge of evaluating *in vitro* method performance over time and across laboratories is
3147 not unique but is often not addressed with *in vivo* test system. *In vitro* test systems have a
3148 large number of dependent variables. The operator is wholly responsible for maintaining test
3149 system homeostasis and consistency.

3150 The importance of positive controls is exemplified by a Draize eye irritation study "Study of
3151 intra- and inter-laboratory variability in the results of rabbit eye and skin irritation tests" in
3152 (Weil and Scala, 1971). [Figure 8](#) shows the 24-hour Draize Maximum Average
3153 Scores (MAS) for two of the test materials: Compound F (46% aqueous triethanolamine
3154 lauryl sulfate) and Compound A (95% ethanol). The data are arrayed according to the MAS
3155 reported for Compound F and A. Without a positive control, there is no means to compare
3156 data from one laboratory to another or one laboratory over time.

Formatted: Font: Times New Roman,
12 pt, Font color: Text 1

3157
 3158 **Figure 8: Sample data from the Draize Eye Irritation test on two chemicals at 24 hours**
 3159 **after instillation**




3160
 3161 The data are arrayed in order of increasing maximum average score (MAS) for compound F.
 3162 The corresponding MAS for compound A is paired with the MAS for compound F from that laboratory.

3163
 3164 [Figure 9](#)[Figure 10](#) shows three concentration response curves from a keratinocyte-based
 3165 cytotoxicity assay (neutral red uptake endpoint) treated with sodium lauryl sulfate. The
 3166 difference between these curves shows the information that can be gained from a concurrent
 3167 positive control. They also illustrate the importance of having multiple concentration points
 3168 in the active range so as to narrowly define the 50% relative survival point. The 50% survival
 3169 point (concentration) for the positive control compound is often the basis for establishing the
 3170 acceptance criteria for a dilution-based cytotoxicity assay. This approach allows increased
 3171 and decreased sensitivity to be readily identified.


Formatted: Font: Times New Roman, 12 pt, Font color: Text 1

3172 **Figure 9: A 96-well based cytotoxicity assay using human keratinocytes**


A. An example of normal positive control concentration response curve for this assay.



B. This is an example of poor cell seeding leading to inconsistent neutral red uptake in replicate wells.



C. An increase in cell number per well has induced a right shift of the SDS cytotoxicity curve.



3173


3174 In selecting a positive control, it is desirable but not always possible to select a single positive
 3175 control to address all endpoints or exposure conditions. In the genetic toxicity *in vitro*
 3176 methods, such as the bacterial reverse mutation (Ames) assay, two positive controls are used
 3177 for each bacterial strain to address direct mutagenic activity and metabolic activation of a
 3178 promutagen with rat liver S9 (Zeiger *et al.*, 1988). In the BCOP, one positive control is used
 3179 for the liquids exposure testing method and another for the solids exposure testing method
 3180 (OECD, 2009).

3181 Finally, the importance of testing the positive control concurrently with the unknown test
 3182 materials each time the *in vitro* method is run is illustrated in [Figure 10](#) (Ulrey et
 3183 al., 2005). This Figure shows the quality control chart for the BCOP using the ethanol control
 3184 for each test performed over a period of 2.5 years. The acceptable upper limit is shown with a
 3185 light line between 60 and 70. There is a cluster of points that extend above the upper limit (in
 3186 mid-2002). These were failed assays which were repeated. The basis of the failures was not
 3187 immediately clear since the eyes looked quite normal upon arrival in the laboratory.
 3188 However, the patterns became clear and the cause was traced to improper handling of the
 3189 eyes in the abattoir. Without the concurrent positive control data, it would not have been
 3190 possible to identify the problem and prevent inappropriate data from being reported. Isolated
 3191 tissues or tissue constructs as test systems can be difficult to properly evaluate visually and so

Formatted: Font: Times New Roman, 12 pt, Font color: Text 1, Do not check spelling or grammar

3192 the functional test provided by the concurrent controls is often the only way to measure their
3193 integrity.

3194 **Figure 10: A quality control chart for the ethanol positive control for the BCOP showing**
3195 **the acceptable range and deviations above and below (March 2002 to November 2005-**
3196 **courtesy of IIVS)**



3197
3198
3199 In the same vein, establishing an expected range for the negative control is important to
3200 assure that the test system performs normally and to determine the detection limit of the
3201 assay.
3202 Finally, it is also important to establish the cut-off value of the acceptance criteria for
3203 quantitative endpoints, i.e. whether it will be accepted not less to or higher than a specified
3204 value, including the number of significant figures. Having a result of 50.4, and an acceptance
3205 criteria of >50 is not acceptable as the number of significant figures are not comparable
3206 (using an acceptance criteria of ≥ 50 would be acceptable), however the preferred approach is
3207 to specify the same number of significant figures both for the acceptance criteria and the
3208 measured result e.g. acceptance criteria of > 50.0 (one significant figure) and a measured
3209 result of 50.4 (would be accepted).

3210 **8.5 Dynamic range/range of application**

3211 The response of the instrument and the *in vitro* method with regard to the readouts of interest
3212 should be known, and should be evaluated over a specified concentration range. An important
3213 characteristic of the method performance, for both quantitative and qualitative *in vitro*
3214 methods is often referred to as Limit Of Detection (LOD). It determines the lowest actual
3215 concentration or signal that can be consistently detected with acceptable precision, but not
3216 necessarily quantified. In case of normally distributed data, the LOD is often determined as

3217 the concentrations at the average response + 3* SD of the negative control range, as this gives
3218 only 1% chance of a false positive.

3219 For quantitative measurements, the boundaries of the dynamic range are determined by the
3220 lowest and highest analyte concentrations that generate results that are reliably produced by
3221 an *in vitro* method without dilution of the sample. To ensure the ability of the method to yield
3222 test results that are directly proportional to the analyte concentration (cell number, enzyme
3223 activity) the linearity of the method should be ensured within the given range. Subsequently,
3224 to facilitate efficient method transfer, the linear regression parameters (correlation
3225 coefficient, slope, y-intercept as well as residual sum of squares) should be submitted along
3226 with a plot of the data. When the upper limit is exceeded, samples fall outside of the linear
3227 range, they may need to be diluted (when the upper limit is exceeded), if possible. In case the
3228 samples give a result below the lower limit of the linear range, it may be useful to change to a
3229 different apparatus with a higher sensitivity or adapt the sample preparation to obtain higher
3230 concentrations, if possible.

3231 The lower limit of linearity is frequently referred to as the lower limit of quantification
3232 (LLOQ) and the upper limit of linearity as the upper limit of quantification (ULOQ). The
3233 upper limit of linearity may be restricted by the highest available concentration in a sample or
3234 by the saturation of the signal generated by the instrument. The lower limit of quantification
3235 may also be based on the negative control range.

3236 For certain assays/methodologies, equations other than the linear can be fit as a standard
3237 curve, provided that the user is operating within the range of the assay/equipment (see 4.1
3238 and above). However, it is recommended that the simplest model that adequately describes
3239 the concentration-response relationship is used. Selection of weighting and use of a complex
3240 regression equation should be justified. (Burd, 2010; EMEA, 2011; FDA, 2001; Viswanathan
3241 *et al.*, 2007).

3242 **8.6 Signal intensity**

3243 The signal to noise (S/N) ratio is frequently applied for methods, which exhibit background
3244 noise (observed as the variation of the blanks) as baseline. It is calculated by comparing
3245 measured signals from samples with positive control item with those of blank samples. Note:
3246 concentrations of positive controls should be chosen in a low to medium range (giving a clear
3247 signal, but rather in the lower range of the linear range). While the S/N ratio is useful for
3248 tests, in which the blank (negative control) and the positive control exhibit comparable levels
3249 of variation, measures that take into account the levels of variability of both the blank and the
3250 sample may be more appropriate for biochemical and cell biology research.

3251 They include:

- 3252 • Signal windows (SW) (Sittampalam, 1997): these allow for the calculation of the
3253 separation between the sample (positive control) and the blank (negative
3254 control/solvent) including the variability of both measurements:

$$3255 \quad SW = \frac{Mean_{sample} - Mean_{control} - 3*(SD_{sample} + SD_{control})}{SD_{sample}}$$

3256 A common SW acceptance criterion is $SW \geq 2$ for biological assays.

- 3257 • Z-factors (Zhang *et al.*, 1999): they represent both the assay signal dynamic range and
3258 the data variation associated with the signal measurements, and therefore are suitable

3261 for assay quality assessment. Z-factors ≥ 0.4 indicate acceptable separation between
3262 the signal and the blank.
3263

$$Z = 1 - \frac{3 * (SD_{sample} + SD_{control})}{Mean_{sample} - Mean_{control}}$$

3264

3265 **8.7 Signal variability and plate uniformity assessment**

3266 The Assay Guidance Manual (HTS Assay Validation, Section 3 (Iversen *et al.*, 2004))
3267 proposes a signal variability and plate uniformity assessment when developing new *in vitro*
3268 methods. This assessment serves to quantify the variability of assay output for the minimum
3269 (base, background), middle (close to EC₅₀) and maximum (highest/maximal) assay readings.
3270 By employing an interleaved format one can calculate the coefficient of variation of the min,
3271 mid and max signals, Z-factors and signal windows. Exemplary acceptance criteria have been
3272 proposed as coefficient variation CV $\leq 20\%$, Z-factor ≥ 0.4 and SW ≥ 2 . In addition, the data
3273 from this assessment can be used to derive the number of required replicates, using power
3274 analysis (Crawley, 2015; Krzywinski and Altman, 2013).

3275 Furthermore, seeding density variation, plate edge effects and drift can also be examined in
3276 this experimental setup using heat maps (e.g. created with Excel) or graphical plots.

3277 Edge effects are differences in growth in outer wells compared to inner wells. They are often
3278 due to uneven evaporation rates or plate stacking and can be a source of variation, as outer
3279 wells can often present as outliers compared to inner wells. Edge effects can be detected by
3280 using a heat map of the readouts from the plate uniformity assessment, or by plotting the
3281 signal readouts as a function of well position.

3282 The plate uniformity assessment can also detect drift effects, i.e. trends in signal from left-
3283 right or top-bottom. Drifts can be due to seeding density variation during the process of initial
3284 cell seeding in plates. For example, cells may be settling down in the master vessel which is
3285 used to store a cell suspension used to seed a particular plate. Additionally, using the same set
3286 of tips on a multichannel pipette while pipetting cells in media compositions prone to
3287 foaming, may compromise the accuracy of the seeding. Higher variability, which cannot be
3288 resolved via technique optimisation may require increased numbers of replicates, more drug
3289 concentrations used to calculate dose-response or a higher numbers of independent
3290 experiments (Iversen *et al.*, 2004).

3291 **8.8 Reliability of endpoint calculations**

3292 The coefficient of variation, signal windows and Z-factors are poorly suited to determine the
3293 variation and reproducibility for potency estimates in dose-response experiments. Comparing
3294 inflection point parameters such as IC₅₀ or EC₅₀ is further complicated by the fact that they
3295 are not normally distributed and require logarithmic transformation before averaging or
3296 statistical comparison (Christopoulos, 1998; Motulsky and Christopoulos, 2004).

3297 **8.8.1 Confidence intervals for EC₅₀s**

3298 While the parameters derived from non-linear regression analysis, like IC₅₀ and EC₅₀ values,
3299 are not normally distributed, their logarithms are (the potency estimates are log-normally
3300 distributed). Therefore, when multiple independent experiments are run, EC₅₀ values from
3301 multiple runs can be combined by using their geometric mean, calculated by averaging the

3302 logEC₅₀ values. The standard error of the mean of the logEC₅₀ values can be further used to
3303 calculate 95%-confidence intervals for the EC₅₀ of a number of independent experiments.
3304 Another way of quantifying dose-response curve reproducibility is by using the minimum
3305 significant ratio (MSR) (Eastwood *et al.*, 2006).

3306 **8.8.2 Minimum significant ratio**

3307 The MSR is a parameter that can be used to quantify assay reproducibility and resolution (the
3308 smallest ratio between compound potency which can be detected in the *in vitro* method). The
3309 so-called 'replicate experiment study' (Assay Guidance Manual, HTS Assay Validation
3310 Section 4 (Iversen *et al.*, 2004)) uses the differences between the potency estimates for 20-30
3311 compounds in two independent assay runs to estimate the within-run assay variability over a
3312 range of potencies. The SD of the difference between logIC₅₀s in run 1 and run 2 is
3313 computed. The $MSR=10^{2 \cdot SD}$ is computed and example acceptance criteria have been given as
3314 $MSR < 3$ (meaning up to three-times potency difference between runs is judged acceptable).

3315 The control compound MSR can be calculated for controls where data for at least six runs is
3316 available and gives an estimate of between-run variability. The control compound MSR is
3317 calculated as $MSR=10^{2 \cdot s}$, where s is the standard deviation of the log₁₀EC₅₀ values across
3318 runs, assuming one EC₅₀ result per run (Haas *et al.*, 2004).

3319 **8.8.3 Variability and outliers**

3320 The degree of variability judged acceptable should be given for the critical parameters of the
3321 assay- CV, Z-factors, SW, MSR or 95%-CI for IC₅₀ measurements.

3322 Criteria to detect/remove outliers should be stated and the reasoning behind should be given
3323 (Motulsky and Brown, 2006; Pincus, 1995).

3324 Variability should be determined for within runs (= repeatability), for within the laboratory (=
3325 intra-laboratory reproducibility), and for between laboratories (inter-laboratory
3326 reproducibility) during validation (Burd, 2010; EMEA, 2011; FDA, 2001; Viswanathan *et al.*,
3327 2007).

3328 **8.8.4 Transferability: in between lab transfer**

3329 The transferability of an *in vitro* method from the developer laboratory to a second laboratory
3330 is a crucial step for demonstrating the robustness of the established procedures and/or SOPs.
3331 This step is necessary to evaluate the practicability of the test and to identify possible sources
3332 of within- and between-laboratory variability. Moreover it provides also an estimation of the
3333 amount of training that will be necessary to successfully transfer the test to an inexperienced
3334 laboratory. A good inter-laboratory reproducibility proves transferability.

3335 **8.9 Accuracy, reliability and uncertainty**

3336 The usefulness of an *in vitro* method depends on its accuracy and reliability to correctly
3337 classify chemicals according to its stated purpose (e.g., sensitivity, specificity, positive and
3338 negative predictivity, false positive and false negative rates). These values often are obtained
3339 from continuous data and categorised accordingly (e.g. as strong, weak). In such case cut off
3340 values are used and their impact on the accuracy and reliability should be taken into account.
3341 The use of confidence bounds taking into account the distance to these cut off values should
3342 be considered.

3343 There are a number of measures that can be used to assess method performance and the
 3344 confidence of the test predictions (Gaddis and Gaddis, 1990), based on Table 6. Below is an
 3345 example that has been used in the validation process of different *in vitro* methods that are
 3346 now OECD test guideline methods.

3347
 3348 **Table 62: Possible outcomes of a test result of a chemical in a validation**

		Gold Standard Test	
	Total number	Condition positive	Condition negative
Test outcome	Test outcome positive	True positive (TP)	False positive (FP)
	Test outcome negative	False negative (FN)	True Negative (TN)

3349 The sensitivity is the ability of a test to reliably classify positive substances:
 3350
 3351

$$\text{Sensitivity (\%)} = 100 \cdot \frac{TP}{(TP + FN)}$$

3352

$$\text{False negative rate } (\beta) = 1 - \text{sensitivity} = \frac{FN}{(TP + FN)}$$

3353 The specificity is the ability to reliably classify negative substances:
 3354
 3355

$$\text{Specificity (\%)} = 100 \cdot \frac{TN}{(TN + FP)}$$

$$\text{False positive rate } (\alpha) = 1 - \text{specificity} = \frac{FP}{(TN + FP)}$$

3356 The positive predictive value (PPV) is the proportion of correct positive responses among
 3357 materials testing positive by an *in vitro* method:
 3358
 3359

$$\text{PPV (\%)} = 100 \cdot \frac{TP}{(TP + FP)}$$

3360 The negative predictive value (NPV) is the proportion of correct negative responses among
 3361 materials testing negative by an *in vitro* method:
 3362
 3363

$$\text{NPV (\%)} = 100 \cdot \frac{TN}{(TN + FN)}$$

3364
 3365 What level of sensitivity, specificity, etc. is acceptable is not standardised. Such levels are
 3366 also very dependent on the list of chemicals with which they are determined, therefore strict
 3367 boundaries in acceptable levels for these accuracy parameters are not realistic. Generally
 3368 though, sensitivities below 75% should not be accepted.

3369 A discussion on what is a gold standard for an *in vitro* test, is very important here, as this is
 3370 where many validations of *in vitro* tests fail currently. Comparison to an *in vivo* test is

3371 problematic, because of several reasons such as: (i) species differences e.g. in case of rat *in*
3372 *vivo* and human cells *in vitro*; (ii) *in vivo* tests may include more pathways to a certain effect
3373 (endpoint) than the *in vitro* test. In this case the *in vitro* test would detect fewer *in vivo*
3374 (“true”) positives. The situation may also be reverse, e.g. due to compensations *in vivo* not
3375 existing *in vitro*; (iii) the test endpoints *in vivo* and *in vitro* may not be 100% identical; (iv)
3376 the *in vitro* test does not reflect *in vivo* toxicokinetics and metabolism, which can lead to
3377 fewer (e.g. in case of toxification by metabolic enzymes) or increased positives (e.g. in case
3378 of detoxification or an absorption barrier) being found by an *in vitro* test without any
3379 metabolic competence, compared to an *in vivo* test; (iv) lack of mechanistic understanding of
3380 the underlying biological processes resulting in a positive or negative effect *in vivo*. Hereby it
3381 should be mentioned that many routinely used animal-based methods correlate poorly to the
3382 effects observed in humans. It might happen, therefore, that an *in vitro* test with human
3383 tissues shows very poor sensitivity and specificity when compared to the animal test and
3384 would therefore seem to be a poor test. Comparison to human data as the golden standard
3385 could in that case show whether the *in vitro* test is really poor, or actually a better test than
3386 the animal test.

3387 The chemicals used for determining accuracy and reliability should therefore be selected with
3388 care (only the ones working through the pathways included in the *in vitro* method) and this
3389 assessment should include toxicokinetic information on these chemicals.

3390 EFSA’s Scientific Committee defines uncertainty as referring to “all types of limitations in
3391 the knowledge available to assessors at the time an assessment is conducted and within the
3392 time and resources available for the assessment”. Examples include:

- 3393 • Possible limitations in the quality and representativeness of data.
- 3394 • Comparing non-standardised data across countries or categories.
- 3395 • Choosing one predictive modelling technique over another.
- 3396 • Using default factors (such as the weight of an average adult).

3397 Since 2013, EFSA’s Scientific Committee has been developing guidance²⁰ on how to assess
3398 uncertainty in a structured and systematic way. The aim is to offer a tool-box of
3399 methodologies – both quantitative and qualitative – and related training for EFSA’s Scientific
3400 Panels and staff as well as for other organisations (e.g. researchers, national authorities) that
3401 carry out scientific work on EFSA’s behalf.

3402 In addition, IPCS has published a guidance document on "Evaluating and expressing
3403 uncertainty in hazard characterization"(IPCS, 2014), which is focussed on animal test results.
3404 However, described concepts such as the benchmark dose (BMD), including its uncertainty
3405 characterisation, can also be used for *in vitro* assays.
3406

3407

²⁰ <https://www.efsa.europa.eu/sites/default/files/160321DraftGDUncertaintyInScientificAssessment.pdf>

Field Code Changed

3408 9 Reporting of Results

3409 *In vitro* methods must be fully documented by good recording and reporting practices and
3410 contain all pertinent details to allow subsequent and adequate analysis and reporting of
3411 results. For example, lot numbers, catalogue numbers, supplier details, and expiry dates for
3412 chemicals and reagents must be listed, as well as temperatures and times (e.g. storage of
3413 chemicals, incubation steps in the *in vitro* method), specific identification of critical
3414 equipment used and, perhaps most importantly, any deviations from SOPs. All this
3415 information must be directly and accurately recorded, signed and dated by the person
3416 performing the activity, as these recordings are important for the correct interpretation of the
3417 results obtained.

3418 Experimental details and results should be easily located; a log page at the front of a
3419 notebook may help track the recordings and observations. Any reference to computer files
3420 containing data should also be catalogued in the notebook. Data files should always be
3421 backed-up in case of computer failure, corruption, or deletion.

3422 Reporting requirements depend on the different development phases of the *in vitro* method.
3423 For regulatory use, requirements for reporting are described in the GLP Principles. Reporting
3424 adequate information and results of all developmental phases will increase the confidence in
3425 the *in vitro* method and would allow for general acceptability by regulatory authorities.

3426 Classical use of *in vitro* method information such as results derived and reported from OECD
3427 *in vitro* test guideline methods (e.g. Ames, *in vitro* skin irritation, *in vitro* eye irritation)
3428 conducted at GLP test facilities are accepted by regulatory authorities of OECD member
3429 states and MAD-adhering economies. In contrast, results derived from non-standard *in vitro*
3430 methods and non-testing methods (QSARs, etc.) are not necessarily accepted by regulators
3431 but might be submitted as supporting information.

3432 There is an increasing tendency towards more transparency when publishing work which
3433 may lead to better reproducibility of published data (Guidelines for Transparency and
3434 Openness Promotion (TOP) Open Science Framework²¹). The EU Competitiveness Council
3435 has also announced their target that all scientific publications resulting from publicly funded
3436 research should be publicly available by 2020²².

3437 Sharing of data in public repositories is also being encouraged and best principles regarding
3438 the publication of scientific data have also been addressed by others, such as the FAIR
3439 (Findable, Accessible, Interoperable, and Reusable) Guiding Principles for scientific data
3440 management and stewardship, by the Nature Publishing Group²³. This initiative not only
3441 promotes more transparency and openness but also promotes the use of computer readable
3442 datasets and data mining so that computers have the ability to access the data autonomously,
3443 unaided by their human operators, which is core to the FAIR Principles.

3444 In addition to the increasing openness and transparency, the publication of negative results is
3445 also gaining more ground e.g. Journal of Negative Results in BioMedicine is an open access,
3446 peer reviewed journal that provides a platform for the publication and discussion of non-
3447 confirmatory and "negative" data.

²¹ <https://osf.io/ud578/>

²² <http://english.eu2016.nl/documents/press-releases/2016/05/27/all-european-scientific-articles-to-be-freely-accessible-by-2020>

²³ <http://www.nature.com/articles/sdata201618>

Field Code Changed

Field Code Changed

Field Code Changed

3448 The information provided by individual *in vitro* methods, as well as *in silico* predictions, can
3449 be combined, interpreted and used for regulatory decision making by means of Integrated
3450 Approaches to Testing and Assessment (IATA) (OECD 2016c). An IATA is an approach that
3451 integrates and weighs all relevant existing evidence and guides the targeted generation of new
3452 data, where required, to build up a hazard or risk assessment acceptable in regulatory
3453 decision-making. Ideally, IATA should be informed by mechanistic understanding of the
3454 underlying toxicokinetics and toxicodynamics. A framework for capturing the toxicodynamic
3455 information is provided by Adverse Outcome Pathways²⁴.

3456 Structured integration of different data types can be performed at different levels, including
3457 raw data and summarised level data. Different levels of data integration can then be used
3458 including Boolean combinations of categorised results, scoring approaches, decision trees,
3459 deterministic and probabilistic approaches. As experience is gained, approaches to data
3460 integration can become standardised. Such approaches, called “defined approaches,” can thus
3461 become core elements of IATA. A defined approach is a formalised decision-making
3462 approach consisting of a fixed data interpretation procedure used to interpret data from a
3463 defined set of information elements (OECD 2006a,b).

3464 When submitting *in vitro* data to a receiving authority, the use of OECD endorsed templates
3465 (e.g. OECD Harmonised Template 201 – Intermediate Effects) is encouraged but is not yet
3466 obligatory. This is facilitated by the use of IUCLID (International Uniform Chemical
3467 Information Database), a software application used to record, store, maintain and exchange
3468 data on intrinsic and hazard properties of chemical substances.

3469 9.1 Publishing

3470 It is essential to have all the results reported in a uniform manner to facilitate their use in the
3471 IATA framework, where the same dataset can be used in many different ways. It is good
3472 practice to publish scientific results in a timely manner. The results will be used and re-used
3473 by other scientists, competitors, modellers or validation study statisticians. Moreover, for any
3474 systemic endpoint the prediction is/will be based on the results of many different studies,
3475 using different methods performed in different facilities.

3476 Therefore, data sharing is encouraged by default, unless there is reason for confidentiality. It
3477 is recommended to not only publish the results, but also the method/SOP (e.g. Nature testing
3478 methods, JOVE, Testing method Exchange, Springer Testing methods, DB-ALM, JRC-
3479 QSAR DB) and using public data-sharing standards and repositories for raw data such as
3480 ISA-TAB (<http://isa-tools.org/>, Dryad, Figshare, and Scientific data
3481 (<http://www.nature.com/sdata/>). In the same vein, assay modifications and further
3482 developments should be published, while tests are still in a development or pre-validation
3483 phase. Such publications should include the changes leading to improvement, the rationale
3484 for them, and this should also entail information on which changes reduce assay performance,
3485 or that do not result in an improvement.

3486 9.2 Reporting of method validation

3487 Validation is at the interface between *in vitro* method development/optimisation and
3488 regulatory acceptance/international recognition and ensures a science-based and
3489 conscientious evaluation of *in vitro* methods and approaches (e.g. Integrated Testing

²⁴ <http://www.oecd.org/chemicalsafety/testing/adverse-outcome-pathways-molecular-screening-and-toxicogenomics.htm>

Field Code Changed

Field Code Changed

Field Code Changed

3490 Strategies, ITS), independent of specific interests, establishing their overall performance and
3491 fitness for a given purpose, i.e. their scientific validity²⁵.
3492

3493 The approach taken by the validation bodies may vary according to the needs of that body as
3494 to whether they will coordinate the validation or whether a validation study should be
3495 submitted to that body for assessment. In general, once a validation study has been
3496 completed, the validation report will be submitted to the validation body for assessment of
3497 the outcome of the validation. While the individual study raw data may not need to be
3498 reported to the validation body, the data analysis and validation report will be required.
3499 Completed validation studies will then undergo independent scientific peer review e.g. by the
3500 EURL ECVAM Scientific Advisory Committee (ESAC), whose outcome will also be made
3501 public.

3502 On completion of the validation study assessment, the recommendations of the validation
3503 body plus the validation report will be made public, e.g. EURL-ECVAM publishes its
3504 EURL-ECVAM Recommendations on its website. Before finalising and publishing these
3505 recommendations, EURL-ECVAM, as a matter of routine, invites comments from the general
3506 public. The aim of an EURL-ECVAM Recommendation is to provide EURL-ECVAM views
3507 on the validity of the *in vitro* method in question, to advise on possible regulatory
3508 applicability, limitations and proper scientific use of the *in vitro* methods, and to suggest
3509 possible follow-up activities in view of addressing knowledge gaps.

3510 The validation study project plan should outline the

- 3511 • Organisations or individuals responsible for data collection
- 3512 • The means of data collection, back-up and archiving
- 3513 • Procedures for data collection (to be established in collaboration with the participating
3514 laboratories
- 3515 • A consistent system of paper-based or electronic labelling of files and folders
3516 including provisions to clearly label final versions.

3517 **9.2.1 Reporting Templates**

3518 To assist in the reporting of data generated in the various laboratories participating in a
3519 validation study, templates are prepared either by the test submitter, the validation body or a
3520 combination of both. These reporting templates should be designed to support later analysis
3521 of the 'high-level' data of validation studies (e.g. by the statistician) and should be validated
3522 before being introduced into the validation process. For a fully GLP-compliant study where
3523 each laboratory is a test facility, the test facility's SOPs and report templates might need to be
3524 used and the validation body subsequently introduces specific templates for analysing the
3525 multi-study validation data.

3526 In case of laboratories operating under GLP, there should be adequate measures put in place
3527 to assure the quality and reliability of the data being reported and recorded. However, in
3528 cases of laboratories that do not operate under quality schemes (e.g. university laboratories) it
3529 is highly advisable to reflect on simple procedures that allow monitoring the adequacy of data
3530 reporting, to ensure that the data used by the statistician for final analysis are correct and
3531 hence relevant. The simplest possible measure would be that data files are checked, at
3532 random, by either a laboratory member or a member of the validation management (the
3533 responsibility should be outlined in the validation project plan). Typical and avoidable

²⁵ <https://eurl-ecvam.jrc.ec.europa.eu/validation-regulatory-acceptance>

3534 mistakes are the incorrect handling (copy / pasting) of numerical values (e.g. wrong column
3535 in an Excel sheet with consequences on automated data analyses), accidental deletion of
3536 embedded formulas or mistakes when normalising values.

3537 Collection of data generated by participating laboratories can be done in a secure
3538 environment by using web-based information exchange tools such as the European
3539 Commission's CIRCABC system.

3540 **9.3 Data reporting for regulatory purposes**

3541 Data captured in GLP studies and results derived thereof will play an important role in
3542 increasing the relevance of *in vitro* data in regulatory contexts. Consideration and ultimately
3543 acceptance of *in vitro* GLP data can be promoted by using a standardised data format. IATA
3544 and AOP knowledge, if properly captured and presented, leads to a better understanding of
3545 toxicity mechanisms, and ultimately the AOP knowledge - derived from testing several
3546 chemicals - may be extrapolated to predict the toxicity of all chemicals that trigger the same
3547 Molecular Initiating Event (MIE) or Key Event (KE). Until recently, the absence of a
3548 template to report Intermediate Effects (like MIEs and KEs) was a limiting factor

3549 The OECD had already designed and published 114 OECD Harmonised Templates (OHTs)
3550 to report test results concerning:

- 3551 • physical/chemical properties (e.g. boiling point, density, flammability, ...),
- 3552 • human toxicity (e.g. carcinogenicity, acute toxicity, ...)
- 3553 • environmental toxicity (e.g. aquatic toxicity, terrestrial toxicity, ...)
- 3554 • other properties describing degradation, accumulation etc.

3555 These templates are geared towards results derived from classical (mostly OECD guideline)
3556 studies, focusing on apical endpoints, i.e. Adverse Outcomes.

3557 However, reporting MIEs or KEs with such a classical OHT would tie them inseparably to
3558 the one Adverse Outcome the one template covers, which is undesirable, as the (*in-vitro*, *in-*
3559 *silico* mechanistic) information is then not easily accessible for building AOPs leading to
3560 other Adverse Outcomes: A Key Event can be relevant not only for one AOP, but several.
3561 Reporting the Intermediate Effect in an "AO-neutral" template makes the data available for
3562 all kinds of AOPs.

3563 A new, AO-neutral OHT was therefore needed that would allow reporting observations from
3564 mechanistic (*in-vitro* and *in-silico*) tests, without immediately locking into one of several
3565 AOs the Intermediate Effect could lead to.

3566 Knowing not only about results of animal tests (classical OHTs), but being able to cross-
3567 reference these test results with the intermediate effect observations (new OHT) has the
3568 potential to lead the way towards a less animal-centred hazard assessment.

3569 The OECD therefore started an initiative to come up with a stable, stakeholder-endorsed
3570 OHT for reporting on "intermediate effects" being observed via *in vitro* assays and possibly
3571 other non-animal test methods (computational predictions etc.). The template was titled
3572 "OHT 201 - Intermediate effects".

3573 OHT 201 was endorsed by the OECD Joint Meeting in 2015 and was finally published in
3574 August 2016, see [http://www.oecd.org/ehs/templates/harmonised-templates-intermediate-](http://www.oecd.org/ehs/templates/harmonised-templates-intermediate-effects.htm)
3575 [effects.htm](http://www.oecd.org/ehs/templates/harmonised-templates-intermediate-effects.htm) for more details.

3576 The basic principle of OHT 201 is that:

- 3577 • one or several objective observation(s) (= results from non-classical test methods)
3578 • lead(s) to one subjective conclusion (= Intermediate Effect present, yes or no).

3579 A properly filled in OHT 201 template therefore conveys a clear statement:

- 3580 • Based on observations $O_1, O_2, \dots O_n$
3581 • a certain chemical
3582 • triggers/does not trigger
3583 • a certain intermediate effect
3584 • on a certain biological level
3585 • at a certain effect concentration.

3586 With OHT 201 being implemented in IUCLID, the ICT system used by industry to fulfil
3587 reporting obligations under more and more legislative programmes (e.g. REACH), the notion
3588 of Intermediate Effects (and implicitly AOPs and predictive toxicology) has started to get
3589 attention in the regulatory world. This is a first step towards the acceptance of results from
3590 alternative tests for regulatory purposes, with the ultimate goal of replacing *in-vivo*-centred
3591 Adverse Outcome observations with alternative-methods-centred IATA/AOP considerations
3592 as the basis for risk assessment.

3593 For *in vitro* methods without a guideline, the Office of Pesticide Programs U.S.
3594 Environmental Protection Agency recommends following OECD Guidance Document 211
3595 (OECD, 2014) for describing non-guideline *in vitro* methods (EPA, 2016).

3596 **10 Storage and retention of records and materials**

3597 As the ultimate goal is to develop an *in vitro* method which will be formally validated for its
3598 future use in a regulatory environment following a quality system (e.g. GLP), it is essential to
3599 have some knowledge of the regulatory requirements specifically relating to the storage and
3600 retention of data, records and materials as the *in vitro* method should be designed so as to be
3601 easily transferrable into a GLP facility.

3602 Data, records and material (specimen) integrity arrangements must be in place throughout the
3603 *in vitro* method lifecycle to ensure that the accuracy and completeness of data, records and
3604 materials are retained. The lifecycle includes all phases in the life of the data, records and
3605 materials, from their initial creation or purchase through processing, use, retention, archival
3606 and retrieval, and eventual destruction. It is vital that formal records used to confirm the
3607 results and how they were obtained are held in a stable/secure form and location which is
3608 documented and traceable and for which there is a minimum storage period. Disposal after
3609 such storage periods should be recorded and a summary report of the stored data and the
3610 means of destruction prepared and held.

3611 Retention arrangements (classified as either archive or backup in case of electronic data)
3612 must be designed to protect data, records and materials from deliberate or accidental changes,
3613 manipulations or deletions thus ensuring integrity throughout the retention period. Archiving
3614 is defined as the long term, permanent retention of completed data and relevant metadata,
3615 records or materials in its final form. Archived data, records or materials may need to be
3616 stored for many years and must be permanently locked such that no changes can be made
3617 without detection. In the case of paper records, storage must be in dry and fire- and light-
3618 proof conditions. In addition to this, they should be easily retrieved for regulatory
3619 inspections.

3620 Storage, retention and archiving requirements vary considerably throughout the lifecycle of
3621 the *in vitro* method. In the early stages, e.g. method development, there are less formal
3622 requirements for storage and retention of records and materials than in the later stages and in
3623 general will follow the university/institute/company policy regarding the storage and
3624 retention of data, records or materials (where such a policy exists). The development phase
3625 should be used to define the raw data, preferably described in the *in vitro* method itself, and
3626 any data (e.g. metadata), records or materials, to be retained when used in a regulatory
3627 environment (e.g. GLP).

3628 As the use of human cells and tissues is critical for the development and implementation of *in*
3629 *vitro* methods for regulatory use in human safety assessment, already in the early stages care
3630 has to be taken regarding their ownership, their identity and who can control their fate. A
3631 number of treaties, laws, and regulations help to guide the ethical collection of human-
3632 derived specimens (Clinchem, n.d.)(<http://www.clinchem.org/content/56/11/1675.full>).

3633 **10.1 Archiving of materials**

3634 An effective archiving system must also provide for the archiving of study samples and
3635 materials, e.g. slides, specimens, test items and reference material under suitable storage
3636 conditions (OECD, 2007a). Specimens and materials should be stored only as long as they
3637 are stable. The Principles of GLP state that: “a sample for analytical purposes from each
3638 batch of test item should be retained for all studies except short-term studies”. The same rules
3639 apply to these archives as apply to the paper based archive, i.e. access restrictions, retrieval
3640 and removal of items, etc.

3641 The storage conditions should be optimal for these samples and often these archives will
3642 require dedicated storage facilities, e.g. low temperature storage such as -20°C, liquid
3643 nitrogen storage or storage of items under inert conditions. Where special storage equipment
3644 is required, the rules governing the control and maintenance of this equipment must be
3645 applied. Where computerised systems are used, these systems must also follow the facility's
3646 policy regarding the use of computerised system, including qualification and validation of
3647 said systems (OECD, 2016b).

3648 Samples of test and reference items or specimen may however be discarded when the quality
3649 of the material no longer permits evaluation. Obviously, the storage conditions should be
3650 optimal for these samples. When samples of test and reference items or specimens are
3651 disposed of before the end of the required retention period, the reason for disposal should be
3652 justified and documented (e.g. the reason might be perishable specimens such as blood
3653 smears, freeze-dried preparations and wet tissues).

3654 **10.2 Document and record management**

3655 Good Document Practices should be followed for both paper records and electronic records
3656 in order to assure data integrity. These principles require that documentation has the
3657 characteristics of being Attributable, Legible, Contemporaneously recorded, Original and
3658 Accurate (sometimes referred to as ALCOA). These essential characteristics apply equally
3659 for both paper and electronic records (WHO, 2016).

3660 When working under a quality system (e.g. ISO, GLP), the management of documents and
3661 records to be retained will be described in the quality system's documentation such as in the
3662 Quality Manual or various SOPs. In a regulatory environment, e.g. GLP, principles and
3663 national laws detail the archiving requirements (OECD, 2007a). Where no formal quality
3664 system exists the general concepts for storage of records still apply.

3665 When storing electronic records it is critical that the relationships between data and metadata
3666 are maintained intact throughout the data life cycle, as the data alone cannot provide the
3667 necessary information to recreate a study.

3668 Electronic records may be stored either as flat files, e.g. individual files, or in a relational
3669 database. Storing data in flat files is relatively trivial to implement, for example files may be
3670 stored on a local drive or on a network share, however the data stored in this way rarely
3671 contains all the relevant metadata, rendering the data incomplete. Due to the nature of flat
3672 files, the user will have full access to the file(s), which presents a greater data integrity risk
3673 since data can be moved, manipulated or even deleted without tracking.

3674 Conversely, a relational database file structure is much more secure as it stores the data and
3675 metadata in different places but maintains the relationship between them. This makes it
3676 inherently more difficult to selectively delete, amend or recreate the original data and the
3677 metadata trail of actions, and most systems provide an audit trail to see who has done what
3678 and when. Such systems mean an upfront investment in software and time (training, etc.), and
3679 the life cycle of the relational database software must also be taken into consideration (e.g.
3680 upgrading, etc.). Once a system has been put in place it is difficult to migrate to another
3681 system and retain all metadata links, specifically the audit trail data. Migrating to a new
3682 system for whatever reason (e.g. current system becomes obsolete) will prove expensive in
3683 cost, time and retraining of personnel.

3684 **10.2.1 Documents and records to be retained**

3685 It is imperative that the historical data, paper-based or in the form of electronic data, are
3686 effectively managed so as to prevent any data integrity issues as this data may be requested
3687 when submitting the method for formal validation.

3688 As compliance with the principles of GLP is required by law for safety studies on chemical
3689 products around the world, it is important that newly developed *in vitro* methods can be
3690 suitable to be performed in a GLP environment from the very beginning of their
3691 development, avoiding an adaption phase at the very end of the process (Coecke *et al.*, 2016).
3692 Ideally, all data supporting the validity of a test method should have been obtained in
3693 accordance with the principles of GLP (OECD, 2005b). Studies which support validations
3694 may or may not be subject to verification depending on compliance monitoring authorities'
3695 programmes²⁶.

3696 In most cases electronic templates (e.g. Excel spreadsheets) will have been developed and
3697 validated for the processing of the study data. These electronic files must be populated with
3698 the acquired data from each laboratory and forwarded to the validation body upon conclusion
3699 of the study. It is recommended that also these electronic files are stored, at least until the
3700 conclusion of the complete validation process. The validation body also has responsibility for
3701 ensuring the integrity of the transfer of these files and their long term storage. It is highly
3702 recommended that the transfer of these files is not performed via email, but using a secure
3703 file transfer system (e.g. https or sftp or similar). The transfer integrity of the files may be
3704 controlled by creating a checksum for each file and checking the checksum at the end of the
3705 transfer steps.

3706 FDA has previously advised that defining paper records as "raw data" (the so called
3707 typewriter rule) does not satisfy the predicate rules, that the industry has misinterpreted the
3708 2003 21 CFR Part 11 Scope and Applications Guidance (FDA, 2003) and that "the printed
3709 paper copy of the chromatogram would not be considered a true copy" ²⁷. Although this
3710 comment was made about chromatographic data, the principles have much wider
3711 implications.

3712 Many electronic records are important to retain in their dynamic (electronic) format, to
3713 enable interaction with the data. Data must be retained in a dynamic form where this is
3714 critical to its integrity or later verification. This should be justified based on risk.

3715 **10.3 Archiving and retention**

3716 In a regulatory environment, when studies are GLP compliant, the archiving retention time is
3717 usually defined in national legislation, however where there is no retention time specified, the
3718 OECD recommends that records and materials should be retained for at least three inspection
3719 cycles so that inspectors can evaluate the compliance of the test facility with the principles of
3720 GLP (OECD, 2007b).

3721 Data is generated during the experimental phase of studies and during this phase the integrity
3722 of the data must be ensured until final archiving of the study. This data will usually be
3723 required for further analysis and as such will not be formally archived until the completion of
3724 the study. It is important that access to this data, both electronic and hard copies, is controlled
3725 until the final archiving upon completion of the study. It is recommended, where possible or

²⁶ <http://www.oecd.org/env/ehs/testing/glp-frequently-asked-questions.htm>

²⁷ <http://www.fda.gov/Drugs/GuidanceComplianceRegulatoryInformation/Guidances/ucm124787.htm>

Field Code Changed

Field Code Changed

3726 feasible, that the electronic data is set as read-only or that an audit trail is provided, detailing
3727 who did what and when.

3728 **10.3.1 Retrieval**

3729 Each facility should have in place procedures concerning the retrieval of archived records and
3730 materials. The procedures should detail who may retrieve records and materials, for how long
3731 and the return of records and materials to the archive. All steps mentioned above need to be
3732 documented and traceable.

3733 In the case of electronic records, viewing the records without the possibility of alteration or
3734 deletion of the archived version does not constitute “retrieval” of a record. Most systems
3735 available nowadays support read-only access, without the possibility to change or delete the
3736 archived record.

3737 **10.3.2 Backup and restore**

3738 When archiving electronic documents, periodic backups should be performed of the
3739 electronic archives. These backups do not constitute archived records, however as they may
3740 be required to be restored in the case a system failure, the same rules regarding access to the
3741 archived electronic records should be applied to access to the backup(s). In general backups
3742 are foreseen for short term storage and not long term storage or archiving and therefore the
3743 long term readability of these archives is usually not an issue, however the restoration of the
3744 backups should also be checked on a regular basis.

3745 Data generated during the experimental phase of the study should also be covered by the
3746 backup and restore policy of the facility.

3747 11 References

- 3748
3749 Allen, M., Bjerke, M., Edlund, H., Nelander, S., Westermark, B., 2016. Origin of the U87MG glioma cell line:
3750 Good news and bad news. *Sci. Transl. Med.* 8, 354re3. doi:10.1126/scitranslmed.aaf6853
- 3751 Anderson, R., O'Hare, M., Balls, M., Brady, M., Brahams, D., Burt, A., Chesné, C., Combes, R., Dennison, A.,
3752 Garthoff, B., Hawksworth, G., Kalter, E., Lechat, A., Mayer, D., Rogiers, V., Sladowski, D., Southee, J.,
3753 Trafford, J., van der Valk, J., van Zeller, A.-M., 1998. The Availability of Human Tissue for Biomedical
3754 Research: The Report and Recommendations of the ECVAM Workshop 32. *Altern. Lab. Anim.* 26, 763–
3755 77.
- 3756 Andrews, P.W., Baker, D., Benvenisty, N., Miranda, B., Bruce, K., Brüstle, O., Choi, M., Choi, Y.-M., Crook,
3757 J.M., de Sousa, P.A., Dvorak, P., Freund, C., Firpo, M., Furue, M.K., Gokhale, P., Ha, H.-Y., Han, E.,
3758 Haupt, S., Healy, L., Hei, D.J., Hovatta, O., Hunt, C., Hwang, S.-M., Inamdar, M.S., Isasi, R.M., Jaconi,
3759 M., Jekerle, V., Kamthorn, P., Kibbey, M.C., Knezevic, I., Knowles, B.B., Koo, S.-K., Laabi, Y.,
3760 Leopoldo, L., Liu, P., Lomax, G.P., Loring, J.F., Ludwig, T.E., Montgomery, K., Mummery, C., Nagy,
3761 A., Nakamura, Y., Nakatsuji, N., Oh, S., Oh, S.-K., Otonkoski, T., Pera, M., Peschanski, M., Pranke, P.,
3762 Rajala, K.M., Rao, M., Ruttachuk, R., Reubinoff, B., Ricco, L., Rooke, H., Sipp, D., Stacey, G.N.,
3763 Suemori, H., Takahashi, T.A., Takada, K., Talib, S., Tannenbaum, S., Yuan, B.-Z., Zeng, F., Zhou, Q.,
3764 2015. Points to consider in the development of seed stocks of pluripotent stem cells for clinical
3765 applications: International Stem Cell Banking Initiative (ISCB). *Regen. Med.* 10, 1–44.
3766 doi:10.2217/rme.14.93
- 3767 Aschner, M., Ceccatelli, S., Daneshian, M., Fritsche, E., Hasiwa, N., Hartung, T., Hogberg, H.T., Leist, M., Li,
3768 A., Mundi, W.R., Padilla, S., Piersma, A.H., Bal-Price, A., Seiler, A., Westerink, R., Zimmer, B., Lein, P.,
3769 2016. Reference compounds for alternative test methods to indicate developmental neurotoxicity (DNT)
3770 potential of chemicals: example lists and criteria for their selection and use. *ALTEX.*
3771 doi:10.14573/altex.1604201
- 3772 ATCC, 2014. *Animal Cell Culture Guide*, Atcc. doi:10.1093/chemse/bjt099
- 3773 Austin, R.P., Barton, P., Cockroft, S.L., Wenlock, M.C., Riley, R.J., 2002. The influence of nonspecific
3774 microsomal binding on apparent intrinsic clearance, and its prediction from physicochemical properties.
3775 *Drug Metab. Dispos.* 30, 1497–503.
- 3776 Bellwon, P., Truisi, G.L., Bois, F.Y., Wilmes, A., Schmidt, T., Savary, C.C., Parmentier, C., Hewitt, P.G.,
3777 Schmal, O., Josse, R., Richert, L., Guillouzo, A., Mueller, S.O., Jennings, P., Testai, E., Dekant, W., 2015.
3778 Kinetics and dynamics of cyclosporine A in three hepatic cell culture systems. *Toxicol. In Vitro* 30, 62–
3779 78. doi:10.1016/j.tiv.2015.07.016
- 3780 Bergström, C., Norinder, U., Luthman, K., Artursson, P., 2002. Experimental and computational screening
3781 models for prediction of aqueous drug solubility. *Pharm. Res.* 19, 182–8.
- 3782 Bessems, J.G., Loizou, G., Krishnan, K., Clewell, H.J., Bernasconi, C., Bois, F., Coecke, S., Collnot, E.-M.,
3783 Diembeck, W., Farcal, L.R., Geraets, L., Gundert-Remy, U., Kramer, N., Küsters, G., Leite, S.B.,
3784 Pelkonen, O.R., Schröder, K., Testai, E., Wilk-Zasadna, I., Zaldivar-Comenges, J.-M., 2014. PBTK
3785 modelling platforms and parameter estimation tools to enable animal-free risk assessment:
3786 recommendations from a joint EPAA--EURL ECVAM ADME workshop. *Regul. Toxicol. Pharmacol.* 68,
3787 119–39. doi:10.1016/j.yrtph.2013.11.008
- 3788 Bilgen, B., Orsini, E., Aaron, R.K., Ciombor, D.M., 2007. FBS suppresses TGF-beta1-induced chondrogenesis
3789 in synoviocyte pellet cultures while dexamethasone and dynamic stimuli are beneficial. *J. Tissue Eng.*
3790 *Regen. Med.* 1, 436–42. doi:10.1002/term.56
- 3791 Bosgra, S., Westerhout, J., 2015. Interpreting in vitro developmental toxicity test battery results: The
3792 consideration of toxicokinetics. *Reprod. Toxicol.* 55, 73–80. doi:10.1016/j.reprotox.2014.11.001
- 3793 Bouhifd, M., Beger, R., Flynn, T., Guo, L., Harris, G., Hogberg, H., Kaddurah-Daouk, R., Kamp, H.,
3794 Kleensang, A., Maertens, A., Odwin-DaCosta, S., Pamies, D., Robertson, D., Smirnova, L., Sun, J., Zhao,
3795 L., Hartung, T., 2015. Quality assurance of metabolomics. *ALTEX* 32, 319–26.
- 3796 Box, G.E., Hunter, J., Hunter, W.G., 2005. *Introduction to fractional factorial experimentation*, in: *Statistics for*
3797 *Experimenters*. Wiley, New York.
- 3798 Brindley, D.A., Davie, N.L., Culme-Seymour, E.J., Mason, C., Smith, D.W., Rowley, J.A., 2012. Peak serum:
3799 implications of serum supply for cell therapy manufacturing. *Regen. Med.* 7, 7–13.
3800 doi:10.2217/rme.11.112
- 3801 Broeders, J.J.W., Blaauboer, B.J., Hermens, J.L.M., 2011. Development of a negligible depletion-solid phase
3802 microextraction method to determine the free concentration of chlorpromazine in aqueous samples
3803 containing albumin. *J. Chromatogr. A* 1218, 8529–35. doi:10.1016/j.chroma.2011.09.064
- 3804 Burd, E.M., 2010. Validation of laboratory-developed molecular assays for infectious diseases. *Clin. Microbiol.*
3805 *Rev.* 23, 550–76. doi:10.1128/CMR.00074-09

3806 Cadena-Herrera, D., Esparza-De Lara, J.E., Ramírez-Ibañez, N.D., López-Morales, C.A., Pérez, N.O., Flores-
3807 Ortiz, L.F., Medina-Rivero, E., 2015. Validation of three viable-cell counting methods: Manual, semi-
3808 automated, and automated. *Biotechnol. Reports* 7, 9–16. doi:10.1016/j.btre.2015.04.004
3809 Chen, Y., Stevens, B., Chang, J., Milbrandt, J., Barres, B.A., Hell, J.W., 2008. NS21: Re-defined and modified
3810 supplement B27 for neuronal cultures. *J. Neurosci. Methods* 171, 239–247.
3811 doi:10.1016/j.jneumeth.2008.03.013
3812 Christopoulos, A., 1998. Assessing the distribution of parameters in models of ligand-receptor interaction: to
3813 log or not to log. *Trends Pharmacol. Sci.* 19, 351–7.
3814 Clinchem, n.d. Clinchem [WWW Document]. URL <http://www.clinchem.org/content/56/11/1675.full>
3815 Coecke, S., Balls, M., Bowe, G., Davis, J., Gstraunthaler, G., Hartung, T., Hay, R., Merten, O.-W., Price, A.,
3816 Schechtman, L., Stacey, G., Stokes, W., 2005. Guidance on good cell culture practice. a report of the
3817 second ECVAM task force on good cell culture practice. *Altern. Lab. Anim.* 33, 261–87.
3818 Coecke, S., Bernasconi, C., Bowe, G., Bostroem, A.-C., Burton, J., Cole, T., Fortaner, S., Gouliarmou, V., Gray,
3819 A., Griesinger, C., Louhimies, S., Gyves, E.M., Joossens, E., Prinz, M.-J., Milcamps, A., Parissis, N.,
3820 Wilk-Zasadna, I., Barroso, J., Desprez, B., Langezaal, I., Liska, R., Morath, S., Reina, V., Zorzoli, C.,
3821 Zuang, V., 2016. Practical Aspects of Designing and Conducting Validation Studies Involving Multi-
3822 study Trials. *Adv. Exp. Med. Biol.* 856, 133–163. doi:10.1007/978-3-319-33826-2_5
3823 Coecke, S., Bowe, G., Milcamps, A., Bernasconi, C., Bostroem, A.-C., Bories, G., Fortaner, S., Gineste, J.-M.,
3824 Gouliarmou, V., Langezaal, I., Liska, R., Mendoza, E., Morath, S., Reina, V., Wilk-Zasadna, I., Whelan,
3825 M., 2014. Considerations in the development of in vitro toxicity testing methods intended for regulatory
3826 use, in: Bal-Price, A., Jennings, P. (Eds.), *In Vitro Toxicology Systems, Methods in Pharmacology and*
3827 *Toxicology*. Springer New York. doi:10.1007/978-1-4939-0521-8_25
3828 Crawley, M.J., 2015. Chapter 1. Fundamentals, in: *Statistics: An Introduction Using R*, 2nd Edition. Wiley, pp.
3829 8–10.
3830 Crean, D., Bellwon, P., Aschauer, L., Limonciel, A., Moenks, K., Hewitt, P., Schmidt, T., Herrgen, K., Dekant,
3831 W., Lukas, A., Bois, F., Wilmes, A., Jennings, P., Leonard, M.O., 2015. Development of an in vitro renal
3832 epithelial disease state model for xenobiotic toxicity testing. *Toxicol. In Vitro* 30, 128–37.
3833 doi:10.1016/j.tiv.2014.11.015
3834 Eastwood, B.J., Farmen, M.W., Iversen, P.W., Craft, T.J., Smallwood, J.K., Garbison, K.E., Delapp, N.W.,
3835 Smith, G.F., 2006. The minimum significant ratio: a statistical parameter to characterize the
3836 reproducibility of potency estimates from concentration-response assays and estimation by replicate-
3837 experiment studies. *J. Biomol. Screen.* 11, 253–61. doi:10.1177/1087057105285611
3838 ECACC, 2010. *Fundamentals Techniques in Cell Culture*.
3839 Ekins, S., Olechno, J., Williams, A.J., 2013. Dispensing processes impact apparent biological activity as
3840 determined by computational and statistical analyses. *PLoS One* 8, e62325.
3841 doi:10.1371/journal.pone.0062325
3842 EMA, 2014. Concept paper on review and update of EMA guidelines to implement best practice with regard to
3843 3Rs (replacement, reduction and refinement) in regulatory testing of medicinal products.
3844 EMEA, 2011. Guideline on bioanalytical method validation EMEA/CHMP/EWP/192217/2009.
3845 EPA, 2016. Process for Evaluating & Implementing Alternative Approaches To Traditional in Vivo Acute
3846 Toxicity Studies for FIFRA Regulatory Use 1–9.
3847 EPA, 2007. *Guidance for Preparing Standard Operating Procedures*.
3848 ESAC, 2008. ESAC statement on the use of FCS and other animal-derived supplements.
3849 Esch, P.M., Moor, C., Schmid, B., Albertini, S., Hassler, S., Donzé, G., Saxer, H.P., 2010. Good Laboratory
3850 Practice (GLP) - Guidelines for the Development and Validation of Spreadsheets. *Qual. Assur. J.* 13, 41–
3851 56. doi:10.1002/qaj.466
3852 EU, 2009. DIRECTIVE 2009/41/EC OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL on the
3853 contained use of genetically modified micro-organisms.
3854 FDA, 2003. Part 11, Electronic Records; Electronic Signatures — Scope and Application. *FDA Guid. Ind.* 12.
3855 FDA, 2001. *Guidance for Industry Bioanalytical Method Validation*.
3856 Frattini, A., Fabbri, M., Valli, R., De Paoli, E., Montalbano, G., Gribaldo, L., Pasquali, F., Maserati, E., 2015.
3857 High variability of genomic instability and gene expression profiling in different HeLa clones. *Sci. Rep.* 5,
3858 15377. doi:10.1038/srep15377
3859 Gaddis, G.M., Gaddis, M.L., 1990. Introduction to biostatistics: Part 3, Sensitivity, specificity, predictive value,
3860 and hypothesis testing. *Ann. Emerg. Med.* 19, 591–7. doi:10.1016/S0196-0644(05)82198-5
3861 Geraghty, R.J., Capes-Davis, A., Davis, J.M., Downward, J., Freshney, R.I., Knezevic, I., Lovell-Badge, R.,
3862 Masters, J.R.W., Meredith, J., Stacey, G.N., Thraves, P., Vias, M., 2014. Guidelines for the use of cell
3863 lines in biomedical research. *Br. J. Cancer* 111, 1021–46. doi:10.1038/bjc.2014.166
3864 Grant, R.J., Roberts, K., Pointon, C., Hodgson, C., Womersley, L., Jones, D.C., Tang, E., 2009. Achieving
3865 accurate compound concentration in cell-based screening: validation of acoustic droplet ejection

3866 technology. *J. Biomol. Screen.* 14, 452–9. doi:10.1177/1087057109336588

3867 Groothuis, F.A., Heringa, M.B., Nicol, B., Hermens, J.L.M., Blaauboer, B.J., Kramer, N.L., 2015. Dose metric

3868 considerations in in vitro assays to improve quantitative in vitro-in vivo dose extrapolations. *Toxicology*

3869 332, 30–40. doi:10.1016/j.tox.2013.08.012

3870 Groten, J.P., Schoen, E.D., van Bladeren, P.J., Kuper, C.F., van Zorge, J.A., Feron, V.J., 1997. Subacute toxicity

3871 of a mixture of nine chemicals in rats: detecting interactive effects with a fractionated two-level factorial

3872 design. *Fundam. Appl. Toxicol.* 36, 15–29. doi:10.1006/faat.1996.2281

3873 Gülден, M., Jess, A., Kammann, J., Maser, E., Seibert, H., 2010. Cytotoxic potency of H₂O₂ in cell cultures:

3874 impact of cell concentration and exposure time. *Free Radic. Biol. Med.* 49, 1298–305.

3875 doi:10.1016/j.freeradbiomed.2010.07.015

3876 Gülден, M., Kähler, D., Seibert, H., 2015. Incipient cytotoxicity: A time-independent measure of cytotoxic

3877 potency in vitro. *Toxicology* 335, 35–45. doi:10.1016/j.tox.2015.07.002

3878 Gülден, M., Mörchel, S., Seibert, H., 2001. Factors influencing nominal effective concentrations of chemical

3879 compounds in vitro: cell concentration. *Toxicol. In Vitro* 15, 233–43.

3880 Gülден, M., Seibert, H., 2003. In vitro-in vivo extrapolation: estimation of human serum concentrations of

3881 chemicals equivalent to cytotoxic concentrations in vitro. *Toxicology* 189, 211–22.

3882 Gunetti, M., Castiglia, S., Rustichelli, D., Mareschi, K., Sanavio, F., Muraro, M., Signorino, E., Castello, L.,

3883 Ferrero, I., Fagioli, F., 2012. Validation of analytical methods in GMP: the disposable Fast Read 102®

3884 device, an alternative practical approach for cell counting. *J. Transl. Med.* 10, 112. doi:10.1186/1479-

3885 5876-10-112

3886 Gupta, K., Rispin, A., Stitzel, K., Coecke, S., Harbell, J., 2005. Ensuring quality of in vitro alternative test

3887 methods: issues and answers. *Regul. Toxicol. Pharmacol.* 43, 219–24. doi:10.1016/j.yrtph.2005.03.010

3888 Haas, J. V., Eastwood, B.J., Iversen, P.W., Weidner, J.R., 2004. Minimum Significant Ratio – A Statistic to

3889 Assess Assay Variability, Assay Guidance Manual. Eli Lilly & Company and the National Center for

3890 Advancing Translational Sciences.

3891 Hall, M.D., Telma, K.A., Chang, K.E., Lee, T.D., Madigan, J.P., Lloyd, J.R., Goldlust, I.S., Hoeschele, J.D.,

3892 Gottesman, M.M., 2014. Say no to DMSO: Dimethylsulfoxide inactivates cisplatin, carboplatin, and other

3893 platinum complexes. *Cancer Res.* 74, 3913–3922. doi:10.1158/0008-5472.CAN-14-0247

3894 Harbell, J., Raabe, H., J. Harbell and H. Raabe, 2014. In Vitro Methods for the Prediction of Ocular and Dermal

3895 Toxicity, in: *Handbook of Toxicology, Third Edition.* CRC Press, pp. 197–231. doi:doi:10.1201/b16632-6

3896 Hartung, T., Bremer, S., Casati, S., Coecke, S., Corvi, R., Fortaner, S., Gribaldo, L., Halder, M., Hoffmann, S.,

3897 Roi, A.J., Prieto, P., Sabbioni, E., Scott, L., Worth, A., Zuang, V., 2004. A modular approach to the

3898 ECVAM principles on test validity. *Altern. Lab. Anim.* 32, 467–72.

3899 Hayashi, M., Dearfield, K., Kasper, P., Lovell, D., Martus, H.-J., Thybaud, V., 2011. Compilation and use of

3900 genetic toxicity historical control data. *Mutat. Res.* 723, 87–90. doi:10.1016/j.mrgentox.2010.09.007

3901 Heringa, M.B., de Wit-Bos, L., Bos, P.M.J., Hakkert, B., 2014. Do current EU regulations for the safety

3902 assessment of chemical substances pose legal barriers for the use of alternatives to animal testing?, RIVM

3903 Report.

3904 Heringa, M.B., Hermens, J.L.M., 2003. Measurement of free concentrations using negligible depletion-solid

3905 phase microextraction (nd-SPME). *TrAC Trends Anal. Chem.* 22, 575–587. doi:10.1016/S0165-

3906 9936(03)01006-9

3907 Heringa, M.B., Hogevoender, C., Busser, F., Hermens, J.L.M., 2006. Measurement of the free concentration of

3908 octylphenol in biological samples with negligible depletion-solid phase microextraction (nd-SPME):

3909 analysis of matrix effects. *J. Chromatogr. B. Analyt. Technol. Biomed. Life Sci.* 834, 35–41.

3910 doi:10.1016/j.jchromb.2006.02.009

3911 Heringa, M.B., Schreurs, R.H.M.M., Busser, F., van der Saag, P.T., van der Burg, B., Hermens, J.L.M., 2004.

3912 Toward more useful in vitro toxicity data with measured free concentrations. *Environ. Sci. Technol.* 38,

3913 6263–70.

3914 Hoelke, B., Gieringer, S., Arlt, M., Saal, C., 2009. Comparison of nephelometric, UV-spectroscopic, and HPLC

3915 methods for high-throughput determination of aqueous drug solubility in microtiter plates. *Anal. Chem.*

3916 81, 3165–72. doi:10.1021/ac9000089

3917 ISCB, 2009. Consensus guidance for banking and supply of human embryonic stem cell lines for research

3918 purposes. *Stem Cell Rev.* 5, 301–14. doi:10.1007/s12015-009-9085-x

3919 ISO, 2015. ISO 15519-2:2015 - Specifications for diagrams for process industry -- Part 2: Measurement and

3920 control.

3921 ISO, 2005. ISO/IEC 17025:2005 - General requirements for the competence of testing and calibration

3922 laboratories.

3923 ISPE, 2012. ISPE GAMP Good Practice Guide. A Risk-Based Approach to GxP Compliant Laboratory

3924 Computerized Systems, Second. ed. ISPE Publication.

3925 Iversen, P.W., Beck, B., Chen, Y.-F., Dere, W., Devanarayan, V., Eastwood, B.J., Farmen, M.W., Iturria, S.J.,

3926 Montrose, C., Moore, R.A., Weidner, J.R., Sittampalam, G.S., 2004. HTS Assay Validation, Assay
3927 Guidance Manual. Eli Lilly & Company and the National Center for Advancing Translational Sciences.

3928 Jochems, C.E.A., van der Valk, J.B.F., Stafleu, F.R., Baumans, V., n.d. The use of fetal bovine serum: ethical or
3929 scientific problem? *Altern. Lab. Anim.* 30, 219–27.

3930 Jouyban, A., Fakhree, M.A., 2012. Experimental and Computational Methods Pertaining to Drug Solubility, in:
3931 Toxicity and Drug Testing. InTech. doi:10.5772/19206

3932 Kleensang, A., Vantangoli, M.M., Odwin-DaCosta, S., Andersen, M.E., Boekelheide, K., Bouhifd, M., Fornace,
3933 A.J., Li, H.-H., Livi, C.B., Madnick, S., Maertens, A., Rosenberg, M., Yager, J.D., Zhaog, L., Hartung, T.,
3934 2016. Genetic variability in a frozen batch of MCF-7 cells invisible in routine authentication affecting cell
3935 function. *Sci. Rep.* 6, 28994. doi:10.1038/srep28994

3936 Knöbel, M., Busser, F.J.M., Rico-Rico, A., Kramer, N.I., Hermens, J.L.M., Hafner, C., Tanneberger, K.,
3937 Schirmer, K., Scholz, S., 2012. Predicting adult fish acute lethality with the zebrafish embryo: relevance
3938 of test duration, endpoints, compound properties, and exposure concentration analysis. *Environ. Sci.*
3939 *Technol.* 46, 9690–700. doi:10.1021/es301729q

3940 Kramer, J.R., Schmidt, N.W., Mayle, K.M., Kamei, D.T., Wong, G.C.L., Deming, T.J., 2015. Reinventing Cell
3941 Penetrating Peptides Using Glycosylated Methionine Sulfonium Ion Sequences. *ACS Cent. Sci.* 1, 83–8.
3942 doi:10.1021/acscentsci.5b00054

3943 Kramer, N.I., Krismartina, M., Rico-Rico, A., Blaauboer, B.J., Hermens, J.L.M., 2012. Quantifying processes
3944 determining the free concentration of phenanthrene in Basal cytotoxicity assays. *Chem. Res. Toxicol.* 25,
3945 436–45. doi:10.1021/tx200479k

3946 Krzywinski, M., Altman, N., 2013. Points of significance: Power and sample size. *Nat. Methods* 10, 1139–1140.
3947 doi:10.1038/nmeth.2738

3948 Leist, M., Efreanova, L., Karreman, C., 2010. Food for thought ... considerations and guidelines for basic test
3949 method descriptions in toxicology. *ALTEX* 27, 309–17.

3950 Loizou, G., Spendiff, M., Barton, H.A., Bessems, J., Bois, F.Y., d'Yvoire, M.B., Buist, H., Clewell, H.J., Meek,
3951 B., Gundert-Remy, U., Goerlitz, G., Schmitt, W., 2008. Development of good modelling practice for
3952 physiologically based pharmacokinetic models for use in risk assessment: the first steps. *Regul. Toxicol.*
3953 *Pharmacol.* 50, 400–11. doi:10.1016/j.yrtph.2008.01.011

3954 Louise, J., de Jong, E., van de Sandt, J.J.M., Blaauboer, B.J., Woutersen, R.A., Piersma, A.H., Rietjens,
3955 I.M.C.M., Verwei, M., 2010. The use of in vitro toxicity data and physiologically based kinetic modeling
3956 to predict dose-response curves for in vivo developmental toxicity of glycol ethers in rat and man.
3957 *Toxicol. Sci.* 118, 470–84. doi:10.1093/toxsci/kfq270

3958 Motulsky, H.J., Brown, R.E., 2006. Detecting outliers when fitting data with nonlinear regression - a new
3959 method based on robust nonlinear regression and the false discovery rate. *BMC Bioinformatics* 7, 123.
3960 doi:10.1186/1471-2105-7-123

3961 Motulsky, H.J., Christopoulos, A., 2004. Fitting models to biological data using linear and nonlinear regression.
3962 National Institutes of Health, 2001. C & D of NIH Publication No: 01-4500. Guidance Document on Using In
3963 Vitro Data to Estimate In Vivo Starting Doses for Acute Toxicity 1–102.

3964 OECD, 2016a. OECD principles No 18 - OECD Position Paper Regarding the Relationship between the OECD
3965 Principles of GLP and ISO/IEC 17025.

3966 OECD, 2016b. OECD principles No 17 - Application of GLP Principles to Computerised Systems, OECD.
3967 doi:ENV/JM/MONO(2007)10

3968 OECD, 2014. OECD Series on Testing and Assessment. Number 211. Guidance Document for describing non-
3969 guideline in vitro test methods. doi:ENV/JM/MONO(2014)35

3970 OECD, 2011. OECD guideline for testing of chemicals. Section 2: effects on biotic systems. In: Guideline 201:
3971 Freshwater Alga and Cyanobacteria, Growth Inhibition Test. Organisation for Economic Cooperation and
3972 Development, Paris, France.

3973 OECD, 2009. OECD Guidelines for the Testing of Chemicals, Section 4, Test No. 437: Bovine Corneal Opacity
3974 and Permeability Test Method for Identifying Ocular Corrosives and Severe Irritants.

3975 OECD, 2007a. OECD series on Principles of Good Laboratory Practice and Compliance Monitoring Number
3976 15. Establishment and Control of Archives that Operate in Compliance with the Principles of GLP. Paris :
3977 OECD Publishing.

3978 OECD, 2007b. OECD GLP 15 - Establishment and Control of Archives that Operate in Compliance with the
3979 Principles of GLP. Guideline. doi:ENV/JM/MONO(2007)10

3980 OECD, 2006a. OECD Guideline for testing of chemicals. Draft proposal for a new guideline. In: Fish Embryo
3981 Toxicity (FET) Test. Organisation for Economic Cooperation and Development, Paris, France.

3982 OECD, 2006b. OECD guideline for testing of chemicals. Section 2: Effects on biotic systems. In: Guideline
3983 221: Lemna sp. Growth Inhibition Test. Organisation for Economic Cooperation and Development, Paris,
3984 France.

3985 OECD, 2005a. OECD Series on Testing and Assessment. Number 34. Guidance Document on the Validation

3986 and International Acceptance of New or Updated Test Methods for Hazard Assessment.
3987 ENV/JM/MONO(2005)14.
3988 OECD, 2005b. OECD Guidance Document on the Validation and International Acceptance of new or updated
3989 test methods for Hazard Assessment No. 34. Guidance 33, 1–96. doi:ENV/JM/MONO(2005)14
3990 OECD, 2004a. The Application of the Principles of GLP to in vitro Studies. OECD Publishing.
3991 OECD, 2004b. OECD series on Principles of Good Laboratory Practice and Compliance Monitoring. Number
3992 14. Advisory Document of the Working Group on Good Laboratory Practice: The Application of the
3993 Principles of GLP to in vitro Studies. ENV/JM/MONO(2004)26, OECD Series on Principles of Good
3994 Laboratory Practice and Compliance Monitoring. OECD Publishing. doi:10.1787/9789264084971-en
3995 OECD, 2004c. OECD GLP 14 - Application of GLP to in vitro Studies. Guideline.
3996 doi:ENV/JM/MONO(2004)26
3997 OECD, 2000. Guidance document on aquatic toxicity testing of difficult substances and mixtures. Environ.
3998 Heal. Saf. Publ. Ser. Test. Assess. - N°23 53. doi:ENV/JM/MONO(2007)10
3999 OECD, 1999. OECD Series on Principles of GLP and Compliance Monitoring Number 4 (Revised) -
4000 QUALITY ASSURANCE AND GLP. Guideline 4. doi:ENV/JM/MONO(99)20
4001 OECD, 1998a. OECD series on principles of good laboratory practice and compliance monitoring. Number 1.
4002 OECD principles on good laboratory practice. ENV/MC/CHEM(98)17. Guideline 33, 1–172.
4003 OECD, 1998b. OECD Principles on Good Laboratory Practice.
4004 OMCL, 2011. OMCL Network of the Council of Europe QUALITY ASSURANCE DOCUMENT.
4005 QUALIFICATION OF EQUIPMENT CORE DOCUMENT PA/PH/OMCL (08) 73 2R.
4006 OMCL, 2009. Validation of Computerised Systems Annex 1: Validation of computerised calculation systems:
4007 example of validation of in-house software PA/PH/OMCL (08) 87 2R.
4008 Ono, K., Satoh, M., Yoshida, T., Ozawa, Y., Kohara, A., Takeuchi, M., Mizusawa, H., Sawada, H., 2007.
4009 Species identification of animal cells by nested PCR targeted to mitochondrial DNA. *In Vitro Cell. Dev.*
4010 *Biol. Anim.* 43, 168–75. doi:10.1007/s11626-007-9033-5
4011 Pamies, D., Bal-Price, A., Simeonov, A., Tagle, D., Allen, D., Gerhold, D., Yin, D., Pistollato, F., Inutsuka, T.,
4012 Sullivan, K., Stacey, G., Salem, H., Leist, M., Daneshian, M., Vemuri, M.C., McFarland, R., Coecke, S.,
4013 Fitzpatrick, S.C., Lakshmipathy, U., Mack, A., Wang, W.B., Daiju, Y., Sekino, Y., Kanda, Y., Smirnova,
4014 L., Hartung, T., 2016. Good Cell Culture Practice for stem cells and stem-cell-derived models. *ALTEX*.
4015 doi:10.14573/altex.1607121
4016 Persson, L.C., Porter, C.J.H., Charman, W.N., Bergström, C.A.S., 2013. Computational prediction of drug
4017 solubility in lipid based formulation excipients. *Pharm. Res.* 30, 3225–3237. doi:10.1007/s11095-013-
4018 1083-7
4019 Phelan, M., Lawler, G., 2001. Cell counting. *Curr. Protoc. Cytom.* 00:3A:A.3A.
4020 Pincus, R., 1995. Barnett, V., and Lewis T.: *Outliers in Statistical Data*. 3rd edition. J. Wiley & Sons 1994,
4021 XVII. 582 pp. *Biometrical J.* 37, 256–256. doi:10.1002/bimj.4710370219
4022 Pistollato, F., Bremer-Hoffmann, S., Healy, L., Young, L., Stacey, G., 2012. Standardization of pluripotent stem
4023 cell cultures for toxicity testing. *Expert Opin. Drug Metab. Toxicol.* 8, 239–57.
4024 doi:10.1517/17425255.2012.639763
4025 Pistollato, F., Louise, J., Scelfo, B., Mennecozzi, M., Accordi, B., Basso, G., Gaspar, J.A., Zagoura, D.,
4026 Barilari, M., Palosaari, T., Sachinidis, A., Bremer-Hoffmann, S., 2014. Development of a pluripotent stem
4027 cell derived neuronal model to identify chemically induced pathway perturbations in relation to
4028 neurotoxicity: effects of CREB pathway inhibition. *Toxicol. Appl. Pharmacol.* 280, 378–88.
4029 doi:10.1016/j.taap.2014.08.007
4030 Pomponio, G., Savary, C.C., Parmentier, C., Bois, F., Guillouzo, A., Romanelli, L., Richert, L., Di Consiglio,
4031 E., Testai, E., 2015. In vitro kinetics of amiodarone and its major metabolite in two human liver cell
4032 models after acute and repeated treatments. *Toxicol Vitro.* 30, 36–51. doi:10.1016/j.tiv.2014.12.012
4033 Reinert, K.H., Giddings, J.M., Judd, L., 2002. Effects analysis of time-varying or repeated exposures in aquatic
4034 ecological risk assessment of agrochemicals. *Environ. Toxicol. Chem.* 21, 1977–1992.
4035 doi:10.1002/etc.5620210928
4036 Riedl, J., Altenburger, R., 2007. Physicochemical substance properties as indicators for unreliable exposure in
4037 microplate-based bioassays. *Chemosphere* 67, 2210–20. doi:10.1016/j.chemosphere.2006.12.022
4038 Rispin, A., Stitzel, K., Harbell, J., Klausner, M., 2006. Ensuring quality of in vitro alternative test methods:
4039 Current practice. *Regul. Toxicol. Pharmacol.* 45, 97–103. doi:10.1016/j.yrtph.2005.03.005
4040 Scenhir, 2015. Final Opinion on the Guidance on the Determination of Potential Health Effects of
4041 Nanomaterials Used in Medical Devices. doi:10.2772/41391
4042 Schmidt, B.Z., Lehmann, M., Gutbier, S., Nembo, E., Noel, S., Smirnova, L., Forsby, A., Hescheler, J., Avci,
4043 H.X., Hartung, T., Leist, M., Kobilák, J., Dinnyés, A., 2016. In vitro acute and developmental
4044 neurotoxicity screening: an overview of cellular platforms and high-throughput technical possibilities.
4045 *Arch. Toxicol.* doi:10.1007/s00204-016-1805-9

4046 Shahdadfar, A., Frønsdal, K., Haug, T., Reinholt, F.P., Brinchmann, J.E., 2005. In vitro expansion of human
4047 mesenchymal stem cells: choice of serum is a determinant of cell proliferation, differentiation, gene
4048 expression, and transcriptome stability. *Stem Cells* 23, 1357–66. doi:10.1634/stemcells.2005-0094
4049 Sittampalam, G.S., 1997. Design of Signal Windows in High Throughput Screening Assays for Drug Discovery.
4050 *J. Biomol. Screen.* 2, 159–169. doi:10.1177/108705719700200306
4051 Smith, K.E.C., Heringa, M.B., Uytewaal, M., Mayer, P., 2013. The dosing determines mutagenicity of
4052 hydrophobic compounds in the Ames II assay with metabolic transformation: passive dosing versus
4053 solvent spiking. *Mutat. Res.* 750, 12–8. doi:10.1016/j.mrgentox.2012.07.006
4054 Smith, K.E.C., Oostingh, G.J., Mayer, P., 2010. Passive dosing for producing defined and constant exposure of
4055 hydrophobic organic compounds during in vitro toxicity tests. *Chem. Res. Toxicol.* 23, 55–65.
4056 doi:10.1021/tx900274j
4057 Stacey, G. (Glyn), Davis, J., 2007. *Medicines from animal cell culture*. Wiley.
4058 Stacey, G., Coecke, S., Price, A., Healy, L., Jennings, P., Wilmes, A., Pinset, C., Sundstrom, M., Myatt, G.,
4059 2016. Ensuring the Quality of Stem Cell-Derived In Vitro Models for Toxicity Testing, in: Eskes, C.,
4060 Whelan, M. (Eds.), *Validation of Alternative Methods for Toxicity Testing*. Springer, pp. 259–297.
4061 Stacey, G.N., 2011. Cell culture contamination. *Methods Mol. Biol.* 731, 79–91. doi:10.1007/978-1-61779-080-
4062 5_7
4063 Stacey, G.N., Hartung, T., 2006. The use of human tissue in drug discovery, scientific, ethical, legal and
4064 regulatory environments, Chapter 9., in: V, M.U. and S. (Ed.), *Drug Testing In Vitro – Breakthroughs and*
4065 *Trends in Cell Culture Technology*. Wiley-VCH, pp. 231–250.
4066 Sullivan, B., 2001. Assay Development in High Density MicroWell® Plates: Use of Well Geometries, Format,
4067 Surface Modification and Optical Properties to Achieve Optimal Assay Performance. *J. Assoc. Lab.*
4068 *Autom.* 6, 47–52. doi:10.1016/S1535-5535(04)00125-X
4069 Tanneberger, K., Rico-Rico, A., Kramer, N.I., Busser, F.J.M., Hermens, J.L.M., Schirmer, K., 2010. Effects of
4070 solvents and dosing procedure on chemical toxicity in cell-based in vitro assays. *Environ. Sci. Technol.*
4071 44, 4775–81. doi:10.1021/es100045y
4072 Ulrey, A.K., Curren, R.D., Harbell, J.W., Mun, G., A, H., Raabe, 2005. Applying Good Laboratory Practices
4073 (GLPs) to In Vitro Studies: One Laboratory’s Perspective.
4074 Usta, S.N., Scharer, C.D., Xu, J., Frey, T.K., Nash, R.J., 2014. Chemically defined serum-free and xeno-free
4075 media for multiple cell lineages. *Ann. Transl. Med.* 2, 97. doi:10.3978/j.issn.2305-5839.2014.09.05
4076 van der Valk, J., Mellor, D., Brands, R., Fischer, R., Gruber, F., Gstraunthaler, G., Hellebrekers, L., Hyllner, J.,
4077 Jonker, F.H., Prieto, P., Thalén, M., Baumans, V., 2004. The humane collection of fetal bovine serum and
4078 possibilities for serum-free cell and tissue culture. *Toxicol. In Vitro* 18, 1–12.
4079 Viswanathan, C.T., Bansal, S., Booth, B., DeStefano, A.J., Rose, M.J., Sailstad, J., Shah, V.P., Skelly, J.P.,
4080 Swann, P.G., Weiner, R., 2007. Quantitative bioanalytical methods validation and implementation: Best
4081 practices for chromatographic and ligand binding assays. *Pharm. Res.* doi:10.1007/s11095-007-9291-7
4082 Vogel, G., 2010. Cell biology. To scientists’ dismay, mixed-up cell lines strike again. *Science* 329, 1004.
4083 doi:10.1126/science.329.5995.1004
4084 Weil, C.S., Scala, R.A., 1971. Study of intra- and interlaboratory variability in the results of rabbit eye and skin
4085 irritation tests. *Toxicol. Appl. Pharmacol.* 19, 276–360.
4086 WHO, 2016. Annex 5 Guidance on good data and record management practices, in: WHO Technical Report
4087 Series, No. 996, 2016, Annex 5. pp. 165–210.
4088 WHO, 2013. WHO Handbook: Quality practices in basic biomedical research, WHO. World Health
4089 Organization.
4090 WHO, 2010. Recommendations for the evaluation of animal cell cultures as substrates for the manufacture of
4091 biological medicinal products and for the characterization of cell banks Table of contents, Who.
4092 Wilson, H.K., Canfield, S.G., Hjortness, M.K., Palecek, S.P., Shusta, E. V., 2015. Exploring the effects of cell
4093 seeding density on the differentiation of human pluripotent stem cells to brain microvascular endothelial
4094 cells. *Fluids Barriers CNS* 12, 13. doi:10.1186/s12987-015-0007-9
4095 Yamasaki, S., Taguchi, Y., Shimamoto, A., Mukasa, H., Tahara, H., Okamoto, T., 2014. Generation of human
4096 induced pluripotent stem (Ips) cells in serum- and feeder-free defined culture and TGF-B1 regulation of
4097 pluripotency. *PLoS One* 9, e87151. doi:10.1371/journal.pone.0087151
4098 Yoshimura, I., Matsumoto, K., 1994. Notes on the use of historical controls. *Environ. Health Perspect.* 19–23.
4099 Zeiger, E., Anderson, B., Haworth, S., Lawlor, T., Mortelmans, K., 1988. Salmonella mutagenicity tests: IV.
4100 Results from the testing of 300 chemicals. *Environ. Mol. Mutagen.* 11 Suppl 1, 1–157.
4101 Zhang, J.J.-H., Chung, T., Oldenburg, K., 1999. A Simple Statistical Parameter for Use in Evaluation and
4102 Validation of High Throughput Screening Assays. *J. Biomol. Screen.* 4, 67–73.
4103 doi:10.1177/108705719900400206
4104 Zimmer, B., Pallocca, G., Dreser, N., Foerster, S., Waldmann, T., Westerhout, J., Julien, S., Krause, K.H., van
4105 Thriel, C., Hengstler, J.G., Sachinidis, A., Bosgra, S., Leist, M., 2014. Profiling of drugs and

4106
4107
4108

environmental chemicals for functional impairment of neural crest migration in a novel stem cell-based test battery. Arch. Toxicol. 88, 1109–26. doi:10.1007/s00204-014-1231-9

DRAFT

4109

List of Figures

4110 | [Figure 1: Hierarchy of governance and management documents](#) [3231](#)

4111 | [Figure 2: A typical "onion ring" structure used to separate different operational areas.](#)

4112 | [3635](#)

4113 | [Figure 3: Growth curve for cells grown in culture. Cells should be subcultured while still](#)

4114 | [in the exponential phase \(ATCC, 2014\)](#) [5150](#)

4115 | [Figure 4: Schematic representation of the processes that can cause the final target](#)

4116 | [concentration to be different than the nominal concentration in an *in vitro* test \(Kramer](#)

4117 | [et al., 2012\)](#) [6968](#)

4118 | [Figure 5: Flow chart to aid in choosing an appropriate dose metric for a specific *in vitro*](#)

4119 | [toxicity test \(Groothuis et al., 2015\)](#) [7271](#)

4120 | [Figure 6: Evolution of a Standard Operating Procedure](#) [7776](#)

4121 | [Figure 7: Example of plate layout. PC positive control, RI reference item, NC negative](#)

4122 | [control, UC untreated control, VC vehicle control, TI test item](#) [8180](#)

4123 | [Figure 8: Sample data from the Draize Eye Irritation test on two chemicals at 24 hours](#)

4124 | [after instillation](#)..... [8483](#)

4125 | [Figure 9: A 96-well based cytotoxicity assay using human keratinocytes](#) [8483](#)

4126 | [Figure 10: A quality control chart for the ethanol positive control for the BCOP showing](#)

4127 | [the acceptable range and deviations above and below \(March 2002 to November 2005-](#)

4128 | [courtesy of IIVS\)](#)..... [8685](#)

4129

Field Code Changed

DRAFT

4130

List of Tables

4131 Table 1: Examples of requirements for documentation concerning the origins of cells and
 4132 tissues [2120](#)

4133 Table 2: Examples of differences between quality assurance and quality control [2827](#)

4134 Table 3: Key criteria to be addressed for data generated under a GLP environment
 4135 (ISPE, 2012) [3332](#)

4136 Table 4: Applicability of Integrity Checks on Cultures [5453](#)

4137 Table 5: Viability testing of cell cultures [6261](#)

4138 Table 2: Possible outcomes of a test result of a chemical in a validation [9089](#)

Field Code Changed

4139

4140

DRAFT

4141 **Annex 1 Good Cell Culture Practice (GCCP)**

4142 See pdf

DRAFT

4143 **Annex 2 Good Cell Culture Practice for stem cells and stem-**
4144 **cell-derived models**

4145 http://www.altex.ch/resources/Pamies_of_160823_v4.pdf

Field Code Changed

DRAFT

4146
4147
4148

Annex 3 Experts participating at the GIVIMP meeting 24-25 February 2015

Maria Astridou	EFSA, Parma, Italy
Camilla Bernasconi	EURL ECVAM, Ispra, Italy
Gerard Bowe	EURL ECVAM, Ispra, Italy
Ann-Charlotte Bostroem	EURL ECVAM, Ispra, Italy
Sandra Coecke	EURL ECVAM, Ispra, Italy
Thomas Cole	EURL ECVAM, Ispra, Italy
Rita Cortvrindt	Vito, Mol, Belgium
Laureano Cuevas	Centro Nacional de Microbiologia, Majadahonda, Spain
Agnieszka Dudra	Bureau for Chemical Substances, Lodz, Poland
Eugene Elmore	University of California Irvine, Trabuco Canyon, Unites States
Chantra Eskes	SeCAM, Magliaso, Switzerland
Varvara Gouliarmou	EURL ECVAM, Ispra, Italy
Andrew Gray	United Kingdom GLP Monitoring Authority, The MHRA Medicines and Healthcare, London, United Kingdom
John Harbell	JHarbell Consulting LCC, Dallas, United States
Minne Heringa	RIVM, The National Institute for Public Health and the Environment, Centre for Safety of Substances and Products, Bilthoven, The Netherlands
Delyan Ivanov	University of Nottingham, Cancer Biology Queens Medical Centre, Nottingham, United Kingdom
Rob Jaspers	Health Care Inspectorate, The Netherlands
Ingrid Langezaal	EURL ECVAM, Ispra, Italy
Marcel Leist	University of Konstanz, Konstanz, Germany
Thomas Lucotte	ANSM, Saint-Denis, France
Marisa Meloni	Vitroscreen SRL, Milano, Italy
Maurits-Jan Prinz	European Commission, Directorate General Internal Market, Industry, Entrepreneurship and SMEs, Brussels, Belgium
Vittorio Reina,	EURL ECVAM, Ispra, Italy
Glyn N Stacey	United Kingdom Stem Cell Bank, Advanced Therapies Division, NIBSC, United Kingdom
Emanuela Testai	Istituto Superiore di Sanita, Rome, Italy
Frederique van Acker	TNO Triskelion, Zeist, The Netherlands
Iwona Wilk-Zasadna	Novartis Vaccines, Basel, Switzerland.
Maria-Chiara Zorzoli	EURL ECVAM, Ispra, Italy

Formatted: French (France)

4149
4150